

Molecular Studies

Equipment and Accessories

Dry bath incubator, DBI-100 Series

A combination of mixing and accurate temperature control, to make your laboratory work efficient and productive



Cat No	Functions	Temperature Range	Mixing frequency
DBI-RS100	Heating, cooling & mixing	$0.1 - 100^{\circ}C$	200rpm-1500rpm
DBI-S100	Heating & mixing	$15 - 100^{\circ}C$	200rpm-1500rpm
DBI-100III	Heating	$15 - 100^{\circ}C$	-

Heating Block Accessories

