

OBTURATION SYSTEM



C-FILL Endodontic obturation system

C-FILL a pack



Feature:

Wireless design, Reduce wrist fatigue
 Holding comfortable ,easy to use
 Heat plunger has easy -lock structure. and can rotate 360 degree
 Can be choose heating method according type of heat plunger.
 Fast heating. quick operation
 Large capacity battery. double battery charging and spare.

Technical:

Li-ion chargeable Battery: DC3.7v 2200mAh
 Heating time: 5s to 200°C
 Working temperature: 150°C,180°C. 200°C, 230°C
 Adaptor Input: AC100-240V
 Output: DV5V,1.5A
 Heat plunger : F,FM,M,ML

C-FILL β back



Features :

Wireless handpiece , easy to use
 Holding and operating very comfortable
 Exactly temperature-control and thermal protector cap can be avoiding scald
 Gun needle can be 360 degree rotation, more easy to get the fill position
 Gun needle screw design, effectively prevent the glue;
 It can be at 200 degree in 30S , can be melting all kinds of I gutta percha
 Large capacity for battery ,which can use long time
 This product has been applied layout-design patent,Coutereit must be investigated;NO ZL2014304851457

Technical:

Li-ion chargeable Battery: DC3.7v 2200mAh
 Heating time: 5s 30s to 200°C
 Working temperature: 150°C,180°C. 200°C,230°C
 Adaptor Input: AC100-240V Output: DC5V,1.5A
 Heat plunger: 23G,25G
 C-Fill set:

Volume(cm): 20.5X18X10

Weight(kg)/PCS: 2.0

Packing material: carton

QTY per master carton: 10pc/ctn



AJJ

HEALTHCARE MANAGEMENT PTE LTD