



INTRA-ORAL SCANNER



INTRA-ORAL SCANNER

intraoral scanner has fast and easy image capturing technology, which can provide reliable accurate & complete data for digital dental impression. Providing patients with a comfortable and efficient treatment experience, our scanners make it easier for dentists to diagnose and communicate with patients Runyes scanner is an open system scanner, which provides the access to the latest technology equipment available in the market. It has a range of applications which includes digital restoration, digital orthodontics navigation template of digital implant, etc.



Powderless, Real color, Light Weight

Easier to operate with net weight of 210g. Simple and elegant with a colour combination of gray and white. Looks trendy and very user friendly. Support multi-angle image capture. Hand-held design makes operation more convenient.



Constant Temperature Protection Design

Dual heating options from internal & external heater, maintain the temperature of the scanner to prevent it from fogging.



Detachable Scanner Tip

Easy to change scanner tip. Scanner tips are autoclavable and multiple times usable. Long life narrow heads for easy access to posteriors.



LED Light Source Without Radiation

High quality light source generates no radiation. Real color scans. Real imaging provided.



High Accuracy

Upto 20um accuracy.



Faster

Moves smoothly while scanning or collecting data. Quick reset to save time. No need to focus from specific distance or specific angle.



Simple Operation

One click installation. Flexible handle operation.



Open Data

Exported in multi-format, STL, PLY, OBJ.



AJJ

HEALTHCARE MANAGEMENT PTE LTD