



General Equipment Mixer Systems

Vortex Mixer, Adjustable Speed, VMX-S

The Vortex Mixers are ideal instruments for mixing liquid components in tubes, resuspension of cells, using an eccentric mechanism.

- Touch operation or continuous mode
- Variable speed control from 0 to 2500rpm
- Used for various mixing applications with optional adapters
- Specially designed vacuum suction feet for body stability
- Robust aluminum-cast construction



Enquire about your required tube adapter accessories



Roller Mixer, Roller Mixer Series

Roller mixers provide a gentle but highly efficient rocking and rolling motion, ideal for mixing blood samples, viscous substances, and liquid-solid suspensions.

- Speed range: 10rpm-80rpm
- Rocking motion: Rolling/Swinging
- Numbers of rollers: 6
- Digital speed control and LED display
- Prevent blood clotting and immunoprecipitation



RMX-TR-S

Cat No	Speed Range (rpm)	Motion	Roller Length(cm)	Dimensions (W×D×H) (cm)
RMX-TR-S	10 – 80	Rolling	28	23 x 44 x 12
RMX-TR-SR	10 - 80	Rolling, Swinging	28	23 x 44 x 12