

CHISON

EBit Series.

Portable ultrasound



EBit 60

- Lightweight at 7.5kg
- 15" LED monitor
- Rotatable from 0° - 30°
- Removable battery, 120 minutes in active mode
- Dual transducer ports (Built-in)
- Probe holders
- Q-image, Q-flow, X-contrast, Q-beam and FHI
- Advanced CW, TDI, Color m mode, ECG, PISA
- Ports: USB, LAN, VGA, DVI, Video, Remote
- Ergonomic trolley (accessory box, printer socket & probe holders, optional: theft-proof lock)
- Short-cut key function
- Multi-language

Tranducer selection for EBit 60



2.0-6.8MHz
Convex C3-E



4.0-15.0MHz
Linear L7-E



7.0-18.0MHz(With FHI)
Linear L12-E



1.5-5.3MHz
Phased Array P3-E



2.0-8.0MHz
Phased Array P6-E



2.0MHz-6.8MHz
Volume V4-EV



4.0-12.0MHz
Transvaginal V6-E



4.0-15.0MHz
Transvaginal V7-E



4.0-15.0MHz
Trans-Rectal L7R-E



2.0-6.8MHz
Micro-Convex MC3-E



4.0-12.0MHz
Micro-Convex MC6-E



4.0-10.7MHz
Micro-Convex MC5-E

EBit 50



EBit 30



*Tranducer selection differs for EBit 50 and EBit 30.