

General Consumables

Centrifuge Tubes 15ml | 50mL

Made of transparent polypropylene, the versatile tubes provide a leak proof storage of liquid samples.

- RNase- and DNase- free
- Non-pyrogenic
- Sterile
- Graduations
- Tight sealing screw cap





Cat No	Volume storage (mL)	Specifications	Max RCF (xg)	Packaging
CT-002-15A	15	Bulk Packed	15,000	25pc/bag, 20bag/ctn
CT-012-15A	15	Racked	15,000	50pc/rack, 10bag/ctn
CT-002-50A	50	Bulk Packed	15,000	25pc/bag, 20bag/ctn
CT-012-50A	50	Racked	15,000	25pc/rack, 00bag/ctn

Screw Cap Tubes

Made of transparent polypropylene, these colored-cap tubes are suitable for a wide range of temperatures from -80°C to 120°C.

- RNase- and DNase- free
- Non-pyrogenic
- · Available with conical or skirted bottom

Cat No	Volume storage (mL)	Specifications	Packaging
801205	0.5	O-ring seal, sterile	100pc/bag, 10bag/ctn
801207	1.5	O-ring seal, sterile	100pc/bag, 10bag/ctn
801209	2.0	O-ring seal, sterile	100pc/bag, 10bag/ctn
801205A	0.5	O-ring seal, sterile	100pc/bag, 10bag/ctn
801207A	1.5	O-ring seal, sterile	100pc/bag, 10bag/ctn
801209A	2.0	O-ring seal, sterile	100pc/bag, 10bag/ctn

