



# General Consumables

## Micro-centrifuge Tubes

Made of transparent polypropylene, the versatile tubes provide a leak proof storage of liquid samples.

- RNase- and DNase- free
- Nonpyrogenic
- Nonsterile tubes are autoclavable
- Graduations
- Low retention specifications available



Cat No	Volume storage(mL)	Specifications	Max RCF (xg)	Packaging
MCT-001-060	0.6	Non-Sterile	12,000	1000pc/box, 10box/ctn
MCT-001-150	1.5	Non-Sterile	12,000	500pc/box, 10box/ctn
MCT-001-150-S	1.5	Sterile	12,000	50pc/bag, 5bag/ctn
MCT-001-150X	1.5	Non-Sterile, Brown	12,000	500pc/box, 10box/ctn
MCT-001-175	1.7	Non-Sterile	12,000	500pc/box, 10box/ctn
MCT-001-200	2.0	Non-Sterile	12,000	500pc/box, 10box/ctn
MCT-001-200S	2.0	Sterile	12,000	50pc/bag, 5bag/ctn
MCT-001-500	5.0	Non-Sterile	25,000	1050pc/bag, 10bag/ctn

Enquire for safe lock designed microtubes

Friendly hinge  
for finger



Crystal clear for  
improved  
visibility



High contrast  
graduations



Large writing area  
for easy  
identification

## Low Retention Tubes

- RNase- and DNase- free
- Nonsterile
- Autoclavable (121°C, 15min)

Cat No	Volume storage(mL)	Specifications	Max RCF (xg)	Packaging
MCT-001-060LR	0.6	Non-Sterile	12,000	500pc/bag, 10bag/ctn
MCT-001-150LR	1.5	Non-Sterile	12,000	500pc/bag, 10bag/ctn
MCT-001-200LR	2.0	Sterile	12,000	500pc/bag, 10bag/ctn

