



AJJ
HEALTHCARE MANAGEMENT PTE LTD

MEDICAL PRODUCTS.

AJJ Healthcare Management Pte. Ltd. ("AJJ")
is a wholly-owned subsidiary of
OEL (Holdings) Limited (SGX:584)

Office
8 Commonwealth Lane
#02-04

Head Quarters
8 Temasek Boulevard
#23-02

Email: marketing@ajjhealthcare.com

24-Hour Hotline: +65 6908 4252

Website: www.ajjhealthcare.com

LinkedIn: [linkedin.com/company/ajj-healthcare-management-pte-ltd](https://www.linkedin.com/company/ajj-healthcare-management-pte-ltd)

Facebook: [business.facebook.com/ajjhealthcare](https://www.facebook.com/ajjhealthcare)



All for the commitment to life.





About AJJ Healthcare Management.

AJJ Healthcare Management Private Limited (“AJJ”) is a wholly-owned subsidiary of OEL (Holdings) Limited (SGX:584), and under the same values of “Fairness”, “Transparency”, “Professional” and “Integrity”, strives to create a brand of its own. AJJ’s core business focuses on high-tech, medical devices, and products of superior quality. AJJ is also customer-focused and seeks to bring the business to greater heights with its excellent services. With the attention on such key focus areas, AJJ aims to achieve the strategic goal of building a Pan-Asian, high-tech, Healthcare company, which includes establishing a medical manufacturing plant and inventory base for medical products in Singapore.

AJJ’s product portfolio consists of High-Intensity Ultrasound For Tumour, Ultrasound Diagnostic, Radiography Diagnostic, Ventilators, Respiratory Support, Anesthesia Machine, Electrocardiogram, Continuous Glucose Monitoring, Negative and Positive Pressure Isolation Booths, Consumables and Accessories

Contents

Breathcare Anesthesia.....001~016

Urinary.....017~029

Hypodermic..... 027~045


Examination.....046~050

Wound Dressing..... 051~063

Medical Non-woven & Accessories.....064~072

Others 073~079

BREATHCARE ANESTHESIA



LB1110 Simple Oxygen Mask

Made of clear, medical grade PVC
Adjustable nose clip
Anti-crush tube
2m oxygen tubing

Cat. NO.	Size
LB111000	XS, Neonate
LB111001	S, Standard child
LB111002	M, Elongated child
LB111003	L, Adult Asia
LB111004	L+, Adult Europe
LB111005	XL, Elongated adult



LB1210 Venturi Mask

Made of clear, medical grade PVC
Two color coded diluters for seven concentrations
Green diluter for 24%,26%,28%,30% White diluter for 35%,40%,50%
Adjustable nose clip
Anti-crush tube

Cat. NO.	Size
LB121000	XS, Neonate
LB121001	S, Standard child
LB121002	M, Elongated child
LB121003	L, Adult Asia
LB121004	L+, Adult Europe
LB121005	XL, Elongated adult



LB1220 Adjustable Venturi Mask

Made of clear, medical grade PVC
Single diluter for one concentration, six diluters can be selected
Adjustable nose clip
Adaptor for high humidity
Anti-crush tube
With six diluters 24% (blue) 28% (yellow)
31% (white) 35% (green) 40% (pink)
50% (orange) 60% (red)

Cat. NO.	Size
LB122000A	LB122000B XS, Neonate
LB122001A	LB122001B S, Standard child
LB122002A	LB122002B M, Elongated child
LB122003A	LB122003B L, Adult Asia
LB122004A	LB122004B L+, Adult Europe
LB122005A	LB122005B XL, Elongated adult



LB1310 Oxygen Mask with Reservoir Bag

Made of clear, medical grade PVC
Standard connector, non-rebreathing
With 1000ml/600ml reservoir bag
Anti-crush tube
2m oxygen tubing

Cat. NO.	Size
LB131001	S, Standard child
LB131002	M, Elongated child
LB131003	L, Adult Asia
LB131004	L+, Adult Europe
LB131005	XL, Elongated adult



LB1411/LB1421/LB1422 Nasal Oxygen Cannula

Made of clear, non-toxic PVC
Two tips, over ear style
Anti-crush tube
With 2m oxygen tubing

Cat. NO.	Size	Cat. NO.	Size
LB141101	Super soft, neonate	LB142103	Flared tip, pediatric
LB141102	Super soft, infant	LB142104	Flared tip, adult
LB141103	Super soft, pediatric	LB142201	Standard tip, neonate
LB141104	Super soft, adult	LB142202	Standard tip, infant
		LB142203	Standard tip, pediatric
		LB142204	Standard tip, adult



LB1610 Oxygen Connecting Tube

Made of Medical grade PVC
The length of the tube can be customized.



LB1510 Tracheostomy Mask

Specialized mask for tracheostomy
Swivel 360° adaptor with 18F/22M port
Size: M / L

Nebulizer Jar



LB20A

Volume: 6cc
Under the condition of gas flow:
7l/min, atomizing volume is a
200mg mg/min, in every 200mg,
The percentage of particulate
with the diameter less than 5um
is more than 70%.



LB20B

Volume: 6cc
Under the condition of gas flow:
7l/min, atomizing volume is a
280mg mg/min, in every 280mg,
The percentage of particulate
with the diameter less than 5um
is more than 70%.

LB20B Patent Number :
ZL 2014 2 0766313.3



LB20C

Volume: 10cc
Under the condition of gas flow:
7l/min, atomizing volume is a about
250mg mg/min, in every 250mg,
The percentage of particulate with
the diameter less than 5um is more
than 70%.

LB20C Patent Number :
ZL 201520637654.5



LB20D

Volume: 6cc
Under the condition of gas flow:
7l/min, atomizing volume is a
about 150-450mg mg/min,
The percentage of particulate
with the diameter less than 5um
is more than 70%.

LB20D Patent Number :
ZL 201520289256.9



LB2011 Nebulizer Mask

Operates even at horizontal position Capacity: 6cc
disinfected with boiling water (only for jar)
2m oxygen tubing

Different type of nebulizer jars are available.

Cat. NO.	Size
LB201100	XS, Neonate
LB201101	Standard child
LB201102	Elongated child
LB201103	Adult Asia
LB201104	Adult Europe
LB201105	Elongated adult



LB2012 Nebulizer Kit

With 1pc adult mask and 1pc child mask
T connector, mouthpiece, 2m oxygen tubing
Different type of nebulizer jars are available.



LB2013 Nebulizer with Mouth Piece

Hand-held, with mouthpiece
T connector, mouthpiece, 2m oxygen tubing
Transparent
With/without corrugated tube
Different type of nebulizer jars are available.





LB60 Suction Connecting Tubing

- 1) Made of transparent material for better visualization
- 2) Striated walls of tube provide superior strength and anti-kinking
- 3) Multiple choices of length
- 4) Tube Length Can be Customization

Cat. NO.	Length (cm)	Size (inch)	Type
LB601118	180	1/4	Long connector
LB601218	180	1/4	Short connector
LB601318	180	1/4	Trumpet connector
LB602118	180	3/16	Long connector
LB602218	180	3/16	Short connector
LB602318	180	3/16	Trumpet connector



LB61 Suction Connecting Tubing with yankauer handle

- 1) Made of transparent material
- 2) Can be combined with different types of yankauer handle and connecting tube to suit various requirement
- 3) Multiple choices of length
- 4) Can choose different connector

Cat. NO.	Length (cm)	Size (inch)	Type
LB6111118	180	1/4	Crown tip with vent
LB6112118	180	1/4	Crown tip without vent
LB6113118	180	1/4	Flat tip with vent
LB6114118	180	1/4	Flat tip without vent
LB6121118	180	3/16	Crown tip with vent
LB6122118	180	3/16	Crown tip without vent
LB6123118	180	3/16	Flat tip with vent
LB6124118	180	3/16	Flat tip without vent

LB3011

Anesthesia Mask

Made from Latex/DEHP free PVC material.
Air cushion assures comfortable face fitting.
Tail Check valve/Upright Check valve
Without check valve type are available.
Single Patient use.
Color-Coded hook options.

Cat.No.	Size	Color	Connector
LB301100	size0# - Neonate	N/A	15mm OD
LB301101	size1# - Infant	N/A	15mm OD
LB301102	size2# - Pediatric	Yellow	22mm ID
LB301103	size 3# - Adult-small	Green	22mm ID
LB301104	size 4# - Adult-medium	Red	22mm ID
LB301105	size 5# - Adult-large	Blue	22mm ID
LB301106	size 6# - Adult-Extra large	Orange	22mm ID



LB3012

PVC free Anesthesia Mask

Made from TPE and PP Material, contain no PVC and are DEHP FREE
The products provides a comfortable superior airtight fit with a clear shell to enhance patient face visibility.
This products ancillary use with anesthesia breathing circuit.
Non-reusable, single use only

Cat.No.	Size	Color	Connector
LB301201	Size1# - Neonate	Blue	15mm OD
LB301202	Size2# - Infant	Orange	15mm OD
LB301203	Size3# - Pediatric	Pink	22mm ID
LB301204	Size4# - Adult-S	White	22mm ID
LB301205	size5# - Adult-M	Yellow	22mm ID
LB301206	size6# - Adult-L	Green	22mm ID



LB3021

Soft Anesthesia Mask

Made from Latex/DEHP free PVC material.
Soft adjustable face cushion assures comfortable face fitting.
With Upright Check valve.
Single Patient use.
No scents or with Vanilla scents are available.

Cat.No.	Size	Color	Connector
LB3021000	Size 000# - Premature	White	15mm OD
LB302100	Size0# - Neonate	Pink	15mm OD
LB302101	Size1# - Infant	Yellow	15mm OD
LB302102	Size2# - Pediatric	Green	22mm ID
LB302103	Size3# - Adult-Small	Red	22mm ID
LB302104	size4# - Adult-Medium	Blue	22mm ID
LB302105	size5# - Adult-Large	Orange	22mm ID



LB3030

Silicon Anesthesia Mask (One-piece)

All-in-one piece soft silicone cushion for better seal and comfort.
Ergonomic shape for lower dead space and better air circulation.
Made of 100% clear, medical grade silicon, latex free
Reusable, all parts are fully autoclavable up to 134°C
Packing: bulk or pe pack available

Cat.No.	Size	Connector
LB303000	size0# - Neonate	15mm OD
LB303001	size1# - Infant	15mm OD
LB303002	size2# - Child	22mm ID
LB303003	size 3# - Adult-small	22mm ID
LB303004	size 4# - Adult-medium	22mm ID
LB303005	size 5# - Adult-large	22mm ID



Anesthesia Breathing Circuit

Anesthesia breathing circuit is bridge of anesthesia machine with Patient, input anesthesia gases and output expiratory air from the patient, to accomplish the normal exchange of oxygen and carbon dioxide gas.

Disposable anesthesia circuit can prevent patients infecting, establish manual airway match with anesthesia machine and respirator. Non-toxic, insipid, made of medical degree material, flexible transparent, air tightness.

LB4511



Anesthesia Circuit-Expandable (Adult/Pediatric)

LB4512



Anesthesia Circuit-Reinforced (Adult/Pediatric)

LB431R



Anesthesia Circuit-Reinforced (Adult/Pediatric)

LB4330



Anesthesia Circuit-Coaxial (Adult)

LB432



Breathing Circuit-Corrugated (Adult/Pediatric)

LB432R



Breathing Circuit-Reinforced (Adult/Pediatric)

LB4410E

LB4410

LB4410R



Catheter Mount (Normal/Expandable/Reinforced)

LB4610



Water Trap

Anesthesia Breathing Bag

LB4910



Anesthesia Breathing bag-Latex Free
Specification: 0.5L, 1L, 2L, 3L

LB4920



Anesthesia Breathing bag-Latex
Specification: 0.5L, 1L, 2L, 3L

Filter

Disposal breathing filter made of medical degree material, mainly used for filtering bacteria and increase gasses humidity degree, effectively prevent patient infecting and mainly used as medical accessory. Efficiency filtration, humidity exchange function.



LB411001 BV Filter



LB411002 BV Filter



LB421001 HME Filter



LB421002 HME Filter



LB421003 HME Filter



LB421004 HME Filter



LB421005 HME Filter



LB421006 HME Filter



LB3110

LB3110 Silicone Resuscitator

Built-in 40-60cm H₂O pressure relief valve.
The set includes oxygen reservoir bag, silicone mask and oxygen tubing.
100% Medical level silicone material.
Latex-free components.
All-in-one intake valve connects to oxygen reservoir bag directly.
Additional accessories and customized label are available.
Built-in 30MM exhale port for PEEP or filter attachment.
All parts (excluding the oxygen tubing & oxygen reservoir) are fully autoclavable up to 134°C.



LB3110L

LB3110L Liquid Silicone Resuscitator

Built-in 40-60cm H₂O pressure relief valve
The set includes oxygen reservoir bag, silicone mask and oxygen tubing
100% Medical level silicone material
Latex-free components
All-in-one intake valve connects to oxygen reservoir bag directly
Additional accessories and customized label are available
Built-in 30mm exhale port for PEEP or filter attachment
All parts (excluding the oxygen tubing & oxygen reservoir) are fully autoclavable up to 134°C
Transparency and resilience is better



LB3120

LB3120 SEBS Resuscitator

Single patient use only.
Built-in 60cm H₂O pressure relief valve.
The set includes oxygen reservoir bag, PVC and oxygen tubing.
100% Medical level SEBS material.
Latex-free components.
All-in-one intake valve connects to oxygen reservoir bag directly.
Additional accessories and customized label are available.
Built-in 30MM exhale port for PEEP or filter attachment.



LB3130

LB3130 PVC Resuscitator

Single patient use only.
Built-in 60cm H₂O pressure relief valve.
The set includes oxygen reservoir bag, PVC mask and oxygen tubing.
100% Medical level PVC material.
Latex-free components.
All-in-one intake valve connects to oxygen reservoir bag directly.
Additional accessories and customized label are available.
Built-in 30M exhale port for PEEP or filter attachment.

LB3150 CPR Mask

Features:

Non-reusable, single use only.
Foldable and non-foldable mask design.
High transparency for better visualization.
100% Medical level PVC material.
Latex-free components.

Specification:

Foldable cushion mask with oxygen port.
One way valve with happy efficiency filter.
Elastic head trap.
Individual colorful hard case package



LB3210 Disposable Silicone Laryngeal Mask

Medical grade silicone for super compatibility and patient comfort
Available with seven sizes for all patient use

Cat. NO.	Type	Applicable People
LB321010	1.0	< 5kg
LB321015	1.5	5-10kg
LB321020	2.0	10-20kg
LB321025	2.5	20-30kg
LB321030	3.0	30-50kg
LB321040	4.0	50-70kg
LB321050	5.0	> 70kg



LB3220 Reusable Silicone Laryngeal Mask

Reusable Silicone Laryngeal Mask are designed to be sterilized and used up to 40 times
Durable high temperature resistant material suitable for autoclave sterilization
Additional user card for record of device usage times and patient information is supplied in each pouch

Cat. NO.	Type	Applicable People
LB322010	1.0	< 5kg
LB322015	1.5	5-10kg
LB322020	2.0	10-20kg
LB322025	2.5	20-30kg
LB322030	3.0	30-50kg
LB322040	4.0	50-70kg
LB322050	5.0	> 70kg



LB3230 Disposable PVC Laryngeal Mask

Medical grade PVC with good compatibility
Available with seven sizes for all patient use

Cat. NO.	Type	Applicable People
LB323010	1.0	< 5kg
LB323015	1.5	5-10kg
LB323020	2.0	10-20kg
LB323025	2.5	20-30kg
LB323030	3.0	30-50kg
LB323040	4.0	50-70kg
LB323050	5.0	> 70kg



Oropharyngeal Airways

Smooth integrated design for optimum patient comfort and safety Half hard non-toxic polyethylene material with color code

LB3310 Guedel Type

Cat. NO.	Size	Color Code
LB3310040	000#, 40mm	Pink
LB3310050	00#, 50 mm	Blue
LB3310060	0#, 60 mm	Black
LB3310070	1#, 70 mm	White
LB3310080	2#, 80 mm	Green
LB3310090	3#, 90 mm	Yellow
LB3310100	4#, 100mm	Red
LB3310110	5#, 110 mm	Light Blue
LB3310120	6#, 120mm	Orange



LB3320 Berman Type

Cat. NO.	Size
LB3320040	000#, 40mm
LB3320050	00#, 50 mm
LB3320060	0#, 60 mm
LB3320070	1#, 70 mm
LB3320080	2#, 80 mm
LB3320090	3#, 90 mm
LB3320100	4#, 100mm
LB3320110	5#, 110 mm



LB3410 Nasopharyngeal Airway

Made form non-toxic PVC, soft and comfortable Only for nasal airway

Cat. NO.	Size
LB3410030	3.0 mm
LB3410035	3.5 mm
LB3410040	4.0 mm
LB3410045	4.5 mm
LB3410050	5.0 mm
LB3410055	5.5 mm
LB3410060	6.0 mm
LB3410065	6.5 mm
LB3410070	7.0 mm
LB3410075	7.5 mm
LB3410080	8.0 mm
LB3410085	8.5 mm
LB3410090	9.0 mm
LB3410095	9.5 mm



Endotracheal Tube

X-ray opaque line throughout the length of the tube Non-toxic PVC, transparent, soft and smooth

Standard Endotracheal Tube

LB5010 Uncuffed

Cat. NO.	Size ID (mm)
LB5010020	2.0
LB5010025	2.5
LB5010030	3.0
LB5010035	3.5
LB5010040	4.0
LB5010045	4.5
LB5010050	5.0
LB5010055	5.5
LB5010060	6.0

Cat. NO.	Size ID (mm)
LB5010070	7.0
LB5010075	7.5
LB5010080	8.0
LB5010085	8.5
LB5010090	9.0
LB5010095	9.5
LB5010100	10.0



LB5010C High volume low pressure cuff

Cat. NO.	Size ID (mm)
LB5010025C	2.5
LB5010030C	3.0
LB5010035C	3.5
LB5010040C	4.0
LB5010045C	4.5
LB5010050C	5.0
LB5010055C	5.5
LB5010060C	6.0

Cat. NO.	Size ID (mm)
LB5010065C	6.5
LB5010070C	7.0
LB5010075C	7.5
LB5010080C	8.0
LB5010085C	8.5
LB5010090C	9.0
LB5010095C	9.5
LB5010100C	10.0



Oral Preformed Endotracheal Tube

LB5020 Uncuffed

Cat. NO.	Size ID (mm)
LB5020020	2.0
LB5020025	2.5
LB5020030	3.0
LB5020035	3.5
LB5020040	4.0
LB5020045	4.5
LB5020050	5.0
LB5020055	5.5
LB5020060	6.0

Cat. NO.	Size ID (mm)
LB5020065	6.5
LB5020070	7.0
LB5020075	7.5
LB5020080	8.0
LB5020085	8.5
LB5020090	9.0
LB5020095	9.5
LB5020100	10.0



LB5020C High volume low pressure cuff

Cat. NO.	Size ID (mm)
LB5020030C	3.0
LB5020035C	3.5
LB5020040C	4.0
LB5020045C	4.5
LB5020050C	5.0
LB5020055C	5.5
LB5020060C	6.0
LB5020065C	6.5

Cat. NO.	Size ID (mm)
LB5020070C	7.0
LB5020075C	7.5
LB5020080C	8.0
LB5020085C	8.5
LB5020090C	9.0
LB5020095C	9.5
LB5020100C	10.0



Endotracheal Tube

X-ray opaque line throughout the length of the tube
Non-toxic PVC, transparent, soft and smooth

Nasal Preformed Endotracheal Tube

LB5030 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5030020	2.0	LB5030065	6.5
LB5030030	3.0	LB5030070	7.0
LB5030035	3.5	LB5030075	7.5
LB5030040	4.0	LB5030080	8.0
LB5030045	4.5	LB5030085	8.5
LB5030050	5.0	LB5030090	9.0
LB5030055	5.5	LB5030095	9.5
LB5030060	6.0	LB5030100	10.0



LB5030C High volume low pressure cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5030030C	3.0	LB5030065C	6.5
LB5030035C	3.5	LB5030070C	7.0
LB5030040C	4.0	LB5030075C	7.5
LB5030045C	4.5	LB5030080C	8.0
LB5030050C	5.0	LB5030085C	8.5
LB5030055C	5.5	LB5030090C	9.0
LB5030060C	6.0	LB5030095C	9.5
		LB5030100C	10.0



Reinforced Endotracheal Tube

LB5040 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5040030	3.0	LB5040065	6.5
LB5040035	3.5	LB5040070	7.0
LB5040040	4.0	LB5040075	7.5
LB5040045	4.5	LB5040080	8.0
LB5040050	5.0	LB5040085	8.5
LB5040055	5.5	LB5040090	9.0
LB5040060	6.0	LB5040095	9.5



LB5040C High volume low pressure cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5040030C	3.0	LB5040065C	6.5
LB5040035C	3.5	LB5040070C	7.0
LB5040040C	4.0	LB5040075C	7.5
LB5040045C	4.5	LB5040080C	8.0
LB5040050C	5.0	LB5040085C	8.5
LB5040055C	5.5	LB5040090C	9.0
LB5040060C	6.0	LB5040095C	9.5



Tracheostomy Tube

Made from non-toxic PVC
Used for tracheostomy



Tracheostomy tube

LB5110 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5110030	3.0	LB5110070	7.0
LB5110035	3.5	LB5110075	7.5
LB5110040	4.0	LB5110080	8.0
LB5110045	4.5	LB5110085	8.5
LB5110050	5.0	LB5110090	9.0
LB5110055	5.5	LB5110095	9.5
LB5110060	6.0	LB5110100	10.0
LB5110065	6.5		

LB5110C High volume low pressure cuff


Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5110030C	3.0	LB5110070C	7.0
LB5110035C	3.5	LB5110075C	7.5
LB5110040C	4.0	LB5110080C	8.0
LB5110045C	4.5	LB5110085C	8.5
LB5110050C	5.0	LB5110090C	9.0
LB5110055C	5.5	LB5110095C	9.5
LB5110060C	6.0	LB5110100C	10.0
LB5110065C	6.5		

FY0636 Endotracheal Intubation Set

Specification:
Endotracheal tube with stylet 1pc
Suction catheter FR.14 with finger controller tip 1pc
Suction catheter FR. 16 with finger controller tip 1pc
Gloves 1pair
Guedel airway 1pc
Mouth pad 1pc
Bite block 1pc
Syringe 5cc
Gauze swab 2pcs
Suction connecting tube 1pc
Facemask 1pc
Cotton balls with silicone oil 2pcs

Endotracheal Tube Introducer (Bougie)

Smooth elbow design of tip can assist intubation
decreasing the damage of soft tissue
scale mark indicates distance to the tip
recommend for difficult intubation and replacing
an endotracheal tube
Available with solid type and hollow type.
Hollow type can facilitate delivery of oxygen in emergency



Cat. NO.	Size (Fr)	Type	Applicable ETT
LB521006	6	Solid	< 3.5mm
LB521010	10	Solid	4.0-5.0mm
LB521015	15	Solid	> 5.5mm
LB522010	10	Hollow	4.0-5.0mm
LB522015	15	Hollow	> 5.5mm



LB6200 Suction Catheter

Made from non-toxic PVC, non-irritant soft PVC

Product Description :

Tube is frosted or smooth
With two side hole
The head end is smooth, openings
With X-ray or without X-ray

Choice of connectors :

LB6210 Plain type


LB6220 Cap-cone

LB6230 Y-type

LB6240 T-type

Specification : Fr5-Fr20

Standard Tube Length : 50cm



LB6310 Feeding Tube

Made from non-toxic PVC, non-irritant soft PVC

Product Description :

Tube is frosted or smooth
With two side hole
The head end is smooth, openings or closed
With X-ray or without X-ray
Scale line

Specification: Fr4-Fr8

Standard Tube Length: 40/45cm

Specification: Fr10-Fr22

Standard Tube Length: 120cm



LB6410 Stomach Tube

Made from non-toxic PVC
non-irritant soft PVC or 100% silicone

Product Description :

Tube is Frosted or smooth
With four side hole
The head end is smooth, openings or closed
With X-ray or without X-ray
Scale line

Specification: Fr5-Fr22

Standard Tube Length: 120cm



LB6510 Rectal Tube

Made from non-toxic PVC , non-irritant soft PVC

Product Description :

Tube is frosted or smooth
With two side hole
The head end is smooth, closed

Specification: Fr20-Fr36

Standard Tube Length: 40cm



LB6600 Nelaton Catheter (PVC)

Made from non-toxic PVC , non-irritant soft PVC

Product Description :

Tube is Play frosted or smooth
With two side hole
The head end is smooth, closed

Specification: Fr6-Fr18

Standard Tube Length: **LB6610** Male 40cm

LB6620 Female 20cm



LB6122 Extractor Mucus

Made from non-toxic PVC

Capacity: 25ml

Specification: Fr6-Fr14


Standard Tube Length: 35cm



LB6810 Closed Suction Catheter

Sputum suction pipe sealing film produced by imported resin, strong toughness, feel is exquisite and conducive to health care workers
Sputum suction tube body design scale mark line, facilitate throughout the sputum suction effect
Closed suction phlegm tube adopts full closed type design, can be used to sync with breathing machine, at the time of sputum suction oxygen to sustainable, convenient and safe operation, to avoid cross-infection
Specification : Fr6-Fr18





Urethral Catheter

- Made from Latex, silicone coated
- Color-coded for size identification
- Sterile by EO
- For single use only
- Hard & soft valves are available

FY0101 2-Way Latex Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
FY010106	6	3	Light Red	27
FY010108	8	3-5	Black	27
FY010110	10	3-5	Grey	27
FY010112	12	5-10/30	White	40
FY010114	14	5-10/30	Green	40
FY010116	16	5-10/30	Orange	40
FY010118	18	5-10/30	Red	40
FY010120	20	5-10/30	Yellow	40
FY010122	22	5-10/30	Violet	40
FY010124	24	30	Blue	40
FY010126	26	30	Pink	40



FY0103 2-Way Female Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
FY010312	12	5-10/30	White	27
FY010314	14	5-10/30	Green	27
FY010316	16	5-10/30	Orange	27
FY010318	18	5-10/30	Red	27
FY010320	20	5-10/30	Yellow	27
FY010322	22	5-10/30	Violet	27



FY0104 2-Way Tiemann Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
Fy010412	12	5-10/30	White	40
Fy010414	14	5-10/30	Green	40
Fy010416	16	5-10/30	Orange	40
Fy010418	18	5-10/30	Red	40
Fy010420	20	5-10/30	Yellow	40
Fy010422	22	5-10/30	Violet	40



FY0102 3-Way Latex Foley Catheter


Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
FY010216	16	30	Orange	40
FY010218	18	30	Red	40
FY010220	20	30	Yellow	40
FY010222	22	30	Violet	40
FY010224	24	30	Blue	40
FY010226	26	30	Pink	40





FY0105 2-Way All Silicone Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
FY010506	6	3	Light Red	31
FY010508	8	3-5	Black	31
FY010510	10	3-5	Grey	31
FY010512	12	5-10/30	White	40
FY010514	14	5-10/30	Green	40
FY010516	16	30	Orange	40
FY010518	18	30	Red	40
FY010520	20	30	Yellow	40
FY010522	22	30	Violet	40
FY010524	24	30	Blue	40
FY010526	26	30	Pink	40



FY0106 3-Way Silicone Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
FY010616	16	30	Orange	40
FY010618	18	30	Red	40
FY010620	20	30	Yellow	40
FY010622	22	30	Violet	40
FY010624	24	30	Blue	40
FY010626	26	30	Pink	40

FY0107 Urethral Catheter (Red Latex)

Cat. NO.	Size (Fr/Ch)	Length(cm)
FY010706	6	27
FY010708	8	27
FY010710	10	27
FY010712	12	40
FY010714	14	40
FY010716	16	40
FY010718	18	40
FY010720	20	40
FY010722	22	40
FY010724	24	40

FY0108 Malecot Catheter

Cat. NO.	Size (Fr/Ch)	Length(cm)
FY010812	12	40
FY010814	14	40
FY010816	16	40
FY010818	18	40
FU010820	20	40
FY010822	22	40
FY010824	24	40
FY010826	26	40
FY010828	28	40
FY010830	30	40



FY0110 Latex T Tube

Made of natural latex , and silicone coated
Used for open gallbladder
The stem is often connected to a bile bag
Size : Fr10-Fr26



FY0111 Condom Catheter

Made of natural latex
By connecting pipe, fixture, membrane, membrane set of edge and sticky
Size : 25mm , 30mm , 35mm , 40mm



FY0121 Silicone Condom Catheter

Made of 100% imported medical grade silicone
Size : 25mm , 30mm , 35mm , 40mm
Sterile by EO



**FY0201****Economic Urine Bag**

Capacity: 1000/1500ml/2000ml
without bottom outlet
Adapter with cap
Anti-Reflux Flutter Valve
Inlet tube: O.D 6.5mm; 90cm Length
Sterilization: EO
Packing: Individual PE Bag, 10pcs/Inner Bag, 250pcs/Ctn

**FY0202****Economic Urine Bag**

Capacity: 2000ml
Push Pull Valve
Adapter with cap
Anti-Reflux Flutter Valve
Inlet tube: O.D 6.5mm; 90cm Length
Sterilization: EO
Packing: Individual PE Bag, 10pcs/Inner Bag, 250pcs/Ctn

**FY0203****Economic Urine Bag**

Capacity: 2000ml
Twist Turn Valve
With bottom outlet
Adapter with cap
Anti-Reflux Flutter Valve
Inlet tube: O.D 6.5mm; 90cm Length
Sterilization: EO
Packing: Individual PE Bag, 10pcs/Inner Bag, 250pcs/Ctn

**FY0204****Economic Urine Bag**

Capacity: 2000ml
T-Tap Valve
with bottom outlet
Adapter with cap
Anti-Reflux Flutter Valve
Inlet tube: O.D 6.5mm; 90cm Length
Sterilization: EO
Packing: Individual PE Bag, 10pcs/Inner Bag, 250pcs/Ctn

FY0205**Urinary Drainage Bag**

Capacity: 2000ml
Twist Turn Valve
With ribbon
Anti-Reflux Flutter Valve
Inlet tube: O.D 9mm; 100cm Length
Sterilization: EO
Packing: Individual PE Bag, 100pcs/Ctn

**FY0206****Urinary Drainage Bag**

Capacity: 2000ml
T- Valve
With ribbon
Anti-Reflux Flutter Valve
Inlet tube: O.D 9mm; 100cm Length
Needleless Sampling Port with Cap
Sterilization: EO
Packing: Individual PE Bag, 100pcs/Ctn

**FY0207****Urinary Drainage Bag**

Capacity: 2000ml
T- Valve
Anti-Reflux Flutter Valve
Inlet tube: O.D 9mm; 100cm Length
Needleless Sampling Port with Cap
With bed sheet clamp
With universal double-hook
Sterilization: EO
Packing: Individual PE Bag, 100pcs/Ctn

**FY0208****Urinary Drainage Bag**

Capacity: 2000ml
With Single Hanger
Anti-Reflux Flutter Valve
Inlet tube: O.D 10mm; 100cm Length
Needleless Sampling Port with Cap
With Slide Clamp
Silicone Outlet Tube with Plastic Folding Clamp
Sterilization: EO
Packing: Individual PE Bag, 40pcs/Ctn



**FY0209****Urinary Drainage Bag**

Capacity: 2000ml
 Twist Turn Valve
 With ribbon
 Anti-Reflux Flutter Valve
 Inlet tube: O.D 9mm; 100cm Length
 Adapter with cap
 Sterilization: EO
 Packing: Individual PE Bag, 100pcs/Ctn

**FY0210****Urinary Drainage Bag**

Capacity: 2000ml
 with bed sheet clamp
 Anti-Reflux Flutter Valve, Air-Vent on the Bag
 Inlet tube: O.D 9mm; 100cm Length
 Needleless Sampling Port with Cap
 With Single Hanger and Rope Hanger
 Silicone Outlet Tube with Metal Folding Clamp
 Sterilization: EO
 Packing: Individual blister Packing, 40pcs/Ctn

**FY0211****Urinary Drainage Bag**

Capacity: 2000ml
 T-Tap Outlet Valve, Rigid Dog House
 With bed sheet clamp, with Slide Clamp
 Inlet tube: O.D 10mm; 110cm Length
 Needleless Sampling Port with Cap
 Three Parts Anti-Reflux Drip Chamber
 Air-Vent on the Bag
 Universal Double-Hook and Rope Hanger
 Sterilization: EO
 Packing: Individual blister Packing, 40pcs/Ctn

**FY0212****Urinary Drainage Bag**

Capacity: 2000ml
 T-Tap Outlet Valve
 With plastic cover, with bed sheet clamp, with Slide Clamp
 Anti-Reflux Drip Chamber, air-Vent on the Bag
 Inlet tube: O.D 10mm; 120cm Length
 Needleless Sampling Port with Cap
 Universal Double-Hook
 Sterilization: EO
 Packing: Individual blister Packing, 40pcs/Ctn

FY0213**Urinary Drainage Bag**

Capacity: 2000ml
 T-Tap Outlet Valve, Rigid Dog House
 Inlet tube: O.D 10mm; 120cm Length
 Needleless Sampling Port with Cap
 with bed sheet clamp, with Slide Clamp
 Anti-Reflux Drip Chamber, Air-Vent on the Bag
 Universal Double-Hook and Rope Hanger
 Sterilization: EO
 Packing: Individual blister Packing, 40pcs/Ctn

**FY0214****Urinary Drainage Bag**

Capacity: 2000ml/2500ml
 T-Tap Outlet Valve, Rigid Dog House
 Inlet tube: O.D 10mm; 120cm Length
 Needleless Sampling Port with Cap
 with bed sheet clamp, with Slide Clamp
 Three Parts Anti-Reflux Drip Chamber
 air-Vent on the Bag
 Universal Double-Hook
 Sterilization: EO
 Packing: Individual blister Packing, 40pcs/Ctn

**FY0215****Urinary Drainage Bag**

Capacity: 4000ml/5000ml
 T-Tap Outlet Valve, Rigid Dog House
 Inlet tube: O.D 10mm; 120cm Length
 Needleless Sampling Port with Cap
 with bed sheet clamp, with Slide Clamp
 Two Parts Anti-Reflux Drip Chamber
 air-Vent on the Bag
 Universal Double-Hook and Rope Hanger
 Sterilization: EO
 Packing: Individual blister Packing, 40pcs/Ctn

**FY0216****Urinary Drainage Bag**

Capacity: 2000ml
 T-Tap Outlet Valve
 Inlet tube: O.D 10mm; 120cm Length
 Needleless Sampling Port with Cap
 with bed sheet clamp, with Slide Clamp
 Two Parts Anti-Reflux Drip Chamber
 air-Vent on the Bag
 Sterilization: EO
 Packing: Individual blister Packing, 40pcs/Ctn



**FY0217****Urinary Drainage Urine Meter**

400cc Two-Chamber Urine meter, 2000ml Urine Bag
 Inlet tube: O.D 10mm; 120cm Length
 Needleless Sampling Port with Cap
 with bed sheet clamp, with Slide Clamp
 Universal Double-Hook and Rope Hanger
 Silicone Outlet Tube with Metal Folding Clamp
 and Splash Guard
 Sterilization:EO
 Packing: Individual blister Packing, 10pcs/Ctn

**FY0218****Urinary Drainage Meter**

400cc Two-Chamber Urine meter, 2500ml Urine Bag
 Inlet tube: O.D 10mm; 120cm Length
 Needleless Sampling Port with Cap
 with bed sheet clamp, with Slide Clamp
 Universal Double-Hook
 T-Tap Outlet Valve, with plastic cover
 Sterilization: EO
 Packing: Individual blister Packing, 10pcs/Ctn

**FY0219****Urinary Drainage Urine Meter**

200cc Urinometer with Drip Chamber
 2000ml Urine Bag
 Inlet tube: O.D 10mm; 120cm Length
 Needleless Sampling Port with Cap
 with bed sheet clamp, with Slide Clamp
 Universal Double-Hook and Rope Hanger
 T-Tap Outlet Valve
 Sterilization: EO
 Packing: Individual blister Packing, 10pcs/Ctn

**FY0220****Urinary Drainage Urine Meter**

500cc Urine meter with 2000ml Urine Bag
 Inlet tube: O.D 10mm; 120cm Length
 Drip Chamber with Air-Vent and Anti-Reflux Valve
 Needleless Sampling Port with Cap, Plastic Hanger
 bed Sheet Clamp, with Slide Clamp
 T-Tap Outlet Valve
 Sterilization: EO
 Packing: Individual blister Packing, 10pcs/Ctn

FY0221**Urinary Drainage Leg Bag**

500ml Leg Bag
 Anti-Reflux Flutter Valve
 Twist Valve
 Two Pre-Assembled Straps
 Sterilization: EO
 Packing: Individual PE Bag, 48pcs/Ctn

**FY0222****Urinary Drainage Leg Bag**

500 ml, non-woven fabric on the back
 Small and white T valve plug
 Inlet tube: O.D 9mm; 30cm Length
 Silicone tube 4.5 cm
 Sterilization: EO
 Packing: Individual blister Packing, 100 pcs/Ctn

**FY0223****Urinary Drainage Leg Bag**

750 ml, non-woven fabric on the back, three cavity
 Small and white T valve plug,
 Inlet tube: O.D 9mm; 30cm Length
 Needleless sampling port with cap
 Silicone tube 4.5cm
 Two straps are available
 Sterilization: EO
 Packing: Individual blister Packing, 100 pcs/Ctn

**FY0224****Urinary Drainage Leg Bag**

750 ml, non-woven fabric on the back
 Small and white T valve plug,
 Inlet tube: O.D 9mm; 30cm Length
 Needleless sampling port with cap
 Two Pre-assembled Straps
 Sterilization: EO
 Packing: Individual blister Packing, 100 pcs/Ctn





FY0229 Pediatric Urinary Collection Bag

Safety Adhesive,PVC Bag , Sponge
Capacity: 100ml, 200ml
Sterilization: EO
Packing: Individual PE Bag, 100 pcs/Ininer Box, 5000Pcs/ctn

FY0235 Closed Wound Drainage System Kit

Made of imported medical grade silicone
This product includes 100cc/200cc/300cc/400cc Silicone Reservoir
7mm/10mm Flat drains or other round drains
and stainless steel trocar
Can use separately or using Trocar for insertion and
combining with silicone Reservoir for fluids syction
and collection
Used for various surgeries
Latex-free
Sterilization: EO
Packing: Individual blister Packing, 20pcs/Ctn




FY0230 Closed Wound Drainage System (spring) FY0231 Closed Wound Drainage System (hollow)

Our closed wound drainage system,including 400cc, 3-spring evacuator,PVC tubing and "Y"connector,and drain with trocar.The 3-spring evacuator design is Minimize Collapse-Resistance,Packed in double liner Tyvek, to make sure mdeical cleanness on O.R.operation.
Capacity: Spring: 200ml, 400ml
Hollow: 100ml, 200ml, 400ml
Silicone or PVC tube, 120cm Length
Sterilization: EO
Packing: Individual blister Packing, 20pcs/Ctn



FY0234 Enteral Feeding Bag Set

Volume: 1000ml, 1200ml
Gravity Type/Pump Type
Enteral Feeding Bag has rigid neck
for easy filling and handling.
Strong,dependable hanging ring.
Easy-to-read graduations
and easy-view translucent bag.
Bottom exit port allows complete drainage
Sterilization: EO
Packing: Individual PE bag,24pcs/ctn



FY0241 Disposable use drainage bottle equipment

Single bottle, 1600ml
Double bottle, 2000ml
Three bottles, 2500ml

Constitutiong:

By attachment,drive pipe,bottle cap,bottle body,dip pipe,
support and suspension hook constitution.



FY0301
Simple Type

FY0302
Foam Type

FY0303
Glue Pan Type



FY0304
Two-piece Type

FY0305
Foam Type

FY0306
Two-piece Type

FY030 New Style Ostomy Bag

Our ostomy bag are of excellent quality and advanced usage,making the wound healed in the good condition
Ostomy bag applications:collect excrement and urine after colostomy, ileostomy and urostomy.

Features:

- 1,The hydrocolloid skin barrier has the performance o anti-maceration, low allergy,and warping and leakage. It has excellent exduation-absorbent and long sticking ability(up to 7 days).And it is soft and comfortable.
- 2,Hydrocolloid skin barrier combines with the highly permeable PU film makes it stick for a long time, no leakage of exudates,convenient to use.
- 3,The material of ostomy bag possesses high obstructing ability,it is comfortable and covert.
- 4,The one-piece system desg n discharge interstinal gas easily and filtering effect it excellent.

Simple type Packing reference:

Foam type 10pcs/bag,400pcs/ctn

Glue Pan type Size:35x31x40cm

Two-piece type G.W.:6kg N.W.:4.5kg

EVA TPN (total parenteral nutrition) Bag

Features:

- Non-PVC material, DEHP free, can be used to fat emulsion and other fat soluble liquid and medicine.
- No adverse effects on infants and young children and pregnant women.
- Flexible tube and bag to avoid kinking and break.



Item	Volume	Details
	500ml	Eva bag body Eva inlet tube
	1000ml	
	2000ml	
	3000ml	
	500ml	Eva bag body
	1000ml	
	2000ml	
	3000ml	

FY0247
Suction Liner

Product Characteristic:
Simplicity
Safety
Economical

Indication for Use:
This product is used for collecting the biohazardous fluids which produced in ICU and operation room.

- Product Features:**
- 1.Double safety design to ensure no liquid is suctioned back to the tubing and suction unit, effectively prevent cross infection and protect the medical staff safe.
 - 2.It's made of environment-friendly polyethylene. It's safe and non-toxic.
 - 3.Easy to use, can be a single use, also can use with the stander.
 - 4.Available in some specifications: with or without shutoff filter, with or without solidifier, 1000mlml,1500ml and 2000ml.

- Product Advantages:**
- 1.Soft and transparent(Convenient for observation liquid waste)
 - 2.Double safety design(To prevent fluid into the suction systems by mistake)
 - 3.It can be connected with relevant suction device or vacuum device, soft bag is seal and reliable inside and can be withstand the pressure of 60kg when fully filled with liquid(Prevent infections liquid overflow caused by medical personnel) .



Specification

Item	Description
Stand	Four tandem stand 4 canisters (a set)
	Twin tandem stand 2 canisters (a set)
Hanger	Twao tandem hanger stand
	Single hanger stand
	Single canister





FY0601 **Disposable Hypodermic Syringe**

1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml
Three-part syringes/Two part syringes (luer slip, special design to prevent plunger slipping out)
Luer slip/luer lock
Latex/latex free piston
Sterilization: EO
PE/Blister individual package, 100pcs/Inner Box
W/needle, W/O needle



FY0602 **Auto Self-destruction Syringe**

0.05ml, 0.1ml, 0.5ml, 1ml
Fixed needle 27G-30G
Sterilization: EO
PE/Blister individual package, 200pcs/Inner Box



FY0603 **Tubercle Bacillus Syringe**

1ml
luer slip
Latex/latex free piston
PE/Blister individual package
Sterilization: EO



FY0604 **Insulin Syringe**

0.3ml, 0.5ml, 1ml
Fixed needle 27G-31G
PE/Blister individual package
Sterilization: EO

FY0617 **Disposable Destruction Syringe**

1ml, 2/3ml, 5ml, 10ml, 20ml
After injection , when you pull the plunger, the needle will be back into the barrel.
It's able to prevent it from hurting people by covering hypodermic needle with barrel.
Luer lock
Sterilization: EO



FY0618 **Safety Clip Syringe**

1ml, 2/3ml, 5ml, 10ml, 20ml
Safety syringe with protective cap
Luer lock
After injection, Push the clip over the needle and lockit onto the clip base
Sterilization: EO



FY0401 **Irrigation Syringe**

50ml, 60ml, 80ml
Flat/Bulb/Ring
PE/Blister individual package
Sterilization: EO



FY0501 **Disposable Infusion Set**

For infusion bags and bottles
20/60 drops = 1±0.1ml
Small Chamber
Latex tube
With Filter
Different types of components and parts can be changed as customer's requirements
For gravity use
EO sterile
PE/Blister individual package





FY0502 Disposable Infusion Set

For infusion bags and bottles
20/60 drops = 1 ± 0.1 ml
Small Chamber
With Y site
Different types of components and parts can be changed as customer's requirements
For gravity use
EO gas sterile
PE/Blister individual package

FY0605 Disposable Hypodermic Needle

Size: 18-30G
Length of needle: 1/2", 5/8", 3/4", 1", etc.
Blister individual package
Sterilization: EO



FY0503 Disposable Infusion Set

For infusion bags and bottles
20/60 drops = 1 ± 0.1 ml
Big Chamber
With Filter
Different types of components and parts can be changed as customer's requirements
For gravity use
EO gas sterile
PE/Blister individual package

FY0606 Scalp Vein Set

Size: 18-30G
Length of needle: 1/2", 5/8", 3/4", 1", etc.
Tube length: 30cm
Luer slip/luer lock
PE/Blister individual package
Sterilization: EO



FY0504 Burette Infusion Set

Burette Capacity: 100ml, 150ml
For infusion bags and bottles
60 drops = 1 ± 0.1 ml
With Filter
With Y Site
Different types of components and parts can be changed as customer's requirements
For gravity use
EO gas sterile
PE/Blister individual package

FY0616 Insulin Pen Needle

Ultra fine cannula for maximum patient comfort, patient point with triple beveling for easy and pain decreases
The needle is fit for Humapen (EliLily, USA)
NovopenIII (Novo Nordisk, Denmark)
Sterile, non-pyrogenic
Sterilization: EO



FY0505 Disposable Blood Transfusion Set

For blood bags
20 drops = 1 ± 0.1 ml
Different types of components and parts can be changed as customer's requirements
For gravity use
EO gas sterile
PE/Blister individual package

FY0701 Disposable Blood Bag

For the collection of a blood volume: 450ml, etc.
Is filled with 63ml of CPDA-1 ANTICOAGULANT SOLUTION
Each 100ml CPDA-1 (U.S.P.) contains:
Citric Acid (anhydrous)-----0.299g
Sodium Citrate (dihydrate)-----2.63g
Monobasic Sodium Phosphate (monohydrate)-----0.222g
Dextrose (monohydrate)-----3.19g
Adenine-----0.0275g
Water for injection, a sufficient quantity to make-----100ml
Needle: 16 gauge
Steam sterilization
Sterile, non-toxic, non-pyrogenic fluid path





FY0508 Heparin Cap

Material: PC
Sterilization: EO
Individual Hard Blister Packing, 100pcs/Box, 5000pcs/ctn



FY0510 Medical Disposable Consumable Spiral Stoppers

Material: ABS Sterilization : EO
Color: Red, White, Blue, etc.
Individual Hard Blister Packing, 100pcs/Box

FY0620 Medical Connector



FY062001 Needle Free Connector

Features:
1. Transparent housing material: Polycarbonate or Copolyester.
2. Metal-free & MRI-compatible.
3. Latex-free.
4. Multiple insertion at least 100 times/day.
5. Priming Volume: 0.09ml.
6. Perfect flowing rate: 350 ml/min under one meter water pressure



FY062002 Y-port Valves

Material: Polycarbonate or Copolyester
Size: ID 4.01mm, 3.6mm



FY062003 Tube End Valves

Material: Polycarbonate or Copolyester
Size: ID 2.0mm, 2.5mm, 2.8mm, 3.7mm, 4.1mm



FY062004 Positive Pressure Needle Free Connector

Features:
1. Patent design of positive pressure avoids back flow of blood during removal of syringe. This will help to prevent blood clot at the catheter tip inside the blood vessel.
2. The housing is impregnated with PC Ag+ which helps to lower infection.
3. Protruding design of the stalk at the injection port follows infection control protocol.
4. High quality spring ensures multiple insertion at the injection site without leaking.
5. Two sealing rings at the upper and lower end of the stalk keeps the connector out air, fluid and outside substance.
6. The straight flow of fluid path creates less turbulence and complies with proper transfusion protocol.
7. The holding capacity of remaining fluid is less than 0.06ml.



FY0511 Micro-flow Regulator

Flow range: 5ml/h-250ml/h
Materials:
Regulator: ABS
Tube: PVC
Other accessories: PC & ABS
Extension Tube Length: 45cm is common, others are also available.
Sterilization: EO
Packing Details: Poly-bag per pcs, 200pcs/inner box, 20boxes/carton



Manifolds



FY0507 Three Way Stopcock

Parts: Body(PC), Plas(PE), Screw Cap(ABS)
Individual Hard Blister Packing, 50pcs/Box, 1000pcs/ctn
Sterilization: EO

FY0509 Extension Tube



FY050901
Disposable Extension Tube

The length of extension tube can be changed as customer's requirements
Different types of components and parts can be changed as customer's requirements
Sterilization: EO
Packing: Individual blister Packing

Tube Length:
10-25cm, 50pcs/ctn
30-50cm, 50pcs/ctn
80cm, 50pcs/ctn
100-120cm, 40pcs/ctn



Screw Extension Tube



Mini-Quantity Infusion
Extension Tube (OD:2.0mm)



T-Extension Tube

FY050902 Needleless Adapters for Infusion Set



- 1.This device is designed to meet the needs of general I.V. Therapy, Anaesthesia Cardiovascular,
- 2.ICU&CCU, Recovery&Oncology.
- 3.Safe for infusion IV medication including most high risk anti-neoplastics, or aspirating blood without risk of hemolysis.

- Catheter Diameter:**
- 1.One-Way Extension Tube with One Normal Needle Free Connector
 - 2.One-Way Extension Tube with One Positive Pressure Needle Free Connector
 - 3.Two-Way Extension Tube with Two Normal Needle Free Connectors
 - 4.Two-Way Extension Tube with Two Positive Pressure Needle Free Connectors
 - 5.Three-Way Extension Tube with Three Normal Needle Free Connectors
 - 6.Three-Way Extension Tube with Three Positive Pressure Needle Free Connectors
 - 7.Four-Way Extension Tube with Four Normal Needle Free Connectors
 - 8.Four-Way Extension Tube with Four Positive Pressure Needle Free Connectors



LB0630
Central Venous Catheter


Advanced automatic drilling technology and equipment, a more accurate location and the pore diameter of the pipe, surface smooth and without cutting residues and injury.

Special custom polyurethane fine materials, delicate and smooth tube, welding head end, ultra smooth surface treatment and X line to show the position, make a catheter inserted into the blood vessels in more easy and safe.



FY3007
Peripherally Inserted Central Catheter Kit

1. The connector is round and smooth with PASV valve.
- A) Valve opens inward for sampling under negative pressure.
- B) Valve opens outward for infusion under positive pressure.
- C) Valve remains closed to prevent air embolism, blood reflux and clot under normal pressure
2. The catheter is made from premium silicone, which has excellent biocompatibility.
3. Adjustable length of catheter.
4. Detectable under X ray, the end tip enhances radiopacity.
5. The style of catheter ensures placement more convenient and smooth.
6. Variety of kit configuration, initial guide wire and positive, pressure needleless connector is available for choice.



FY0631
Anesthetic Products Kit

- | | |
|-------------------------------|----------------------------|
| Standard Kits include: | Optional Compound |
| 1. Central Venous Catheter | Kits include: |
| 2. Guide-wire | 1. Central Venous Catheter |
| 3. Vessel Dilator | Standard Kit Accessories |
| 4. Clamp | 2. 5ml Syringe |
| 5. Fastener: Catheter Clamp | 3. Surgical Gloves |
| 6. Introducer Needle | 4. Surgical Pledget |
| 7. Introducer Syringe | 5. Surgery Sheet |
| 8. Injection Needle | 6. Surgery Towel |
| 9. Injection Cap | 7. Sterile Brush |
| | 8. Gauze Pad |
| | 9. Suture of Needle |
| | 10. Wound Dressing |
| | 11. Scalpel |





FY0632 Anesthetic Products Kit

Standard Kits include:

1. Drainage Catheter
2. Guide-wire
3. Vessel Dilator
4. Introducer Needle
5. Introducer Syringe
6. Injection Needle
7. Injection Cap
8. Connecting tube

Optional Compound

- Kits include:
1. Drainage Catheter
 2. 5ml Syringe
 3. Surgical Gloves
 4. Surgical Pledget
 5. Surgery Sheet
 6. Surgery Towel
 7. Sterile Brush
 8. Gauze Pad
 9. Suture of Needle
 10. Wound Dressing
 11. Scalpel



FY0635 Disposable Pressure Transducer

Disposable Pressure Transducer is made of medical PC and PVC material.

Operating Parameters:

Operating Pressure Range: -30mmHg~+300mmHg

Operating Temperature Range: 15°C ~40°C

Storage Range: -25°C ~+70°C

Atmosphere Pressure Range: 56.7kPa~113kPa

Zero Output: ± 25 mmHg

Sensitivity: $5 \mu V/V/mmHg \pm 2\%$

Input Impedance: $1200 \Omega \sim 3200 \Omega$

Output Impedance: $285 \Omega \sim 315 \Omega$

Zero Thermal Drift: $\pm 0.3mmHg/^{\circ}C$

Zero Drift: $\leq 2mmHg$

Natural Frequency: 40Hz nominal for a standard kit >200Hz for transducer alone

Sensitivity Thermal Drift: $\pm 0.1\%/^{\circ}C$

Leak Current: $< 2 \mu A$ at 230VRMS, 50Hz

Note: rated voltage 6.00VDC; ambient temperature 25 °C.

Specification:

130cm pressure tubing with M/F connector	1pc
30cm pressure tubing with M/F connector	1pc
Disposable pressure transducer	1pc
Flush device 3 cc/h (30cc/h)	1pc
Male non-vented Cap	2pcs
150cm single IV set	1pc
3-way stopcock	2pcs



FY0633 Blood Purification Products Kit

Standard Kits include:

1. Hemodialysis Catheter
2. Vessel Dilator
3. Introducer Needle
4. Syringe
5. Guide-wire
6. Adhesive Wound Dressings
7. Heparin Caps
8. Scalpel
9. Needle with Suture

Optional Compound

- Kits include:
1. Hemodialysis Catheter
 2. 5ml Syringe
 3. Surgical Gloves
 4. Surgical Pledget
 5. Surgery Sheet
 6. Surgery Towel
 7. Sterile Brush
 8. Gauze Pad



FY0634 Blood Purification Products Kit

Characteristics:

1. The Catheter is made from premium silicone, which has excellent biocompatibility.
2. The dacron cuff can protect immigration of catheter and reduce the chance of CRBI (Catheter Related Bloodstream Infections).
3. The silicone material makes catheter more floppy, hence it can minimize the trauma to the vein and reduce the incident of clot.

Standard Kits include:

1. Silicone Catheter
2. Vessel Dilator
3. Valved Peelable Introducer
4. Introducer Needle
5. Guide-wire
6. Tunneler (Tocar)
7. Adhesive Wound Dressings
8. Heparin Caps
9. Gauze Sponges
10. Scalpel
11. Syringe
12. Needle with Suture



FY0615 Epidural Set

Features:

Epidural needle: 16G/18*80/90mm

Epidural catheter: 1.0/0.8mm, clear marked, length longer than 850mm with 3 holes or more holes, closed tip

Epidural filter: 0.22um, Luer lock

LOR syringe: 7ml, 10ml available

Catheter adaptor: TBA, secure the catheter completely

Single use to avoid the cross infection

Use in the clinic of epidural nerve or subarachnoid anesthesia



FY050601
Dental Needle

Size: 27g x 38mm 30g x 38mm, etc.
Sterilization: EO
Packing: 100pcs/Box, 5000pcs/ctn



FY050602
Dental Needle

Size: 18G-25G
Packing: 100pcs/bag, 600bags/ctn



FY0611 Spinal Needle

Size: 18g~28g x 90mm
Sterilization: EO
Packing: 1pc/Bag, 50pcs/Box, 1000pcs/ctn



FY0612 A.V. Needle

Size: 17g, 18g, etc.
Sterilization: EO
Packing: 1pc/Bag, 10pcs/Box, 500pcs/ctn



FY0607



FY0608

I.V. Cannula

Normal/Safety
Pen type/With wings/With injection port
Size: 14G-26G
PTFE/FEP/PU flex catheter
Sterilization: EO
Packing: Hard Blister Packing



FY0609

All the popular sizes are available in both a carbon or stainless steel sterile
Blade shapes #9 to 17 are used in histology units and are all compatible with the associated
#3 handles (Blade shape 11 can be also used for gardening with a plastic or metal handle).
Blade shapes #18 to 36, #60 and #70 are used for major invasive procedures and are all
compatible with the associated #4 handles
Sterilization: Gamma
Packing: 1pcs/Sealed Aluminum Foil Bag,
100pcs/box, 5000pcs/ctn



FY0901
Disposable Surgical Blade

FY2209 Disposable electrosurgical pencil (Tips)



FY2209 Disposable electrosurgical pencil

FY220901



Electrode Shape: Blade
Handle Length: 158mm
Hand control
Cable Length: 3m
Exposed electrode length: 40mm

FY220902



Electrode Shape: Loop
Handle Length: 158mm
Hand control
Cable Length: 3m
Electrode Exposed Length: 40mm

FY220903



Electrode Shape: Blade
Handle Length: 158mm
Foot control
Cable Length: 3m
Exposed electrode length: 40mm



FY0902 Stitch Cutters

All sizes are available in both a carbon or stainless steel sterile
Designed with a rounded bull nose tip, the Index family of Stitch Cutter Blades are ideal for cutting the wide range of modern sutures available today whilst minimising patient discomfort
The range consists of a short standard stitch cutter with a larger grip. The long handle enables a more aseptic technique allowing good visibility whilst working further from the wound site
Size: 65mm, 110mm
Sterilization: Gamma
Packing: 1pcs/Sealed Aluminum Foil Bag, 100pcs/box, 5000pcs/ctn



**FY0907 Twist Lancets
(Exclusive to Roche Diagnostics)**

5 twist lancets (28G) are used with Roche Accu-Chek Softclix Lancing Device and the needle tip is siliconized for easy and smooth motion during lancing
Size: 28G, 30G
Packing: 100pcs/Box, 20000pcs/ctn



FY0903 Surgical Handle

Index stainless steel metal handles are supplied non-sterile and individually wrapped
Sterilization: Gamma
Packing: 1pcs/Sealed Aluminum Foil Bag, 10pcs/box, 500pcs/ctn



FY0910 Safety Lancet

Ergonomic design offers a more secure hold when to activate the device
Tri-bevel needle tip ensures the sharpness and accuracy
The spring buffer the puncture making the procedure a more comfortable experience
The complete range of needle sizes from 18G to 32G ensure the accurate
Blood samples for all popular meters
The needle retracts immediately after use
Sterilization : Gamma
Size: 21G, 23G, 26G



FY0904 Surgical Scalpels

Standard sheath guarded and blade protected sterile disposable scalpels offer convenience in a number of user environments: on the ward, in emergency or trauma units or even out in the field with paramedics
Material: Carbon Steel or Stainless Steel
Sterilization: Gamma
Packing: 10pcs/Box, 500pcs/ctn



FY0911 Pressure Activated Safety Lancets

Needle is fully shielded before and after use
Pre-loaded, pressure activated design
Ensures a good blood sample every time
Ultra sharp tri-bevel needle or bi-bevel blade, high speed penetration
Variety of needle specifications to suit most capillary blood requirements
Sterilization : Gamma
Size: 21G, 23G, 26G



FY0905 Steel Lancet

Low cost: simple, flat, steel lancet
Sample volume: excellent blood flow
High quality: hardened surgical steel
Minimum discomfort: modern asymmetric design
Sterile: sterilized by gamma radiation
Size: 40 × 5.5 × 1.5mm
Packing: 5pcs of packaged, 200pcs/Box, 20000pcs/ctn



FY0912 Button Activated Safety Lancets

Needle tip is sealed with cap through injection moulding
Pre-loaded, button operated design
Single-use only, sterilised by gamma irradiation, cannot be reused
High quality, ultra sharp tri-bevel needle, high speed penetration
Variety of needle specifications to suit most capillary blood requirements
Sterilization : Gamma
Size: 21G, 23G, 26G



FY0906 Twist Lancets

Conventional twist-off cap design
Sterilized by gamma-radiation
Smooth tri-bevel needle tip for comfortable sampling experience
Versatile with most lancing devices
Size: 21G, 23G, 26G, 28G, 30G
Packing: 100pcs/Box, 20000pcs/ctn



**FY0908 Lancing Device (Releasable)
FY0909 Lancing Device**

The user can install the lancet with the device easily. And the device will work with the lancet smoothly without any blockage. The lancet won't fall on its own when to fire the device. Retracting automatically after firing, the launcher is safe and reliable.





FY0913 / FY0914 / FY0915 Medical Razor
Carbon steel or stainless steel Single or double blade All razors come with the highest quality stainless steel blades or carbon steel , Comfort-coated for smooth shaves



Surgical Suture


Cat. NO.	Type
FY1001	Chromic Catgut
FY1002	Plain Catgut
FY1003	Silk Braided
FY1004	Nylon Monofilament
FY1005	Polyester Braided
FY1006	Polypropylene
FY1007	Polyglycolic Acid

Thread Length: 45cm, 60cm, 75cm, 90cm, 150cm etc.
Thread Diameter: 8/0,7/0,6/0,5/0,4/0,3/0,2/0,1/0,1,2,3 etc.
Needle Length: 6mm,8mm,12mm,18mm,35mm,40mm,50mm etc.
Needle Curvature: Straight,1/2circle,1/2circle(double),1/4circle,1/4circle(double),3/8circle,3/8circle(double),5/8circle,loop round
Cross Section: round bodied,round bodied(heavy),curved cutting, curved cutting(heavy),reverse cutting,reverse cutting(heavy) etc.



FY1502 Splints

Aluminum plate Macromolecule foam material Nontoxic and nonbenzenoid
Free of effect of temperature and climate Type available: Fold, Unfold and Roll
Normal size: 5cm×12cm, 5cm×15cm, 7cm×11cm
11cm×44cm, 11cm×91.5cm, 11cm×100cm



FY1501 Thermal Accident Blanket

Material: PET Film
Waterproof, use to offset hypothermic reactions, reflective surface provides high visibility, reflects sun and heat for cooling in hot climates, space-age design retains up to 90% of body heat,Assists in prevention of post-accident shock and first aid application
Size: 132×210cm, 140×220cm, 160×210cm or other size per customer demand
Color: Gold, Sliver, Gold+Sliver



EXAMINATION





FY0802 Surgical Gloves

Size: 6.0#, 6.5#, 7.0#, 7.5#, 8.0#, 8.5#, 9.0#
Powdered/Powder free
1pair/pouch, gamma sterile
AQL: 1.5
Packing: 25pairs/Box, 500 Pairs/ctn



FY0803 Latex Examination Gloves

Size: S,M,L
Powdered/Powder free
100pcs/box, non-sterile
AQL: 2.5
Packing: 100pcs/Box, 1000 Pcs/ctn



FY0804 Nitrile Examination Gloves

Size: S,M,L
Powder free
100pcs/box, non-sterile
AQL: 1.5/2.5
Packing: 100pcs/Box, 1000 Pcs/ctn



FY0805 Vinyl Gloves

Size: S (3.0g,3.5g,4.0g) ; M (3.5g, 4.0g, 4.5g) ;
L (4.0g,4.5g,5.0g)
Powdered/Powder free
100pcs/box, non-sterile
AQL: 4.0
Packing: 100pcs/Box, 1000 Pcs/ctn
Corlor: Blue, White



FY0806 PE Gloves

Size: S,M,L
Thickness: 0.018mm, 0.02mm, 0.022mm etc.
Type: HDPE/LDPE
Packing: 100pcs/Box, 10000 Pcs/ctn



FY0807 Two Fingers Cots

Size: S,M,L
Powdered/Powder free
Packing: 100pcs/Inner Bag, 5000 Pcs/ctn

FY2201 Vaginal Speculum



FY220101 Side Screw Type

Size : S,M,L



FY220102 Fastener Type

Size : S,M,L



FY220103 Spanish Type

Size : S,M,L



FY220104 Middle Screw Type

Size : S,M,L



FY220106 Push Pull Type

Size : S,M,L



FY220107 American Type

Size : S,M,L






FY220201 Cervical Brush




FY220202 Cervical Brush



FY220203 Cervical Brush



FY220204 Cervical Brush



FY220205 Cervical Scraper



FY220206 Cervical Scraper



FY220207 Cervical Scraper



FY220208 Cervical Scraper



**FY2101
Wooden Tongue Depressor**

Size: Adult (150mm × 18mm × 1.6mm)
Child (140mm × 14mm × 1.6mm)
EO Sterile / Non-Sterile
Packing: Single medical paper wrapped/Bulk
100pcs/Box, 5000pcs/Ctn




**FY2102
Cervical Scraper (Wooden)**

EO Sterile / Non-Sterile
Packing : Single medical paper wrapped/Bulk
100pcs/Box, 5000pcs/Ctn



**FY2205
Gynecological Kit**

Configuration :
1pc speculum
1pc Gynaecology pad(40*60cm)
1pair latex glove
1pc cervical spatula



**FY1904 Condom
Specification:**

Surface: smooth, dotted, ribbed lubricated
Size: L, M, S (55mm±2、52mm±2、49mm±2).
Be available with spermicide nonoxynol 9
Be available with various fruit smell,such as apple/ banana/ strawberry/ cherry etc.
Be available with square or rectangle packing ;
If quantity is large enough, the condom is available with various colorsuch as natural/ red/pink/green etc.
Specification is available with En600,ISO4074, ASTM, WHO, SABS, etc.



**FY1901
HCG Test Strip**

Size: 2.5mm, 3.0mm, etc
Response time: 5 minutes
Sensitivity: 25mIU/ml
Specificity: 99%
Packing: 1pc/aluminum foil bag



**FY1902
HCG Test Devise**

Size : 3.0mm , 3.5mm , etc
Response time: 5 minutes
Sensitivity : 25mIU/ml
Specificity : 99%
Packing : 1pc/aluminium foil bag



**FY1903
HCG Midstream**

Size : 4.0mm , 4.5mm , 5.0mm , etc
Response time: 5 minutes
Sensitivity : 25mIU/ml
Specificity : 99%
Packing : 1pc/aluminium foil bag

WOUND DRESSING

FY1101 Gauze Swabs

The pads are folded all by machine, pure 100% cotton yarn ensure the product soft and adherent. Superior absorbency makes the pads perfect for absorbing blood and exudates. In accordance with the customers' requirements we can produce different kinds of pads, such as folded, and folded, with x-ray and non x-ray. The adherent pads are perfect for operation. General specification: 2"×2", 3"×3", 4"×4" PLY: 8ply, 12ply, 16ply Mesh: 19×10, 19×15, 20×12, 26×18 EO Sterile / Non-Sterile



FY1105 Non-woven Swabs

Non-woven swabs from 70% polyethylene and 30% terylene. The pad is designed to wick away fluids and disperse them evenly. Product have been cut like "O" and "Y", so it is easy to use. It is mainly used to absorb blood and exudates during operation and cleaning the wounds. General specification: 30g/m², 40g/m²; 2"×2", 3"×3", 4"×4" PLY: 4ply, 8ply EO Sterile / Non-Sterile




FY1103 Prewashed Laparotomy Sponge

Laparotomy Sponge for absorbing blood and exudates or for cleaning wound. Cotton Yarn: 40's x 40's It can prevent the wound from being infected. Our laparotomy sponges are sewn and cross-stitched. General specification: 45×45cm/30×30cm, PLY: 8ply, 16ply Mesh: 19×15/26×18 EO Sterile / Non-Sterile Prewashed/ Non-Prewashed



FY1102 Gauze Roll

Material : 100% cotton Cotton Yarn: 40's×40's Weaving according to customer requirements, high pressure, high temperature degreasing, bleaching, refined but become, with x-ray or without x-ray Products without fluorescence, soft, has the very strong absorption ability, is the first choice for medical units use products, quality in line with the British pharmacopoeia standards, are widely used in other markets General specification: 36"×100yard-4ply, 36"×100m-4ply, 36"×1000m-4ply; Mesh: 19×11/19×15





FY1104 Gauze Bandage

Material : 100% cotton
Cotton Yarn: 40's×40's
Weaving according to customer requirements,
high pressure, high temperature degreasing,
bleaching, refined but become
Packing : 1roll/white paper, 1doz
blue kraft paper packing
General specification:
2" ,3" ,4" ,6" ×4.5m/5m/10m;
Mesh:19×15,20×12,26×18
EO Sterile / Non-Sterile

FY1108 Cotton Ball

Soft & Pure balls are also highly absorbent, offering
outstanding cleansing properties as they absorb
fluids extremely quickly, making them ideal for use
with toner and nail varnish or make up remover,
without drips or mess.
Made from 100% pure cotton
Gentle to the skin
General specification: 0.2g, 0.3g, 0.5g, 1.0g
EO Sterile / Non-Sterile



FY1106 Cotton Roll

Absorbent cotton is pure 100% cotton fabric,
degreased and bleached by advanced way.
The cotton is perfect for beauty; light drainae
wounds, the alcohol sterile, cleaning the medical
instruments.
Packing become rolls by white and blue paper
or PE bag.
General specification:50g,100g,200g,250g,454g,
500g,1000g
EO Sterile / Non-Sterile

FY1110 Cotton Applicator (Cotton Swab)

Cotton Tip Applicators are available in sterile and
non sterile, different lengths and packet sizes, with
a choice of wooden and plastic stems; sterile cotton
tipped applicator or non-sterile cotton tipped applicators.
Various size can according to the specifications of
an individual purchaser, use for any skin care,
peeling, facial, nail, eyebrows applications.
Size: 3", 4", 6", 8" with wooden, plastic, using 100%
absorbent cotton
EO Sterile / Non-Sterile



FY1109 Gauze Ball

100% pure cotton yarn of our product ensures
a superior absorbency.
There are two types of this product: with x-ray
and non x-ray ball.
It is mainly used to absorbing blood and exudates
during operation and claeing the wounds.
General specification: Dia.:3cm, 4cm, 5cm;
mesh:19×15, 26×18
EO Sterile / Non-Sterile

FY1401 Dental Cotton Roll

The dental cotton roll is made of 100% cotton,
it is widely used on dentistry and dental cleaning.
Size: 0.4g, 0.5g, 0.7g
Size: 8 × 38mm, 10 × 38mm, 12 × 38mm
EO Sterile / Non-Sterile




FY1107 Zig-zag Cotton

Absorbent cotton is pure 100% cotton fabric,
degreased and bleached by advanced way.
The cotton is perfect for beauty; light drainae
wounds, the alcohol sterile, cleaning the medical
Instruments.
Packing: sealed poly bag, drawstring bag
or self-sealing bag
General specification: 25g,50g,60g,100g,250g,500g
EO Sterile / Non-Sterile

FY1406 Oral Cleaning Swab

Material: PP and sponge
Color: green, red, blue, white or others;
Product Specification: 10cm, 15cm, 18cm,
or as your request.
Gently cleans, moistens and refreshes mouth,
teeth, gums and disinfect the wounds.
Raised foam ridges remove particles and
mucous for enhanced oral hygiene





FY1402 Salive Ejector

Material: PVC
Length: 15cm
Packing: 100pcs/bag,
20bags/ctn,
35×29×24cm



FY1419 3 way syringe tip

3 Function Sparayer Covers
185mm×127mm ×170mm,
500pcs/box, 18000pcs/ctn



FY1404 Dental Bib

Item: Disposable Bibs
Material: 2×20g paper Laminated 20g PE film
Color: White/green/pink/blue etc.
Packing: 125pcs/bag, 4bags/ctn 35×25×25cm



FY1408 Disposable Plastic Forceps

Material: plastic
Colour: red, blue, green, yellow,
white or as customers' request.
Description: Our first and tweezer
are widely used in hospitals, clinics
and also in first aid.
Dimension: 10.8cm / pc, 11.5cm/pc,12.5cm / pc ,
13 cm / pc, 14cm/pc,18 cm/pc etc.
Note: Logo and print can be finished
according to customers' request.



FY1409 Haemostatic Forceps

Material: plastic/ABS
Colour: red, blue, green, yellow,
white or as customers' request.
Dimension: 14×7.5cm

FY1201 Spandex Elastic Bandage

Material : cotton and spandex
Color : natural,flesh or bleach
Weight : 60g/m²,70g/m²,75g/m², etc.
Size:5cm×4.5m,7.5cm×4.5m,
10cm×4.5m,15cm×4.5m



FY1202 PBT Elastic Bandage

Material : cotton and polyester
Color : natural or bleach
Weight : 25g/m²,28g/m²,30g/m²
Size:5cm×4.5m,7.5cm×4.5m,
10cm×4.5m,15cm×4.5m



FY1203 Cotton Elastic Bandage

Material : 100% cotton
Color : natural or bleach
Weight : 60g/m²,70g/m²,75g/m², etc.
Size:5cm×4.5m,7.5cm×4.5m,10cm×4.5m,
15cm×4.5m



FY1204 High Elastic Bandage

Material : rubber and polyester
Color : white,flesh or bleach
Weight : 80g/m²,100g/m², etc.
Size:5cm×4.5m,7.5cm×4.5m,
10cm×4.5m,15cm×4.5m





FY1205 Elastic Net Bandage

Material: rubber and polymide
Color: white
Weight: 120g/m², etc.
Size: 0.8cm×25m, 1.7cm×25m, 2cm×25m,
2.2cm×25m, 2.6cm×25m,
3.8cm×25m, 4.4cm×25m, .5cm×25m,
6.8cm×25m, 8cm×25m



FY1206 Tubular Bandage

Material: rubber and cotton
Color: natural
Size: 4.5cm×10m, 6.25cm×10m,
5cm×10m, 6.75cm×10m, 7.5cm×10m, 8.75cm×10m,
10cm×10m, 12cm×10m, 17.5cm×10m,
21.5cm×10m, 32.5cm×10m



FY1220 Stockinette Bandage

Material: 100% Cotton
Color: Natural white, or other color color customer's demand
Introduction: Made from 100% natural cotton, which is used
for the wound and sport. It can be made into roll
Size: 5.0cm×20m, 6.0cm×20m, 8.0cm×20m,
10cm×20m, 15cm×20m



FY1207 Non-woven Self-adhesive Bandage

Material: non-woven
Weight: 35g/m², 40g/m², 45g/m², etc.
Color: 12 colors
Size: 5cm×4.5m, 7.5cm×4.5m,
10cm×4.5m, 15cm×4.5m



FY1208 Cotton Self-adhesive Bandage

Material: cotton
Weight: 75g/m², 80g/m², 85g/m², etc.
Color: 12 colors
Size: 5cm×4.5m, 7.5cm×4.5m,
10cm×4.5m, 15cm×4.5m

FY1209 Plaster Of Paris Bandage

Material: cotton and plaster
Weight: 390g/m², 410g/m², etc.
Color: white
Size: 5cm×2.7m, 7.5cm×2.7m,
10cm×2.7m, 15cm×2.7m



FY1218 Cast Padding

Material: natural cotton , polyester , etc
Color: white
Size: 5cm×2.7m, 7.5cm×2.7m,
10cm×2.7m, 15cm×2.7m



FY1210 Spandex Self-adhesive Bandage

Material: spandex
Weight: 75g/m², 80g/m², 85g/m², etc.
Color: skin
Size: 5cm×4.5m, 7.5cm×4.5m,
10cm×4.5m, 15cm×4.5m



FY1216 Cool Bandage

Cool Bandage is a cold bandage which uses a
revolutionary new technique for treating
injuries, muscle pain and strains. This bandage
reduces recovery time from injuries by combining
both the cooling and compressing process in one
easy to use step.
Color: brown, blue
Size: 7.5cm×3.2m, 10cm×3.2m





FY1217 Fiberglass Casting Bandage

Description: Fiber Glass casting tape consists of a lightweight knitted polyester fabric, impregnated with a water absorbent polyurethane resin.
Setting time: 3-5 minutes
Weight-bearing: 20 minutes
Excellent conformability & mouldability
Tack free
100% transparent to X-rays
Very good water vapour permeability
With elasticity, and flexibility
Blue, Green, White, Pink, Red color
Size: 5cm×3.6m, 7.5cm×3.6m, 10cm×3.6m

FY1214 Silk Surgical Tape

The silk surgical tape mainly uses low-sensitization medical hot-melt adhesive, the back face is coated with an uncoiling agent, the back lining material is mainly made of soft, light and thin silk cloth serving as a base material, and the adhesive is coated in a shape of mesh.
Packing: simple package, plastic can, with cutter
Color: white
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.



FY1211 Zinc Oxide Plaster

Made by using cotton cloth as a base material and coating adhesive cement mainly made of zinc oxide and rubber or medical hot-melt adhesive.
Packing: simple package, plastic can, galvanized can
The product is high in stickiness, and can be cut into a saw-tooth shape
Color: white, flesh
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.

FY1215 Sport Tape

The sports tape uses cotton cloth serving as a base material and adhesive cement mainly made of zinc oxide and rubber, and the adhesive cement is coated in a shape of strip.
Different colors available
Width: 2.5cm, 3.8cm, 5cm, etc.
Length: 10yard, 15yard, etc.



FY1212 Non-woven Surgical Tape

The nonwoven surgical tape mainly uses low-sensitization medical hot-melt adhesive, the back face is coated with an uncoiling agent, and the back lining material is mainly made of soft, light, thin and air-permeated nonwoven cloth serving as a base material.
Packing: simple package, plastic can, with cutter
Color: white, flesh
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.

FY1301 Drilled Plaster

The product is non-toxic, low in sensitization, not irritant to the skin and subjected to EO sterilization; the product is a kind of externally applied dressing and has the effects of stopping bleeding and protecting wounds, and it is applied to affected parts such as small wounds and abrasions; and the product can prevent the intrusion of bacteria and foreign bodies and accelerate the wound healing.
Specification: PVC (natural/colored), PU, PE, elastic cloth, non-woven cloth, cartoon, etc.
Size: φ22mm, φ25mm, 38cm×38cm, 72cm×19cm, 76cm×19cm, 76cm×38cm, etc.



FY1213 PE Surgical Tape

The PE surgical tape mainly uses low-sensitization medical hot-melt adhesive, the back face is coated with an uncoiling agent, and the back lining material is mainly made of soft PE films serving as a base material.
Packing: simple package, plastic can, with cutter
Color: transparent
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.

FY1302 Healing Capsicum Plaster

The product is mainly formed by two parts: a base material (cotton material) and separate paper. The product has the features of air permeation, moisture permeation, tight fixing, high adaptability, convenient use, etc.





FY1303 First Aid Plaster

The product is low in sensitization, excellent in biological compatibility, bacteria-free and high in air permeation, and doesn't irritate or damage the skin nor affect the natural functions of the skin. Specification: PU, spun-laced non-woven cloth, etc. Size: 8cm×6cm, 8cm×10cm, 10cm×10cm, 10cm×15cm, 10cm×20cm, 10cm×25cm, etc.



FY1304 Eye Pad

There are two kind of Eye pad. One is non adhesive 100% Cotton eye pad. The other is made of an elastic non-woven dressing with absorbent pad. Both dressing are Specially indicated for strabismic, weak sighted therapy and nursing after ocular operation. There are different sizes for adults and children available.



FY1305 Non-woven Adhesive Pad

The external layer of the dressing pad is formed by spun-laced nonwoven cloth, the internal layer is formed by needle-punched nonwoven cloth, and the surface is formed by release paper; the base material is soft, and the product should be sterilized and is mainly used for wound dressing in the department of surgery of a hospital.



FY1310 Vaseline Gauze Swab(Paraffin Gauze)

The product is made from medical degreased gauze or non-woven together with Vaseline. It can lubricate the skin and protect the skin from cracks. It is widely used on clinic.



FY130601

FY130602
2% Chlorhexidine

FY130601 / FY130602 Alcohol Prep Pad (Alcohol Swab)

Contents: 70% Isopropyl Alcohol
Non-Woven Size : 65×30mm, 65×58mm, etc
Type: 2ply, 4ply, etc.
Packet Size: 50×50mm



FY1307 Povidone-Iodine Prep Pad

Contents: 10 % Povidone-Iodine Solution
Non-Woven Size : 65×30mm, 65×58mm, etc
Type: 2ply, 4ply, etc.
Packet Size: 50×50mm



FY1309 Povidone-Iodine Swabstick

Contents: 10 % Povidone-Iodine Solution
Swabstick Size: 70mm, 100mm
Packet Size : 110×40mm, 140×40 mm



FY1315 Povidone-Iodine

Cleaning Scrub Swabstick 1pc
Antiseptic Prep Swabstick 1pc



FY1308 Sting Relief Medicated Pad

Contents: Antipruritic solution
Non-Woven Size: 65×30mm, 65×58mm, etc
Type: 2ply, 4ply, etc.
Packet Size: 50×50mm



FY1311 Alcohol Prep Pad with Chlorhexidine

Contents: 70% Isopropyl Alcohol + 2% Chlorhexidine
Non-woven 50g/m²
Size: 65mm×30mm-2ply
Sachet Size: 50mm×50mm
Package: 1pc/Airtight Sachet, 100pcs/box



FY1312 Wet Wipes

Specifications:

Material: Spunlace Nonwoven, ES hot air through Nonwoven, Thermal-bonded Nonwoven, Air-laid paper, Wet Strength Paper, Bamboo Fiber.

Weight: 25-80gsm for your choice

Size: 15x20cm, 16x20cm, 18x20cm, 20 x 30cm, etc.

Wipes count: 1 to 200pcs for your choice.

Folding Size: Z fold, M fold.

Formula: Gentle and safe, Alcohol Free, all kinds of fragrance for your choice.

Fragrance: Lemon, Lavender, Camellia, etc.

Function extracts for your choice, such as Aloe Vera & Vitamin E and so on.

Packaging type: PE/PET film, plastic canister/tube, or plastic box, etc.



FY1313 Disinfectant Wipes

Specifications:

Material: Spunlace Nonwoven, ES hot air through Nonwoven, Thermal-bonded Nonwoven, Air-laid paper, Wet Strength Paper, Bamboo Fiber.

Weight: 25-80gsm for your choice

Size: 15x20cm, 16x20cm, 18x20cm, 20 x 30cm, etc.

Wipes count: 1 to 200pcs for your choice.

Folding Size: Z fold, M fold.

Fragrance: Lemon, Lavender, Camellia, etc.

Function extracts for your choice, such as Aloe Vera & Vitamin E and so on.

Packaging type: PE/PET film, plastic canister/tube, or plastic box, etc.

**MEDICAL
NON-WOVEN
& ACCESSORIES**



FY1601 Round Nurse Cap

Material: non-woven PP
Weight: 10g/m², 12g/m², 14g/m², 16g/m², etc
Size: 18", 21", etc
Color: white, blue, green, etc



FY1603 Doctor Cap

Material: non-woven PP
Weight: 20g/m², 25g/m², 30g/m², etc
Size: 61×14cm, 64×23cm, etc
Type: tie/elastic back
Color: white, blue, green, etc



FY1602 Clip Cap

Material: non-woven PP
Weight: 10g/m², 12g/m², 14g/m², etc
Size: 20", 21", 22", etc
Type: single elastic/double elastic
Color: white, blue, green, etc



FY1604 Net Cap

Made from high quality nylon net fabric material
Weight: 2.7g/m², etc
Size: 19", etc
Color: white, blue, green, etc

FY1609 Surgical Hood

Material: non-woven PP
Weight: 16g/m², 20g/m², 25g/m², 30g/m², etc
Size: 30×124cm, etc
Type: tie
Color: white, blue, green, etc



FY1605 Non-woven Face Mask

Material: non-woven PP or non-woven PP and filter paper
Specification: ear-loop/tie and 2ply (18g/m²+25g/m²)
/3ply (18g/m²+20g/m²+25g/m²) /4ply(22g/m²+25g/m²+23g/m²+25g/m²)
Size: 17.5×9.5cm
Color: white, blue, green, etc



FY1632 Non-woven Face Mask with Charcoal

Description:
Conforms to American NEL SON N95 test standards.
Adjustable nose seal, easy to fold and carry. Mini-fiber high electrostatic filter, filtration efficiency of dust & powder is over 95%.
Repel continually contamination of dangerous blood and saliva continually.
Specification: Ear-loop type: 4 ply
Detailed weight: 4 ply: 25g/m²+20g/m²+40g/m²+25g/m²
Size: 18×9cm
Color: white, green, blue etc.
Packing: 50pcs/box, 40boxes/carton, Dia.: 520×380×300cm



FY1606 Paper Face Mask

Material: high quality wet-strength paper
Specification: ear-loop and 1ply/2ply
Size: 20 × 7cm
Color: white



FY1608 Face Mask With Shield

Composition : non-woven face mask and PET shield (prevent mist)
Specification : ear-loop and 3ply/4ply
Size : 17.5×9.5cm
Color : white, blue, green, etc



FY1629 Face Guard

Material : PET
Type : elastic double-sided anti-fog
Pack : 1pc/pe bag



FY1301 Protective Spectacle

Impact resistant PC lens, Lens thickness: 2mm
Integrated side legs
High strength structure
Integrated wide side shields





FY160701 Simple-Dust-Proof FaceMask

Conforms to FFP1& FFP2
Fiber Glass Free



FY160702 New-Type-Dust-Proof FaceMask

Material: terylene (conforms to American
NEL SON N95 test standards)
Specification: ear-loop
Size: 220g,240g,260g,etc
Color: white,blue,etc



FY160703 N95 Facemask (with/without valve)

Material: polypropylene (conforms to American
NEL SON N95/FFP1/FFP2/FFP3 test standards)
Specification: ear-loop
Size: without valve or with valve
Color: white



FY160704 Folded Facemask

Material: polypropylene (conforms to American
NEL SON N95/FFP1/FFP2/FFP3 test standards)
Specification: ear-loop
Size: without valve or with valve
Color: white



FY1610 Surgical Gown

Standard and Reinforced
Material: Spunlace, SMS, SU/US nonwoven and so on
Gram weight of material: from 35g/m² to 50g/m²
Color can be decided by customer
Whether sterilize is be decided by customer
Packaging way can be bulk or individual
Size: S, M, L, XL, XXL, XXXL



FY1611 Disposable Isolation Gown

Material: SPP or PP+PE non-woven fabrics
or as customers' requirement
Color: Blue, yellow, green and so on.
Size: S/M/L/XL/XXL
Details: knitted or elastic cuff
with velcro or ties on the neck and ties on the waist.



FY1613 Lab Coat

Material: non-woven PP/PE
Color: white/blue/yellow/green
Weight: 30g, 35g, 40g, 45g, etc
Size: M/L/XL/XXL

FY1612 CPE Coat

Material: CPE
Color: Blue
Size: 62g, 125×200cm, etc
Characteristics:
prevent and isolate dust,
particle, alcohol, blood
Bacteria and virus invading



FY1610R Reinforced Surgeon Gown

Specification:
Material: Spunbond Meltdown Spunbond (SMS)
Color: whiteue/yellow/green
Weight: 40gsm/45gsm/48gsm/50gsm/55gsm etc.
Size: S/M/L/XL/XXL
Suited to hospital use & for surgery
With ties at back of the neck and waist
Comfortable sleeve tailoring for maximum
Freedom of movement
Soft knit cuffs for wrist comfort
Secure sleeve seams for protection
Extra protection reinforce on sleeve & chest
And full-body on selected gowns
Spunbond meltdown spunbond(SMS) material gown breathable,
quiet and supple,
High repellency, AAMI standard, sterile with comfortable design
for operating theater
With knitted cuff, Velcro neck back closure, side waist band
Reinforced arm and additional liner at front
1 hand towel is available
1 wrap towel is available
Sterile method: non-sterile or EO sterile
Packing: 10pcs/poly bag or 1 pc/sterile bag



FY1615 Disposable Protective Coverall

Chemical protective coverall includes: liquid-tight protective clothing, spray-tight protective clothing, protective clothing for use against solid particulates, protective clothing against liquid chemicals and so on.
Material: PP/PP+PE/SMS and so on
Weight: 35-60g
Size: M/L/XL/XXL
Color: Blue/white/yellow/ dark blue
Application: provide reliable protective products for industry/lab/ nuclear radiation/ agriculture/ military combat readiness and so on



FY1614 Disposable Scrub Suit

Size: S/M/L/XL/XXL and so on
Material: SS/PP/SMS 20g - 60g
Color: Blue/white/yellow/ dark blue
Packing: 1pc/PE bag, 50bags/carton



FY1618 CPE Shoe Cover

Material: CPE, with comfortable elastic opening
Color: blue, green, etc
Weight: 1.6g, 1.8g, 2.0g, 2.1g, 2.5g
Size: 36×15cm, 40×16cm
Thickness: 0.018mm, 0.02mm, 0.021mm and so on



FY1617 PE Shoe Cover

Material: PE, with comfortable elastic opening
Color: blue, green, etc
Weight: 1.6g, 1.8g, 2.0g, 2.1g, 2.5g
Size: 36×15cm, 40×16cm
Thickness: 0.018mm, 0.02mm, 0.021mm and so on



FY1619 Non-woven Shoe Cover

Material: non-woven PP
Color: white, blue
Regular/Non-skid
Weight: 30g/m², 40g/m², 45g, etc
Size: 36×15cm, 38×18cm, 40×16cm



FY1616 PE Sleeve

Made of PE, with comfortable elastic
Color: white, blue and red, etc
Thickness: 0.018mm
Size: 20×40cm, 22×42cm, 22×46cm



FY1622 PE Apron

Made of polyethylene, with adjustable tie
Color: white, blue, etc
Size: 110×70cm, 135×71cm, 160×90cm
Thickness: 0.018mm, 0.02mm and so on


FY1623 Non-woven Bed Sheet

Soft touch and harmless to your skin
100%biodegradable and environmentally friendly
SPEC: light blue, 60×180cm
30gsm other size and specification are also acceptable



FY1624 Non-woven Bed Cover

Disposable bed cover made of pp non-woven fabric helps you to be clean & hygiene on your travel and regular purpose
SPEC: 30g, 210×90×20cm, other size and specification are also acceptable



FY1625 Nonwoven Pillow cover

Disposable pillow covers are used widely in hotels, hospitals, home, and airline and in all travels

Specification:
80gsm, 50×60cm, with 5cm gusset,
other size and specification are also acceptable
Packing: 100pcs/polybag, 1bag/ctn
Ctn size: 54×32×38 /100pcs





FY1620 Under Pad

Material: non-woven fabrics+2ply medical paper
+ fluff pulp+PE film, non-woven fabrics+4ply
medical paper+PE film, etc
with or without SAP
Size: 40×60cm, 60×90cm, etc
Weight can be decided by customer



FY1621 Adult Diaper

First layer: Non-woven fiber, natural soft
and breathable.
Second layer: fuff pulp+super absorbent polymer
+classic tissue.
Breathable PE film with cloth like back Non-woven
or casting film.
Standard Size:M: 650×800mm
L: 750×900mm, XL:840×1030mm



FY140702 Surgical Brush

Medica Brush has an effective, soft and skin-friendly brush
head for the nails, knuckles and skin.
Attached to the brush is a sponge designed for cleaning the
hands and forearms.
In addition a separate nail cleaner is supplied in the packaging.
The Medica Brush saves on handling costs since brush, soap
and nail cleaner are all supplied in one packaging.
In addition to a dry sterile version, one has a choice of several
impregnated versions:
Medical Brush Dry Sterile
Medical Brush Chlorhexidine
Medical Brush PVP Iodine



FY140701 Surgical Hand Brush

Surgical scrub brush is used to significantly reduce the number
of micro-organisms
on the hands and forearms before performing surgeries.
The nail cleaner feature makes removal of micro-organisms
under the nail quick and easy
Dry Surgical Brush includes a variety of antimicrobial formulations,
including CHG, PCMX,
and Povidone Iodine
All kinds of spray cans are available



FY150301 Sterile Delivery Pack

Non-woven surgical gown-1pc
86×80cm packing clothes 2ply(Non-woven and thin film)-1pc
97×73cm sheets 3ply(Thin film, Non-woven and paper pulp)-1pc
97×73cm Triangle trouser 3ply(Thin film, Non-woven and paper pulp)-2pcs
66×43cm Non-woven cloth 2ply(Thin film and Non-woven)-2pcs
66×43cm Non-woven cloth 1ply-1pc; Umbilical cord clamp-1pc
8cm cotton sticker-2pcs; PE glove-1pairs; Latex examination glove-1 paris
8×10cm Gauze swab-5pcs; Face mask-1pc; Cap-1pc



FY150302 Sterile Delivery Pack

Instrument table cover 150×200CM	1PC
Side drape 100×130CM	1PC
Simple drape 100×130CM	1PC
Hand towel 40×40CM	1PC
Leggings drape 75×120CM	2PCS
Clip	1PC
Reinforced gown L 45gsm SMS	1PC
Gauze 7.5×7.5cm	10PCS
Wrapping cloth 100×100CM	1PC
Side drape 75×90CM	1PC



FY150303 Sterile Delivery Pack

1Pc Instrument Table Cover Size 150 x 200
1Pc Caesarean Drape With Pouch Size 200 x 300,60gsm SPP Bi-lamination.
1Pc Baby Blanket size 75 x 90
3Pc Hand Towels size 40 x 40
1Pc Op Tapes size 10 x 50
3Pc Reinforced Surgical Gown Size L, Individual packing, EO Sterile.



FY4201 Spinal Epidural CVC Sterile Pack

Fenestrated pe Drape 610 × 700mm one,
Luer Lock Syringe Two,3ml and 5ml
Applicator Sponge Three,
Instrument table cover one,120×100cm
Hand Towel two,30×40cm
Green Gauze 100×100mm eight
Filter Needle 19g,two
Medical Plastic Tray one,24×25×4.5cm
1pack/blister packing,30packs/ctn

OTHERS



FY1701
Hot Water Bag

Material: calcium carbonate \geq 30% rubber
Size: 500cc, 1000cc, 1750cc, 2000cc, etc
Packing: 1pc/PE bag



FY1706
Ear Syringe

Material: calcium carbonate \geq 25% rubber or PVC
Size: 30ml, 60ml, 90ml, etc
Packing: 1pc/PE bag



FY1705
Air Cushion

Material: calcium carbonate \geq 30% rubber
Size: ϕ 35cm, ϕ 40cm, ϕ 42.5cm, ϕ 45cm, etc
Packing: 1pc/PE bag



FY1702
Fabric Ice Bag

Made of canvas
Size: 6", 9", 11", etc
Packing: 1pc/PE bag



FY1703
Instant Ice Bag

Material: PE+PA infilling (carbamide + distilled water)
Weight: 160g, 280g
Size: 19.5 \times 12cm, 24 \times 12cm



FY170404
Hot/Cold Pack

Material: nylon taffeta infilling (glycerol)
Weight: 200g
Size: 20.5 \times 15.5cm



FY170403
Hot/Cold Pack

Material: PVC+carbamide + distilled water
Weight: 40g
Size: 9.8×8cm




FY170401
Oval Ice Bag

Material: Rubber
Max Length: 28cm
Max Width: 17cm



FY170402
Hot/Cold Pack

Material: PVC+carbamide + distilled water
Weight: 65g
Size: 10×10cm



FY2301 Heat-Sealing Gusseted (Reel)

Size: (75+25)mm×100mm, 1roll/box, 5rolls/ctn, 22.5×46×23cm
(100+50)mm×100mm, 1roll/box, 4rolls/ctn, 22.5×47×23cm
(150+50)mm×100mm, 1roll/box, 2rolls/ctn, 22.5×34×23cm
(400+75)mm×100mm, non-box, 1roll/ctn, 22×41×22cm



FY2302 Self-sealing flat pouch

Size: 50×130mm 200pcs/box, 14400pcs/ctn, 270.6×37×30cm
90×230mm 200pcs/box, 4800pcs/ctn, 38.2×23×30cm
135×255mm 200pcs/box 3600pcs/ctn, 26.8×42.4×30cm
190×330mm 200pcs/box, 2400pcs/ctn, 34.3×29.6×30cm
250×370mm 200pcs/box, 1200pcs/ctn, 38.3×26.3×30cm



FY2001
Feeding Bottle

Material: PC/Glass+Liquid silicone nipple
Description: 150ml, 225ml, 250ml, 280ml




FY2002
Umbilical Cord Clamp

Material: ABS, PE, PP
Color: blue, white, etc
Size: 5cm, 5.5cm



FY2305 Paper/Paper Gusseted Pouch

Size: 90×50×150mm
90×50×250mm
110×30×190mm
140×75×250mm etc.



FY2304 Heat-Sealing Gusseted Pouch


Size: (75+25)mm×300mm, 100pcs/box, 3000pcs/ctn, 31.3×40×30cm
(100+50)mm×300mm, 100pcs/box, 2400pcs/ctn, 31.3×42.2×30cm
(150+50)mm×400mm, 100pcs/box, 1200pcs/ctn, 41.3×31.6×30cm
(200+55)mm×400mm, 100pcs/box, 600pcs/ctn, 21.3×41.3×30cm
(250+60)mm×450mm, 100pcs/box, 600pcs/ctn, 26.3×46.3×30cm



FY2003
Umbilical Cord Clamp Cutter




FY2004
Amniotomy Hook



FY2303 Heat-sealing flat reel pouch

Size: 55mm×200m 1roll/box, 6rolls/ctn, 22.5×42×23cm
75mm×200m 1roll/box, 5rolls/ctn, 2.5×46×23cm
100mm×200m 1roll/box, 4rolls/ctn, 22.5×47×23cm
150mm×200m 1roll/ctn, 2rolls/ctn, 22.5×34×23cm
400mm×200m non-box, 1roll/ctn, 22×41.5×22cm



FY2307 Medical Crepe Paper

Size: 30×30mm
50×50mm
70×70mm
120×120mm etc.



FY2306
Sterilization Indicator Tape

Model: EO, Autoclave and so on



FY2309
Steam Indicator Strip

Size: 15×100mm, 15×200mm, 200pcs/box, etc.



FY0619 Tourniquet

Easy to wipe clean
Highly flexible, latex free cotton ribbon which is gentle on the skin.
sturdy snap-on clip
quick and safe releasing




FY0621 Non-Latex Tourniquets

Material: TPE



FY2207
Electrosurgical Pad

Size: Child, Adult
Material: PE Foam or Non-woven, Aluminum-foil, Hydro-gel
Packing: 5pcs/pouch, 500pcs/carton
Carton size: 50×30×30cm,
G.W/N.W.: 16/15kg /carton




FY2206
Disposable E.C.G Accessories

Model : non-woven , foam , etc



FY0625 Pre-cut Roll Tourniquet

Material: TPE



FY0624 Latex Tourniquets

Material: TPE



FY2210 E.C.G Paper

Classification : single lead electrocardiograph ,
three lead electrocardiograph ,
six lead electrocardiograph ,
twelve lead electrocardiograph



FY3008 Disposable Oxi-Sensor

9 pin, Neonate/Adult
1pc/blister packing, 40pcs/box. 4boxes/ctn



FY2103 I.D. Bracelet

Material: PVC
Model: write-on, insert-card
Size: infant, adult, adult/infant
Color: blue, pink, white and so on

**FY180101**

Capacity: 0.8L
Product Size:
Top Dia: 125mm,
Bottom Dia: 95mm,
Height: 115mm
Carton Size: 510 x 390 x 570 mm
Quantity in one Carton : 180pcs

**FY180102**

Capacity: 3.0L
Product Size:
Top Dia: 170mm,
Bottom Dia: 150mm,
Height: 150mm
Carton Size: 520 x 520 x 460 mm
Quantity in one Carton : 90pcs

**FY180103**

Capacity: 3.0L
Product Size:
Top Dia: 200mm,
Bottom Dia: 160mm,
Height: 140mm
Carton Size: 620 x 340 x 450 mm
Quantity in one Carton : 70pcs

**FY180104**

Capacity: 2.0L
Product Size:
Top Dia: 200mm,
Bottom Dia: 120mm,
Height: 170mm
Carton Size: 405 x 270 x 450 mm
Quantity in one Carton : 30pcs

**FY180105**

Capacity: 1.0L
Product Size:
Top Dia: 110mm,
Bottom Dia: 90mm,
Height: 160mm
Carton Size: 360 x 320 x 460 mm
Quantity in one Carton : 100pcs

**FY180108**

Capacity: 3.0L
Product Size:
Top Dia: 150mm,
Bottom Dia: 120mm,
Height: 260mm
Carton Size: 480 x 320 x 480 mm
Quantity in one Carton : 30pcs

**FY180112
Needle Cutter**

Item: Needle Cutter
Material: ABS + Blade
Description:
Dia 11.33cmX9.6cm, manual
Packing: 60pcs/ctn,
48X36X49cm

**FY180109
Paper Safety Box**

75pcs/carton,
62x30x258.5cm,
20/18.5kgs

**FY1802
Biohazard Specimen Bag**

Size: 8"x10"
Thickness : 2mil 3ply
Size: 6"x10"
Thickness: 2mil 3ply

**AJJ**

HEALTHCARE MANAGEMENT PTE LTD

Our Promise To You.**Professional**

We set high standards with a high level of work ethic. Always be reliable, responsible, provide value and service promptly

Integrity

We always adhere to the highest level of medical ethics, based on a core of an organisational ecosystem of ethical culture

Transparency

We build and maintain trust by having honest and open communication and feedbacks with all stakeholders

Fairness

We ensure all stakeholders are treated equally and respectfully. Have zero tolerance for racism, sexism and xenophobia

Best Service Best Healthcare

AJJ Healthcare Management's core values of professional, integrity, fairness, transparency aligns with how we perceives all our stakeholders, especially the customers and patient. At AJJ, we believe that we can bring better healthcare for all by providing safe medical devices of superior quality for our customers. As a global medical device provider, AJJ will keep innovating and provide excellent services to all of our customers. As what we always say, it is all for the commitment to life.

Contact Us

24-Hour Customer Service Hotline
+65 6908 4252

