



AJ+

HEALTHCARE MANAGEMENT PTE LTD

Product Catalogue





CONTENTS

Surgical Clothing

Surgical Drapes

Surgical Packs

Surgical Accessories



All the material used for the production of disposable medical devices are certified for medical purpose.

Our products exceed all the requirements of EN 13795 and bear the CE label according to EC directive 94/92EEC medical devices.



SMS Non-woven (Spunbond-Meltblown-Spunbond)

3ply polypropylene composite fabric, comfortable to wear, with a strong barrier to fluids and particles.
With following additional options: Alcohol and fluids repellent/Anti-static

SMMMS Non-woven (Spunbond-Meltblown-Meltblown-Meltblown-Spunbond)

5ply polypropylene material composition, with 3 layers of Meltblown fabric, even more secure.
With following additional options: Alcohol and fluids repellent/Anti-static



Composite Non-woven

2 or more layers of materials are compounded together for different functions such as higher strength, better fluid absorption, barrier properties, etc.

Laminated Non-woven

Non-woven fabrics are coated with a layer of plastic film to achieve a better performance of fluid protection.



Spun-lace Non-woven

Non-woven material made by bonding process of wet or dry fibrous webs. Viscose, cotton, polyester and polypropylene are mainly used as raw materials.
Soft and absorbent, used for surgical hand towels, surgical caps, premium surgical gowns and sometimes as a reinforcement layer.

Material Legend:

The flexibility of our production structure enables us to develop products to meet specific demands of our customers. A variety of products and materials are used to adjust to different medical procedures, we offer a range of best cost-efficient solutions and surgical sets compositions.



Transparent PE

Transparent PE

Fluid repellent and resistant.

Non-woven Fabric

SMS, SMMMS, PP
Ensure high resistance and high protection to microorganisms and liquid penetration.



Fabric



Reinforcement

Reinforcement Layer

Hydrophobic and hydrophilic reinforcement of different materials.

Incise Film

Covering and protecting the surgical area to maintain sterility and provide skin protection.



Incise



Adhesive

Adhesive

Medical used, hypoallergenic and non-toxic.

Elastic material

Used in different surgical drapes to protect extremities.



Elastic material



Blue PE

Blue PE

Fluid repellent and resistant.



SURGICAL CLOTHING

Preventing direct microorganisms contact in the operating field is the fundamental function of the surgical gown: to decrease the transmission of skin flora from health care staff and to protect the staff against the blood-borne pathogens of the patient. Comparing to the traditional reusable gowns, disposable surgical gowns perform a stronger barrier between the medical staff and the patient as presenting a more solid, reliable and reproducible bacterial impermeability.

AJJ offers a wide range of disposable surgical gowns made by high quality non-woven fabrics to provide a more thorough protection as we as enhance the comfortableness and performance for the surgical team.

Protection and comfort

Not only light and soft, AJJ is always looking for:

More breathable and comfortable surgical clothing

-Customisable

A wide range of surgical clothing of different level of protection that can be used in surgeries of various departments.

-Performance

Perfectly block fine dust and liquid.

No flocculent particles.

Collar with velcro makes it convenient to put on and fix.

-Quality

Sterilised label on package to easily confirm if the product is well under non-bacteria condition.


Reinforced technique, durable protection

The production and quality inspection are in accordance with the European standards and American standard AAMI safety guidelines.

Standard Surgical Gown

1. Blue SMS material
2. Knitted cuff, four waist bands
3. Ultrasonic welding
4. Velcro on collar
5. Two hand towels and one SMS wrap
6. Sterilized/Non-sterile


REF:	Size	A	B	Quantity
FJY0001	M	120	140	50
FJY0002	L	120	150	50
FJY0003	XL	130	150	50
FJY0004	XXL	140	160	50



Reinforced Surgical Gown

1. Blue SMS material
2. Reinforcement on the front side and sleeves
3. Knitted cuff, four waist bands
4. Ultrasonic welding
5. Velcro on collar
6. Two hand towels and one SMS wrap
7. Sterilized/Non-sterile

REF:	Size	A	B	Quantity
FJY0006	M	120	140	50
FJY0007	L	120	150	50
FJY0008	XL	130	150	50
FJY0009	XXL	140	160	50




Material Legend: Fabric Reinforcement



Standard Surgical Gown

1. Blue spun-lace material
2. Knitted cuff, four waist bands
3. Knitted sleeve
4. Velcro on collar
5. Two hand towels and a wrap
6. Sterilized/Non-sterile


REF:	Size	A	B	Quantity
FJY0009	M	120	140	50
FJY0010	L	120	150	50
FJY0011	XL	130	150	50
FJY0012	XXL	140	160	50



Standard Surgical Gown with Raglan Sleeve

1. Blue SMS material, Raglan sleeve
2. Knitted cuff, four waist bands
3. Ultrasonic welding
4. Velcro on collar
5. Two hand towels and a wrap
6. Sterilized/Non-sterile


REF:	Size	A	B	Quantity
FJY0017	M	120	140	50
FJY0018	L	120	150	50
FJY0019	XL	130	150	50
FJY0020	XXL	140	160	50



Reinforced Surgical Gown

1. Blue spun-lace material
2. Reinforcement on the front side and sleeves
3. Knitted cuff, four waist bands
4. Knitted sleeve
5. Velcro on collar
6. Two hand towels and a wrap
7. Sterilized/Non-sterile


REF:	Size	A	B	Quantity
FJY0013	M	120	140	50
FJY0014	L	120	150	50
FJY0015	XL	130	150	50
FJY0016	XXL	140	160	50



Reinforced Surgical Gown with Raglan Sleeve

1. Blue SMS material, Raglan sleeve
2. Knitted cuff, four waist bands
3. Ultrasonic welding
4. Velcro on collar
5. Two hand towels and a wrap
6. Sterilized/Non-sterile
7. Reinforcement on the front side and sleeves

REF:	Size	A	B	Quantity
FJY0021	M	120	140	50
FJY0022	L	120	150	50
FJY0023	XL	130	150	50
FJY0024	XXL	140	160	50





Material Legend: Fabric Reinforcement

Material Legend: Fabric Reinforcement



Scrub Suit


1. Blue SMS material
2. Two pockets
3. Tie on waist
4. Packed in PE bag
5. Sterilized/Non-sterile



REF:	Size	A	B	C	D	E	Quantity
FSY0001	S	71	54	49	104	32	100
FSY0002	M	74	57	53	106	34.5	100
FSY0003	L	76	60	57	108	37	100
FSY0004	XL	78	63	64	113	39.5	100
FSY0005	XXL	80	65	64	115	41	100
FSY0006	3XL	82	67	65	118	42	100
FSY0007	4XL	84	70	70	121	43.5	100
FSY0008	5XL	86	73	73	124	45	100

Short Lab Coat

1. Blue SMS material
2. Two pockets
3. Cotton collar
4. Knitted cuffs
5. 4 snaps
6. Blue/green/pink
7. Packed in PE bag
8. Sterilized/Non-sterile




REF:(普通)	Size	A	B	C	D	Quantity
FSY0009	M	20	126	58	78	100
FSY0010	L	20	130	58	79	100
FSY0011	XL	21	135	58	83	100

Material Legend: Fabric

Long Lab Coat


1. Blue SMS material
2. Two pockets
3. Cotton collar
4. Knitted cuffs
5. 4 snaps
6. Packed in PE bag
7. Sterilized/Non-sterile



REF:(加长)	Size	A	B	C	D	Quantity
FSY0012	M	20	126	58	96	100
FSY0013	L	20	130	58	100	100
FSY0014	XL	21	135	58	105	100

Isolation Gown

1. PP material or PP+PE coating
2. Collar with tie
3. 1 waist band
4. Elastic cuffs/knitted cuffs
5. White/blue/green/yellow
6. Packed in PE bag
7. Sterilized/Non-sterile



REF:	Size	A	B	Quantity
FSY0015	L	120	140	100
FSY0016	XL	135	160	100


Material Legend: Fabric



Colonoscopy Pants

1. Dark blue PP
2. Elastic waist band
3. Split in the back with a flexible flap
4. Packed in PE bag
5. Sterilized/Non-sterile


REF:	Size	A	B	Quantity
FSY0017	M	57	66	100
FSY0018	L	60	70	100
FSY0019	XL	69	80	100



Patient Gown

1. Made by PP/SMS
2. 1 waist band
3. Half sleeves
4. Packed in PE bag
5. Sterilized/Non-sterile


REF:	Size	A	B	Quantity
FSY0020	L	120	140	100
FSY0021	XL	135	160	100



Patient Gown

1. Made by PP/SMS
2. 1 waist band
3. Collar with tie
4. Half sleeves
5. Packed in PE bag
6. Sterilized/Non-sterile

REF:	Size	A	B	Quantity
FSY0022	L	120	140	100
FSY0023	XL	135	160	100




Material Legend:  Fabric

Boot Cover

1. Made by PP+PE
2. Entry with elastic
3. Ties in the back
4. Consolidated by tape
5. Packed in PE bag
6. Sterilized/Non-sterile


REF:	Size	A	B	Quantity
FJX0001	L	45	37	500
FJX0002	XL	50	42	500



Boot Cover

1. Made by PP+PE
2. Elastic on top and in the middle part
3. Consolidated by tape
4. Packed in PE bag
5. Sterilized/Non-sterile


REF:	Size	A	B	Quantity
FJX0003	L	45	37	500
FJX0004	XL	50	42	500



Shoes Cover

1. PP material
2. Anti-slip/Normal
3. Packed in PE bag
4. Sterilized/Non-sterile

REF:	Size	A	B	Quantity
FSX0001		38	16	1000(pcs)
FSX0002		40	18	1000(pcs)




Material Legend:  Fabric




Nurse Cap/Bouffant Cap

REF:	Color	Size	Quantity
FSM0001	Blue	53cm	1000
FSM0002	Green	53cm	1000




Nurse Cap

REF:	Color	A	B	Quantity
FSM0003	Blue	33.5	14.5	1000
FSM0004	Green	33.5	14.5	1000




Doctor Cap

REF:	Color	A	B	Quantity
FSM0005	Blue	64	15	1000
FSM0006	Green	64	15	1000



Doctor Cap with Sweat Wicker


REF:	Color	A	B	Quantity
FSM0007	Blue	64	15	1000
FSM0008	Green	64	15	1000



Material Legend:  Fabric


CPE Apron

REF:	Color	A	B	Quantity
FSY0031	Blue	190	115	200



PE Apron


REF:	Color	A	B	Quantity
FSY0032	Blue	72	120	500
FSY0033	White	72	120	500



Material Legend:  Fabric



Coverall/Protective Suit



REF	Size	A	B	C	D	E	F	G	H
FSY0034	S	56	165	86	80	8~16	12~24	73	15
FSY0035	M	61	169	88	81	9~18	12.5~25	76	15
FSY0036	L	66	173	91	82	9~18	12.5~25	79	15
FSY0037	XL	71	180	94	84	9~18	12.5~25	84	15
FSY0038	2XL	76	185	97	85	9~18	12.5~25	88	15
FSY0039	3XL	81	188	100	86	9.5~19	13~26	92	15

Material Legend:  Adhesive tape  Fabric

SURGICAL DRAPE AND PACKS

AJJ provides a wide selection of both standard surgical drapes and customised drapes to meet our customers' needs.

Surgical drapes isolate the surgical area from non-sterile areas of the patient's body and the operating table, reducing the risk of germ transfer in surgeries.

Custom-surgical packs contain everything needed for the surgery, providing the medical staff and the patient with better safety solutions and greater efficiency in the surgery.

All packs are customised to your requirements.

Please contact us at +65 6908 4252 to customise your pack.

Our office address: 8 Commonwealth Lane, #02-04, Singapore 149555



Basic/Universal Drapes


Drape Sheet

Product Code	Size	Packaging unit	Packaging
FJD0001	45x45cm	1 pc.	700pcs
FJD0002	50x60cm	1 pc.	600pcs
FJD0003	60x75cm	1 pc.	400pcs
FJD0004	75x90cm	1 pc.	250pcs
FJD0005	90x120cm	1 pc.	200pcs
FJD0006	120x150cm	1 pc.	100pcs
FJD0007	150x180cm	1 pc.	100pcs




Side Drape with Adhesive

Product Code	Size	Packaging unit	Packaging
FJD0021	45x45cm	1 pc.	600pcs
FJD0022	50x60cm	1 pc.	500pcs
FJD0023	60x75cm	1 pc.	350pcs
FJD0024	75x90cm	1 pc.	200pcs
FJD0025	90x120cm	1 pc.	200pcs
FJD0026	120x150cm	1 pc.	100pcs
FJD0027	150x180cm	1 pc.	100pcs



Reinforced Side Drape

Product Code	Size	Packaging unit	Packaging
FJD0031	75x90cm	1 pc.	150pcs
FJD0032	90x120cm	1 pc.	100pcs
FJD0033	100x150cm	1 pc.	80pcs
FJD0034	150x200cm	1 pc.	60pcs




Material Legend:



Basic/Universal Drapes


Aperture Drape

Product Code	Size	Fenestration	Packaging unit	Packaging
FJD0051	45x45cm	ø6cm	1 pc.	600pcs
FJD0052	50x60cm	ø6cm	1 pc.	500pcs
FJD0053	60x75cm	ø7cm	1 pc.	350pcs
FJD0054	75x90cm	ø7cm	1 pc.	200pcs
FJD0055	90x120cm	ø8cm	1 pc.	200pcs
FJD0056	120x150cm	ø8cm	1 pc.	100pcs




Head Drape

Product Code	Size	Packaging unit	Packaging
FJD0041	150x240cm	1 pc.	50pcs
FJD0042	160x250cm	1 pc.	50pcs



Foot Drape

Product Code	Size	Packaging unit	Packaging
FJD0043	180x180cm	1 pc.	50pcs
FJD0044	180x200cm	1 pc.	50pcs



Material Legend:






Basic/Universal Drapes


Minor Procedure Drape

Product Code	Size	Fenestration	Packaging unit	Packaging
FJD0057	45x75cm	ø8cm	1 pc.	200pcs
FJD0058	75x90cm	7x9cm	1 pc.	150pcs



Reinforced Aperture Drape

Product Code	Size	Fenestration	Packaging unit	Packaging
FJD0061	475x90cm	ø7cm	1 pc.	200pcs
FJD0062	100x100cm	ø7cm	1 pc.	150pcs
FJD0063	100x120cm	ø8cm	1 pc.	120pcs
FJD0064	120x150cm	ø8cm	1 pc.	100pcs




Split Drape

Product Code	Size	Fenestration	Packaging unit	Packaging
FJD0071	120x130cm	10x50cm	1 pc.	80pcs
FJD0072	120x150cm	10x50cm	1 pc.	60pcs




Waterproof Split Drape

Product Code	Size	Fenestration	Packaging unit	Packaging
FJD0073	150x180cm	10x55cm	1 pc.	80pcs
FJD0074	150x210cm	10x55cm	1 pc.	60pcs




Material Legend:



Basic/Universal Drapes

Mayo Stand Cover

Product Code	Size	Packaging unit	Packaging
FJD0091	60x137cm	1 pc.	50pcs
FJD0092	80x145cm	1 pc.	50pcs



Back Table Cover

Product Code	Size	Packaging unit	Packaging
FJD0081	100x150cm	1 pc.	50pcs
FJD0082	140x190cm	1 pc.	50pcs
FJD0083	150x210cm	1 pc.	50pcs
FJD0084	150x230cm	1 pc.	50pcs




Leggings

Product Code	Size	Packaging unit	Packaging
FJD0101	45x60cm	2 pcs	100pcs



Leggings

Product Code	Size	Packaging unit	Packaging
FJD0102	75x120cm	2 pcs	100pcs
FJD0103	80x130cm	2 pcs	80pcs



Material Legend:






Basic/Universal Accessories

Stockinettes

Product Code	Size	Packaging unit	Packaging
FJQ0001	25x75cm	1 pc.	50pcs
FJQ0002	30x90cm	1 pc.	50pcs
FJQ0003	30x120cm	1 pc.	50pcs




OP Tape

Product Code	Size	Colour	Packaging unit	Packaging
FJQ0011	10x50cm	blue	1 pc.	800pcs
FJQ0012	10x50cm	green	1 pc.	800pcs



Hand Towels

Product Code	Size	Colour	Packaging unit	Packaging
FJQ0021	53g/m² 30x40cm	white	1 pc.	500pcs
FJQ0022	68g/m² 30x40cm	white	1 pc.	500pcs
FJQ0023	60g/m² 30x40cm	white	1 pc.	500pcs
FJQ0024	43x66cm	blue	1 pc.	100pcs



Baby Blanket

Product Code	Size	Packaging unit	Packaging
FJQ0025	60x80cm	1 pc.	50pcs
FJQ0026	75x90cm	1 pc.	50pcs



Material Legend:




Basic/Universal Accessories

Instrument Pouch/Diathermy Pouch

- With double tape
- Customizable size
- Made by transparent PE


Product Code	Size	Packaging unit	Packaging
FJD0111	35x40cm	1 pc.	200pcs
FJD0112	40x60cm	1 pc.	150pcs
FJD0113	35x90cm	1 pc.	100pcs



Fluid Collection Pouch

- With drainage ports and screen
- With double tape
- With malleable stripz


Product Code	Size	Packaging unit	Packaging
FJD0115	50x60cm	1 pc.	80pcs
FJD0116	96x57cm	1 pc.	50pcs



Knee Pouch

- With drainage ports and malleable strip
- With double tape
- 5cm circular fenestration in the middle of elastic material


Product Code	Size	Packaging unit	Packaging
FJD0117	107x92cm	1 pc.	50pcs



Knee Pouch

- With drainage ports and malleable strip
- With double tape
- 6cm circular fenestration in the middle of elastic material

Product Code	Size	Packaging unit	Packaging
FJD0118	96x80cm	1 pc.	50pcs



Material Legend:






Basic/Universal Accessories

Shoulder Pouch

- With drainage ports and malleable strip
- With double tape
- Circular fenestration in the middle of 20*25cm elastic material


Product Code	Size	Packaging unit	Packaging
FJD0121	80/90x80cm	1 pc.	50pcs



Shoulder Pouch

- With drainage ports and malleable strip
- With double tape


Product Code	Size	Packaging unit	Packaging
FJD0122	50/115x55cm	1 pc.	50pcs



Caesarean Pouch

- With drainage ports and malleable strip
- With double tape
- With medical incise area

Product Code	Size	Packaging unit	Packaging
FJD0119	72x75cm	1 pc.	50pcs
FJD0120	75x90cm	1 pc.	50pcs



Material Legend:




Basic/Universal Accessories

Circular Protective Cover

- Made by transparent PE
- Latex-free elastic bands


Product Code	Size	Packaging unit	Packaging
FJD0131	ø80cm	1 pc.	150pcs
FJD0132	ø100cm	1 pc.	100pcs
FJD0133	ø120cm	1 pc.	80pcs



Rectangular Protective Cover

- Made by transparent PE
- Latex-free elastic bands

Product Code	Size	Packaging unit	Packaging
FJD0135	75x75cm	1 pc.	150pcs
FJD0136	75x90cm	1 pc.	100pcs
FJD0137	100x100cm	1 pc.	80pcs



Surgical Camera Cover

- Made by transparent PE
- With adhesive



Product Code	Size	Packaging unit	Packaging
FJD0141	13x150cm	1 pc.	80pcs
FJD0142	13x250cm	1 pc.	80pcs


Material Legend:






UNIVERSAL PACK I REF:FJB0035

NO.	Product name	Size	Quantity
1	Head Drape	150x240cm	1
2	Foot Drape	175x180cm	1
3	Side Drape	75x90cm	2
4	Mayo Stand Cover	80x145cm	1
5	OP Tape	10x50cm	2
6	Hand Towels	30x40cm	2
7	Standard Surgical Gown	XL:130x150cm	2
8	Suture Bag	16x28cm	1
9	Back Table Cover	140x190cm	1
10	SMS Wrap	100x100cm	1




UNIVERSAL PACK III REF:FJB0037

NO.	Product name	Size	Quantity
1	Head Drape	150x250cm	1
2	Foot Drape	180x200cm	1
3	Side Drape	100x200cm	2
4	Mayo Stand Cover	80x145cm	1
5	Side Drape	40x65cm	4
6	OP Tape	10x50cm	2
7	Suture Bag	16x28cm	1
8	Hand Towels	30x40cm	6
9	Standard Surgical Gown	L:120x150cm	1
10	Reinforced Surgical Gown XL	XL:130x150cm	2
11	Back Table Cover	150x230cm	1
12	SMS Wrap	100x100cm	1



UNIVERSAL PACK II REF:FJB0036

NO.	Product name	Size	Quantity
1	Head Drape	150x240cm	1
2	Foot Drape	175x180cm	1
3	Side Drape	75x90cm	2
4	Mayo Stand Cover	80x145cm	1
5	OP Tape	10x50cm	2
6	Hand Towels	30x40cm	2
7	Suture Bag	16x28cm	1
8	Back Table Cover	140x190cm	1
9	SMS Wrap	100x100cm	1



Material Legend:



Material Legend:






EENT/Ophthalmology Surgery

Ophthalmic Drape 60x60cm

FDP01


- 6*9cm oval fenestration
- 10*15cm incise area
- 25*30cm fluid collection pouch with malleable strip



Ophthalmic Drape 130x140cm

FDP02


- 10*13cm rectangular fenestration
- 15*18cm incise area
- 25*30cm fluid collection pouch with malleable strip



Ophthalmic Drape 150x200cm

FDP03

- Rectangular incise area, 15*18 cm
- Rectangular fenestration, 9*12cm
- Two fluid collection pouches with malleable strip, 25*30cm




Material Legend:



Ophthalmic Drape 160x180cm

FDP04


- Rectangular incise area, 15*18cm
- Rectangular fenestration, 9*12cm
- Fluid collection pouch with malleable strip, 32*68cm



Head Turban Drape 100x100cm

FDP18


- Adhesive on top, 55cm
- Turban
- Glue on the dotted lines



Head Bar Drape 100x100cm

FDP19

- Adhesive on top, 55cm
- Glue on the dotted lines



Material Legend:






Thyroid T-Drape 200/250x315cm

FDP20


- Absorbent reinforcement: 60*75cm
- 12*12cm diamond fenestration with adhesive
- Tube holders



Craniotomy Drape 180/310x320cm

FDP34


- Absorbent reinforcement, 60*90cm
- Rounded block incise with aperture, 20*30cm
- Attached fluid collection pouch, 58*114cm
- Clear arm board covers, 70*100cm
- Tube holders



Thyroid Drape 200x300cm

FDP21


- Absorbent reinforcement: 60*75cm
- 12*12cm diamond fenestration with adhesive
- Tube holders



Craniotomy T-Drape 240x260cm

FDP35

- Absorbent reinforcement, 50*75cm
- Rounded block incise with aperture, 20*30cm
- Attached fluid collection pouch with screen and port, 58*114cm
- Opening of pouch with 3 malleable strips, 40*44cm



Material Legend:



Material Legend:






Craniotomy Drape 200x300cm


FDP60

- Absorbent reinforcement, 75*100cm
- Rounded block incise with aperture, 20*30cm
- Attached fluid collection pouch with screen and port, 58*114cm
- Tube holders




EENT PACK I REF:FJB0023

NO.	Product name	Size	Quantity
1	Split Drape	200x300cm	1
2	Head Bar Drape	100x100cm	1
3	Mayo Stand Cover	80x145cm	1
4	Hand Towels	30x40cm	4
5	Standard Surgical Gown	XL:130x150cm	2
6	Suture Bag	16x28cm	1
7	OP Tape	10x50cm	1
8	Back Table Cover	150x190cm	1
9	SMS Wrap	100x100cm	1




OPHTHALMIC PACK I REF:FJB0024

NO.	Product name	Size	Quantity
1	Ophthalmic Drape	130x140cm	1
2	Hand Towels	30x40cm	2
3	SMS Wrap	80x80cm	1



OPHTHALMIC PACK II REF:FJB0025

NO.	Product name	Size	Quantity
1	Ophthalmic Drape	100x150cm	1
2	Hand towels	30x40cm	2
3	Standard Surgical Gown	XL:130x150cm	2
4	Back table cover	100x150cm	1
5	SMS Wrap	100x100cm	1



Material Legend:




Material Legend:





OPHTHALMIC PACK III REF:FJB0026

NO.	Product name	Size	Quantity
1	Ophthalmic Drape	160x180cm	1
2	Mayo Stand Cover	60x137cm	1
3	OP Tape	10x50cm	1
4	Suture Bag	16x28cm	1
5	Hand Towels	30x40cm	2
6	Standard Surgical Gown	L:120x150cm	1
7	Standard Surgical Gown	XL:130x150cm	1
8	Back Table Cover	150x190cm	1
9	SMS Wrap	100x100cm	1




Extremity/Orthomax Surgery

Extremity Drape 200/250x315cm

FDP05


- Absorbent reinforcement, 75*100cm
- 6cm circular elastic fenestration
- Tube holders



Bilateral Extremity Drape 240/280x350cm


FDP06

- Absorbent reinforcement, 100*140cm
- Two 6cm elastic fenestration
- Tube holders



OPHTHALMIC PACK IV REF:FJB0027

NO.	Product name	Size	Quantity
1	Ophthalmic Drape	150x200cm	1
2	Hand Towels	30x40cm	2
3	Back Table Cover	100x150cm	1
4	SMS Wrap	100x100cm	1




Material Legend:



Material Legend:






Hand Drape 180/370x285cm

FDP47

- Absorbent reinforcement: 50*100cm
- 3cm circular fenestration
- Tube holders



Orthomax Body Split Drape 220x290cm

FDP15

- Absorbent reinforcement, 60*75cm
- 13*70cm split with 5cm adhesive
- PE reinforcement, 40*50cm


Material Legend:



Bilateral Extremity Drape 220x320cm

FDP08


- Absorbent reinforcement, 75*100cm
- 6cm circular elastic fenestration
- Tube holders



Extremity Drape 200x300cm

FDP59

- Absorbent reinforcement, 50*100cm
- 6cm circular elastic fenestration
- Tube holders



Material Legend:







Orthomax U-Drape 200x300cm

FDP17

- Absorbent reinforcement, 75*90cm
- Split with 5cm adhesive, 7.5*65cm
- Tube holders



Orthomax Body Split Drape 220x290cm

FDP33


- Absorbent reinforcement, 75*100cm
- Split with adhesive, 13*70cm
- Tube holders

Material Legend:




EXTREMITY PACK REF:FJB0012

Quantity	Product name	Size
1 pc.	Extremity T-Drape	200/250x315cm
1 pc.	Split Drape	150x180cm
1 pc.	Stockinette	30x90cm
1 pc.	Mayo Stand Cover	80x145cm
2 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	L:120x150cm
2 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	XL:130x150cm
2 pcs	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Video Camera Cover	13x250cm
2 pcs	Back Table Cover	140x190cm
1 pc.	SMS Wrap	100x100cm




EXTREMITY PACK REF:FJB0013

Quantity	Product name	Size
1 pc.	Extremity Drape	200x300cm
2 pcs	Side Drape	100x150cm
1 pc.	Legging	45x60cm
1 pc.	Mayo Stand Cover	80x145cm
2 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	L:120x150cm
2 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	XL:130x150cm
2 pcs	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Video Camera Cover	13x250cm
2 pcs	Back Table Cover	140x190cm
1 pc.	SMS Wrap	100x100cm




Material Legend:






BILATERAL EXTREMITY PACK REF:FJB0014

Quantity	Product name	Size
1 pc.	Bilateral Extremity Drape	240/280x350cm
2 pcs	Side Drape	100x150cm
2 pcs	Legging	45x60cm
1 pc.	Mayo Stand Cover	80x145cm
2 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	L:120x150cm
2 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	XL:130x150cm
4 pcs	OP Tape	10x50cm
2 pcs	Suture Bag	16x28cm
1 pc.	Skin Marker and Ruler	ST
1 pc.	Video Camera Cover	13x250cm
2 pcs	Back Table Cover	140x190cm
1 pc.	SMS wrap	100x100cm




BILATERAL EXTREMITY PACK REF:FJB0015

Quantity	Product name	Size
1 pc.	Bilateral Extremity Drape	220x320cm
2 pcs	Split Drape	150x180cm
2 pcs	Legging	45x60cm
1 pc.	Mayo Stand Cover	80x145cm
2 pcs	Hand towels	30x40cm
1 pc.	Reinforced Surgical Gown	L:120x150cm
2 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	XL:130x150cm
4 pcs	OP Tape	10x50cm
2 pcs	Suture Bag	16x28cm
1 pc.	Skin Marker and Ruler	ST
1 pc.	Video Camera Cover	13x250cm
2 pcs	Back table cover	140x190cm
1 pc.	SMS wrap	100x100cm




Material Legend:




HAND PACK REF:FJB0016

Quantity	Product name	Size
1 pc.	Hand Drape	180/370x285cm
1 pc.	Drape Sheet	150x150cm
1 pc.	Mayo Stand Cover	80x145cm
4 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	L:120x150cm
1 pc.	Reinforced Surgical Gown	XL:130x150cm
2 pcs	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Back Table Cover	140x190cm
1 pc.	SMS wrap	100x100cm



Material Legend:






Arthroscopy Surgery

Arthroscopy T-Drape with Fluid Control Pouch

FDP09


- Absorbent reinforcement, 75*100cm
- Circular elastic fenestration, 6cm
- Attached fluid collection pouch with malleable strip, 96*82cm
- Tube holders



Arthroscopy Drape with Fluid Control Pouch 200x300cm

FDP10


- Absorbent reinforcement, 75*100cm
- Circular elastic fenestration, 6cm
- Attached fluid collection pouch, 107*92cm
- Tube holders



Arthroscopy T-Drape with Fluid Control Pouch

FDP11


- Absorbent reinforcement, 75*100cm
- Two circular elastic fenestration, 6cm
- Attached fluid collection pouches with malleable strip, 96*82cm
- Tube holders



Arthroscopy T-Drape with Fluid Control Pouch

FDP12


- Absorbent reinforcement, 75*100cm
- Circular elastic fenestration, 6cm
- Attached fluid collection pouch, 107*92cm
- Tube holders



Material Legend:




Material Legend:






KNEE ARTHROSCOPY PACK REF:FJB0007

Quantity	Product name	Size
1 pc.	Knee Arthroscopy T-Drape	200/250x315cm
1 pc.	Split Drape	150x180cm
1 pc.	Stockinette	30x90cm
1 pc.	Mayo Stand Cover	80x145cm
2 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	L:120x150cm
2 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	XL:130x150cm
1 pc.	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Yankauer Suction Tube	23+304cm
1 pc.	Back Table Cover	140x190cm
1 pc.	SMS Wrap	100x100cm




KNEE ARTHROSCOPY PACK REF:FJB0009

Quantity	Product name	Size
1 pc.	Knee Arthroscopy T-Drape	200/250x315cm
1 pc.	Split Drape	150x180cm
1 pc.	Mayo Stand Cover	80x145cm
1 pc.	Stockinette	30x90cm
2 pcs	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Skin Marker and Ruler	ST
1 pc.	Video Camera Cover	13x250cm
1 pc.	Back Table Cover	140x190cm
1 pc.	SMS wrap	100x100cm




KNEE ARTHROSCOPY PACK REF:FJB0008


Quantity	Product name	Size
1 pc.	Knee Arthroscopy Drape	200x300cm
1 pc.	Split Drape	150x180cm
1 pc.	Mayo Stand Cover	60x137cm
1 pc.	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Legging	45x60cm
1 pc.	Sheet Drape	140x190cm
1 pc.	Back Table Cover	140x190cm
1 pc.	SMS Wrap	100x100cm



Material Legend:



Material Legend:






Laparotomy/Laparoscopy Surgery

Chest U-Drape 200x260cm

FDP13


- Absorbent reinforcement, 75*90cm
- 30*75cm U Split
- Tube holders



Nephroscopy Drape with Pouch 180x300cm

FDP53


- Absorbent reinforcement with built-in tubing system
- PE reinforcement, 50*60cm
- Oval adhesive fenestration with 5cm circular aperture, 20*30cm
- Fluid collection pouch with malleable strip, drainage port and screen, 65*110cm
- Four tube holders, 7.5*13cm



Laparotomy Drape 190x240cm

FDP16


- Absorbent reinforcement: 40*80cm
- 10*30cm rectangular fenestration with 2.5cm adhesive



Laparoscopy/Pelviscopy Drape with gussets 100/250x240cm

FDP69

- 18/23*18cm trapezoidal abdominal fenestration with adhesive
- 7.5*15cm perineal fenestration
- Privacy flap, 36*40cm
- Two tube holders
- Absorbent reinforcement, 35*75cm



Material Legend:



Material Legend:






Laparoscopic Perineal Drape 100/250*240cm

FDP70


- 18/23*18cm trapezoidal abdominal fenestration with adhesive
- 7.5*18cm perineal fenestration
- Privacy flap, 36*40cm
- Absorbent reinforcement, 50*75cm
- 62*75cm fluid collection pouch with screen and drainage ports
- 6 tube holders



Laparotomy T-Drape 200/260x315cm

FDP73


- Absorbent reinforcement, 75*75cm
- 10*30cm fenestration with 2.5cm adhesive
- 4 tube holders



Laparotomy T-Drape 200/260x315cm

FDP74


- Absorbent reinforcement, 75*75cm
- 10*30cm fenestration with 2.5cm adhesive
- 4 tube holders



Laparoscopic Cholecystectomy Drape with troughs 200/260x315cm

FDP75

- Absorbent reinforcement, 50*90cm
- Two open-ended, poly troughs, 25*75cm
- 30*30cm fenestration with 2.5*35cm adhesive
- 4 tube holders



Material Legend:




Material Legend:






LAPAROSCOPIC PERINEAL PACK REF:FJB0003

Quantity	Product name	Size
1 pc.	Laparoscopic Perineal Drape	100/250x240cm
2 pcs	Leggings	75x120cm
1 pc.	Mayo Stand Cover	80x145cm
4 pcs	Hand Towels	30x40cm
2 pcs	Standard Surgical Gown	XL:130x150cm
1 pc.	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Back Table Cover	140x190cm
1 pc.	SMS Wrap	100x100cm




LAPAROTOMY PACK REF:FJB0005

Quantity	Product name	Size
1 pc.	Laparotomy Drape	200/260x315cm
2 pcs	Drape Sheet	140x160cm
1 pc.	Mayo Stand Cover	60x138cm
4 pcs	Hand Towel	30x40cm
2 pcs	Reinforced Surgical Gown	XL:130x150cm
1 pc.	Suture Bag	16x28cm
4 pcs	Side Drape	45x60cm
1 pc.	Back Table Cover	150x210cm
1 pc.	SMS Wrap	100x100cm



LAPAROSCOPY / PELVICOPY PACK REF:FJB0004

Quantity	Product name	Size
1 pc.	Laparoscopy / Pelviscopy Drape	100/250x240cm
1 pc.	Under Buttocks Drape	90x100cm
2 pcs	Leggings	75x120cm
1 pc.	Mayo Stand Cover	80x145cm
4 pcs	Hand Towels	30x40cm
2 pcs	Reinforced Surgical Gown	XL:130x150cm
1 pc.	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Drape Sheet	75x90cm
1 pc.	Back Table Cover	150x210cm
1 pc.	SMS Wrap	100x100cm



Material Legend:



Material Legend:






Cardiovascular/Neurology Surgery

Femoral Angiography Drape 220x300cm

FDP22


- Two circular adhesive fenestration, 15cm
- Aperture, 9cm
- Central absorbent reinforcement, 75*100cm
- Clear side panels on both sides of the patient



Femoral Angiography Drape 210x300cm

FDP23


- Two circular adhesive fenestration, 15 cm
- Aperture, 9cm
- Central absorbent reinforcement, 75*100cm
- One clear side panel



Femoral Angiography Drape 220x360cm

FDP24


- Absorbent reinforcement: 75x120cm
- Two 15*15cm rectangular incise area
- Two circular adhesive fenestration, 15cm
- Aperture: 9cm
- Two clear side panels



Femoral Angiography Drape extra long 220x400cm

FDP25


- Absorbent reinforcement, 80*100cm
- Two circular adhesive fenestration, 15cm
- Aperture, 9cm
- 25*75cm fluid collection pouch with malleable strip and velcros
- Two clear side panels



Material Legend:



Material Legend:






Radial Angiography Drape 240x380cm

FDP26


- Two 12cm radial adhesive fenestration, with 6cm aperture
- Two 15cm femoral adhesive fenestration, with 9cm aperture
- Two clear side panels on both sides of patient
- Central absorbent reinforcement, 75*170cm



Pacemaker Drape 220x290cm

FDP27

- Two 15*18 cm rectangular adhesive fenestration, with 10*10cm aperture
- Absorbent reinforcement, 75*90 cm
- Catheter pouch with malleable strip, 30*30cm




Material Legend:



Central Line Drape 180x250cm

FDP28


- Absorbent reinforcement: 50*50cm
- 10cm circular adhesive fenestration
- 25*25cm clear PE window around fenestration



Central Line Drape 200x280cm

FDP29

- Absorbent reinforcement: 50*110cm
- Two 12*15cm rectangular adhesive fenestration
- 30*25cm clear PE window around each fenestration



Material Legend:






Central Line Drape 180x250cm

FDP30


- Circular fenestration with adhesive (15*15cm) around, 10cm
- Clear PE window around fenestration, 25*25cm
- Absorbent reinforcement, 50*50cm



Central Line Drape 200x280cm

FDP31


- Absorbent reinforcement, 50*110cm
- Two 12*15cm rectangular adhesive fenestration
- 2.5cm adhesive around each fenestration
- 30*25cm clear PE window around each fenestration



Cardiovascular Split Drape 217.5/260x340cm

FDP57


- Split with adhesive, 7.5*200cm
- 35*63.5cm top reinforcement with 50*63.5cm flap
- 2 split reinforcement, 22.5*150cm



Bilateral Split Drape 200x200cm

FDP58

- Absorbent reinforcement, 50*120cm
- Two U splits with 2.5cm adhesive, 15*50cm
- Groin towel with adhesive, 90*30cm



Material Legend:



Material Legend:






Cardiovascular Drape 200/300x300cm

FDP61


- Fenestration with incise
- 50*120cm reinforcement
- Two 32*75cm fluid collection pouches



Cardiovascular U Split Drape 215/275x335cm

FDP62

- 15*210cm U split with adhesive
- Split reinforcement with hook-and-loop closures
- Clear anesthesia screen
- Contoured armboard covers




Material Legend:



Universal Spine Drape 240/280x350cm

FDP76

- Absorbent reinforcement, 50*100cm
- 10*30cm fenestration with 2.5*35cm adhesive
- 220*28cm instrument pouches
- 4 tube holders




Material Legend:






ANGIOGRAPHY PACK REF:FJB0019

NO.	Product name	Size	Quantity
1	Angiography Drape	210x300cm	1
2	Hand Towels	30x40cm	2
3	Reinforced Surgical Gown	XL:130x150cm	2
4	Rectangular Banded Bag	75x90cm	1
5	Circular Banded Bag	Φ100cm	1
6	Solution Cup	500cc	1
7	X-Ray Gauze	7.5x7.5cm	10
8	Betadine Sponge	/	1
9	Gallipots	1000cc	1
10	Guid Wire	150cm	1
11	Back Table Cover	150x190cm	1
12	SMS Wrap	100x100cm	1



ANGIOGRAPHY PACK REF:FJB0020

NO.	Product name	Size	Quantity
1	Angiography Drape	240x380cm	1
2	Circular Banded Bag	Φ120cm	1
3	Rectangular Banded Bag	75x90cm	1
4	X-Ray Gauze	7.5x7.5cm	10
5	Hand Towels	30x40cm	4
6	Reinforced Surgical Gown	XL:130x150cm	2
7	Betadine Sponge	/	1
8	Solution Cup	500cc	1
9	Gallipots	1000cc	1
10	Guid Wire	150cm	1
11	Back Table Cover	140x190cm	1
12	SMS Wrap	100x100cm	1




Material Legend:




CARDIOVASCULAR PACK REF:FJB0021

NO.	Product name	Size	Quantity
1	Cardiovascular U Split Drape	215x75x340cm	1
2	Bilateral Split Drape	200x200cm	1
3	Side Drape	75x90cm	1
4	Mayo Stand Cover	80x145cm	1
5	Hand towels	30x40cm	4
6	Reinforced Surgical Gown	XL:130x150cm	1
7	OP Tape	10x50cm	2
8	Suture Bag	16x28cm	2
9	Back table cover	150x230cm	1
10	Back table cover	140x190cm	1
11	SMS wrap	100x100cm	1



CARDIOVASCULAR PACK REF:FJB0022

NO.	Product name	Size	Quantity
1	Cardiovascular U Split Drape	218/260x340cm	1
2	Bilateral Split Drape	200x200cm	1
3	Top Drape	150x240cm	1
4	Mayo Stand Cover	80x145cm	1
5	Hand towels	30x40cm	2
6	Shoe Covers	40x20cm	2
7	PE Shoe Covers	40x20cm	2
8	OP Tape	10x50cm	2
9	Suture Bag	16x28cm	2
10	Side Drape	75x90cm	1
11	Back Table Cover	150x230cm	1
12	Back Table Cover	140x190cm	1
13	SMS wrap	100x100cm	1



Material Legend:






Obstetrics Surgery

Caesarean Section Drape 160/ 230x300cm

FDP36


- 360°fluid collection pouch, 72*72cm
- Block incise area, 30*30cm
- Tube holders



Caesarean Section Drape 200/ 260x320cm

FDP37


- Absorbent reinforcement:25*50cm
- Block incise area:35*35cm
- 360°fluid collection pouch:75*90cm



Caesarean Section Drape 160x300cm

FDP38


- Block incise, 30*30cm
- 360°fluid collection pouch with malleable strip, 75*80cm
- Tube holders



Caesarean Section Drape 200x320cm

FDP39

- Block incise, 30*30cm
- 360°fluid collection pouch with malleable strip, 75*80cm
- Absorbent reinforcement, 25*50cm
- Tube holders



Material Legend:



Material Legend:






Under Buttocks Drape 100x112cm

FDP40


- Absorbent reinforcement, 25*50cm
- Fluid collection pouch, 82*74cm with filtering screen, drainage port and malleable strip



Under Buttocks Drape 85x36cm

FDP41


- Absorbent reinforcement: 25*50cm
- Fluid collection pouch: 82*74cm with filtering screen, drainage port and malleable strip



Under Buttocks Drape 90x100cm

FDP42


- Reinforcement: 35*50cm



Under Buttocks Drape 150x150cm

FDP43

- Reinforcement: 60*50cm
- Fluid collection pouch: 82*74cm with filtering screen, drainage port and malleable strip



Material Legend:



Material Legend:






Caesarean Section Drape 160x300cm

FDP64


- 75*80cm Fluid collection pouch with malleable strip
- 30*30cm incise area
- 18*18cm pear-shaped fenestration
- 3 tube holders



Caesarean Section Drape 160/230x300cm

FDP65

- 75*80cm Fluid collection pouch with malleable strip
- 30*30cm incise area
- 18*18cm pear-shaped fenestration
- 5 tube holders
- Contoured armboard covers




Material Legend:




CAESAREAN SECTION PACK REF:FJB0030

NO.	Product name	Size	Quantity
1	Caesarean Section Drape	160/230x300cm	1
2	Mayo Stand Cover	60x137cm	1
3	Hand Towels	30x40cm	2
4	Umbilical Clamp		1
5	Baby Blanket	60x80cm	1
6	Back Table Cover	150x190cm	1
7	SMS Wrap	100x100cm	1



CAESAREAN SECTION PACK REF:FJB0031

NO.	Product name	Size	Quantity
1	Caesarean Section Drape	200/260x320cm	1
2	Mayo Stand Cover	80x145cm	1
3	Hand Towels	30x40cm	4
4	Reinforced Surgical Gown	XL:130x150cm	2
5	Umbilical Clamp		1
6	Baby Blanket	75x90cm	1
7	Suture Bag	16x28cm	1
8	Back Table Cover	150x190cm	1
9	SMS Wrap	100x100cm	1




Material Legend:






CAESAREAN SECTION PACK REF:FJB0032

NO.	Product name	Size	Quantity
1	Caesarean Section Drape	160x300cm	1
2	Mayo Stand Cover	80x145cm	1
3	Hand Towels	30x40cm	4
4	Reinforced Surgical Gown	XL:130x150cm	2
5	Umbilical Clamp		1
6	Baby Blanket	75x90cm	1
7	Suture Bag	16x28cm	1
8	Back Table Cover	140x190cm	1
9	SMS Wrap	100x100cm	1




DELIVERY PACK REF:FJB0034

NO.	Product name	Size	Quantity
1	Under Buttocks Drape	85x36cm	1
2	Drape Sheet	100x100cm	1
3	Leggings	75x120cm	2
4	Side Drape	75x90cm	1
5	Hand Towels	30x40cm	2
6	Reinforced Surgical Gown	XL:130x150cm	1
7	Baby Blanket	75x90cm	1
8	Umbilical Clamp		1
9	Back Table Cover	140x190cm	1
10	SMS Wrap	100x100cm	1



DELIVERY PACK REF:FJB0033


NO.	Product name	Size	Quantity
1	Under Buttocks Drape	100x112cm	1
2	Leggings	75x120cm	2
3	Reinforced Side Drape	100x150cm	1
4	Hand Towels	30x40cm	2
5	Reinforced Surgical Gown	XL:130x150cm	1
6	Baby Blanket	75x90cm	1
7	Umbilical Clamp		1
8	Back Table Cover	150x190cm	1
9	SMS Wrap	100x100cm	1



Material Legend:



Material Legend:






Hip/Urology Surgery

Hip Drape 240x350cm

FDP50


- Absorbent reinforcement: 75*140cm
- 17*20cm pear-shaped elastic fenestration
- Four tube holders
- 5*18cm adhesive



Hip Drape 220/280x350cm

FDP51


- Absorbent reinforcement: 60*150cm
- Elastic area: 35*40cm
- 17*20cm pear-shaped fenestration
- Fluid collection pouch: 75*150cm with 4 velcros
- Two tube holders



By-pass Drape 160x240cm

FDP63


- 44*90cm U split with 5cm adhesive
- One I-shaped perineal flap with double tape
- Absorbent reinforcement: 75*140cm



TUR Drape 170/240x175cm

FDP66


- 6cm top circular fenestration with incise
- 3.5cm bottom circular fenestration with finger-cot
- 20*20cm top absorbent reinforcement
- 50*60cm bottom absorbent reinforcement
- 96*57cm fluid collection pouch with screen and drainage ports



Material Legend:




Material Legend:






T.U.R. PACK REF:FJB0006

Quantity	Product name	Size
1 pc.	T.U.R. Drape	170/240x175cm
4 pcs	Hand Towels	30x40cm
1 pc.	Standard Surgical Gown	XL:130x150cm
1 pc.	Reinforced Surgical Gown	XL:130x150cm
1 pc.	Back Table Cover	140x190cmv
1 pc.	SMS Wrap	100x100cm



HIP PACK REF:FJB0010

Quantity	Product name	Size
1 pc.	Hip T-Drape	220/280x350cm
1 pc.	Split Drape	150x210cm
1 pc.	Stockinette	30x120cm
1 pc.	Mayo Stand Cover	80x145cm
2 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	L:120x150cm
2 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	XL:130x150cm
1 pc.	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Yankauer Suction Tube	23x304cm
1 pc.	Back Table Cover	150x230cm
1 pc.	SMS Wrap	100x100cm



Material Legend:




Shoulder/Orthopedic Surgery

Shoulder Drape with pouch 160x250cm

FDP44


- Fluid collection pouch with drainage ports, 50/115*55cm
- Shoulder split, 5*50cm
- Adhesive around split



Shoulder Drape with pouch 160x250cm

FDP45

- Fluid collection pouch with drainage ports, 50/115*55cm
- Adhesive fenestration, 28*35.6cm
- Pear shaped aperture, 18*25cm



Material Legend:






Shoulder Drape 240x400cm

FDP46


- Fluid collection pouch 80*90cm with drainage port
- 10*11cm sector shaped fenestration with adhesive in the back
- Two velcros



Shoulder Drape 160x250cm

FDP48


- Absorbent reinforcement, 50*75cm
- Contoured shoulder split with 5cm adhesive, 7.5/15*20cm
- PE area, 30*30cm
- Tube holders



Shoulder Drape 160x280cm

FDP49


- Absorbent reinforcement, 50*120cm
- Two contoured shoulder split with
- 5cm adhesive, 5/9*15cm
- Tube holders



Patient Isolation Drape with pouch 320x200cm

FDP52

- Adhesive rectangular fenestration, 24*50cm
- Fluid collection pouch with malleable strip, 50*60cm
- Two instrument pouches, 18*30cm
- Two tube holders, 15*20cm
- Adhesive strip along top



Material Legend:




Material Legend:






SHOULDER PACK REF:FJB0017

Quantity	Product name	Size
1 pc.	Shoulder Drape with Pouch	240x400cm
1 pc.	Split Drape	150x180cm
1 pc.	Stockinette	25x75cm
1 pc.	Mayo Stand Cover	80x145cm
4 pcs	Hand Towels	30x40cm
2 pcs	Reinforced Surgical Gown	XL:130x150cm
2 pcs	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Back table cover	140x190cm
1 pc.	SMS wrap	100x100cm



SHOULDER PACK REF:FJB0018

Quantity	Product name	Size
1 pc.	Shoulder Drape with Pouch	160x250cm
1 pc.	Orthomax Body Split Drape	220x300cm
1 pc.	Mayo Stand Cover	80x145cm
1 pc.	Stockinette	25x75cm
4 pcs	Hand Towels	30x40cm
2 pcs	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Split Drape	150x180cm
1 pc.	Back table cover	140x190cm
1 pc.	SMS wrap	100x100cm




Material Legend:



ORTHOPEDIC PACK REF:FJB0011

Quantity	Product name	Size
1 pc.	Orthopedic Split Drape	220x300cm
1 pc.	Split Drape	150x210cm
1 pc.	Head Drape	150x240cm
1 pc.	Mayo Stand Cover	80x145cm
4 pcs	Hand Towels	30x40cm
1 pc.	Reinforced Surgical Gown	L:120x150cm
1 pc.	Reinforced Surgical Gown	XL:130x150cm
1 pc.	OP Tape	10x50cm
1 pc.	Suture Bag	16x28cm
1 pc.	Instrument Pouches	35x40cm
1 pc.	Stockinette	30x90cm
1 pc.	Back Table Cover	140x190cm
1 pc.	SMS Wrap	100x100cm



Material Legend:






Dental Surgery

Dental Drape 牙科单 100x150cm

FDP54


- 10.8*13cm fenestration, 40cm to the top



Dental Drape 牙科单 100x150cm

FDP55

- Teardrop fenestration 10*15cm, 40cm to the top
- Absorbent reinforcement: 30*40cm




Material Legend:



Dental Drape 牙科单 100x150cm

FDP56

- U split, 10*10cm
- Adhesive on top, 2.5cm




Material Legend:






DENTAL PACK REF:FJB0028

NO.	Product name	Size	Quantity
1	Dental Drape	100x150cm	1
2	Drape Sheet	50x60cm	2
3	Aperture Drape	75x90cm	1
4	Gauze	5x5-4p	10
5	Video Camera Cover	8x120cm	2
6	OP Tape	5x50cm	2
7	Suture Bag	16x28cm	1
8	Hand Towel	30x40cm	2
9	Standard Surgical Gown	L:120x150cm	2
10	Back Table Cover	140x190cm	1
11	SMS Wrap	100x100cm	1



DENTAL PACK REF:FJB0029

NO.	Product name	Size	Quantity
1	Dental Drape	100x150cm	1
2	Drape Sheet	75x90cm	1
3	Gauze	5x5-4p	10
4	Video Camera Cover	8x120cm	2
5	OP Tape	5x50cm	2
6	Suture Bag	16x28cm	1
7	Hand Towels	30x40cm	4
8	Standard Surgical Gown	L:120x150cm	2
9	Back Table Cover	140x190cm	1
10	SMS Wrap	100x100cm	1



Material Legend:




Lithotomy /Cystoscopy Surgery

Lithotomy Drape 80/160x270cm

FDP67


- Anesthesia screen, 160*100cm
- Perineal fenestration with adhesive, 7.5*15cm
- With leggings
- Absorbent reinforcement, 45*50cm
- Fluid collection pouch with screen and drainage ports, 70*75cm



Lithotomy Drape 75/135x195cm

FDP68

- 7.5*15cm fenestration with 2.5cm adhesive
- Absorbent reinforcement, 30*60cm
- Tube holders



Material Legend:






Cystoscopy Drape 75x180cm

FDP71


- 5.5*12.5cm diamond fenestration
- 2.5*10cm adhesive
- Two 21*31cm U split with four 2.5*3cm adhesive



Cystoscopy Drape 90/150x215cm

FDP72

- 7.5*18cm fenestration with 2.5cm adhesive
- 25*25cm fenestration with mesh screen
- Film reinforcement, 45*100cm
- Two tube holders




Material Legend:



CYSTOSCOPY PACK

REF:FJB0001


Quantity	Product name	Size
1 pc.	Cystoscopy "T" Drape	90/150x215cm
2 pcs	Leggings	75x120cm
4 pcs	Hand Towels	30x40cm
1 pc.	Standard Surgical Gown	XL:130x150cm
1 pc.	Reinforced Surgical Gown	XL:130x150cm
1 pc.	Back Table Cover	140x190cm
1 pc.	SMS wrap	100x100cm



CYSTOSCOPY PACK

REF:FJB0002

Quantity	Product name	Size
1 pc.	Cystoscopy Drape	75x180cm
2 pcs	Leggings	75x120cm
2 pcs	Hand Towels	30x40cm
1 pc.	Standard Surgical Gown	XL:130x150cm
1 pc.	OP Tape	10x50cm
1 pc.	Drape Sheet	90x90cm
1 pc.	Back Table Cover	140x190cm
1 pc.	SMS Wrap	100x100cm



Material Legend:

