


Effectively prevention of cross-infection

The fiber-optical tips and
Stainless handpiece sleeve can be autoclaved.



Large-capacity battery

11.1V 2600mAh x2 (57.7Wh) lithium battery,
which can be used for one week on a single charge.



Complete accessories

Fiber-optical tips,
TMJ therapy tip,
whitening tip,
biostimulation tip,
laser protective glasses



Effective prevention

Three pairs of laser protective glasses can effectively prevent the dentist,
assistant and patient from laser radiation.



Metal handle enables flexible movements.

Laser parameters

- Wavelengths & optical power:
 - 450 + 20 nm : Pmax=3W
 - 650 + 20 nm : Pmax=200mW
 - 976 + 20 nm : Pmax=5W
- Laser system:
 - 450 nm: Class IV;
 - 976 nm: Class IV;
 - 650 nm: Class II
 (according to IEC 60825-1)
- Emission modes:
 - CW (continuous wave), chopped 1 Hz to 20 kHz
- Aiming beam:650 + 20 nm / Pmax<5mW
- Battery:11.1V 2600mAh x2 (57.7Wh)

C-SMART-I PRO ENDO MOTOR



Model: C-SMART-I PRO



Order name: C-SMART-I PRO

Technical:


Adaptor working voltage	AC100-240V,50/60HZ;
Adjustable torque,range	0.6-5.2N.cm(6-52mN.m)
5 setting for adjustable rotational speed ranging:	150rpm-650rpm
Volume(cm)	23.5X20X14.5
Weight(kg)/PCS	2.8
Packing material	carton
QTY per master carton	6 pc/ctn

Touch control, System inside
With LED light
With Apex locator function

Built-in mainstream system



Endo motor & Apex locator
Dual mode



Big color LCD screen



Model: C-SMART-MINI

Order name: C-SMART-MINI

Technical:


Adaptor working voltage	100-240Vac,50/60HZ;
Adjustable torque,range	0.6-3.9N.cm(6-39mN.m)
5 setting for adjustable rotational speed ranging:	125rpm-625rpm
Volume(cm)	19X13X6
Weight(kg)/PCS	0.7
Packing material	carton
QTY per master carton	12 pc/ctn




Order name: C-SMART-MINI (LED)

Technical:

Adaptor working voltage	AC100-240V,50/60HZ;
Adjustable torque,range	0.6-3.9N.cm(6-39mN.m)
5 setting for adjustable rotational speed ranging:	125rpm-625rpm
Volume(cm)	19X13X6
Weight(kg)/PCS	0.7
Packing material	carton
QTY per master carton	12 pc/ctn



LED light



Reciprocation



Reciprocation



Non-contact recharge design



LED screen



Non-contact recharge design