

Supragingival scaling

Implant maintenance









Subgingival scaling


















Root canal treatment













Root canal retreatment Removal of foreign body



ED6 ED7 Special material (ED6, ED7) Thin, flexible, good conductivity to ultrasound



For removal of foreign body at bottom 1/3 of root canal.



For removal of 1For removal of foreign body foreign body and at top 1/3 of instrument fractures root canal (stainless steel) in root canal



GD30

ED15 For removal of temporary cement and instrument fracture (stainless steel) in root canal



Diamond-coated ball head Locate root canal orifice with security, accuracy and efficiencye

Cavity preparation







For the preparation of mesiocclusion cavity



For the preparation of distal cavity











ED12D

136

PULP TESTER





HEALTHCARE MANAGEMENT PTE LTD

Order name:C-PULSE



Feature:

- 1-Preset speed mode(high-mid-low speed)
- 2-Agentle,pulsed stimulus begins to increase at a rate specifically designed for patient comfort.3-If the patient indicates perception,simly release the button, the stimulus stops immediately,but the numbers remain frozen on the screen for about
- 3-minutes.
- 4-Specifically designed for patient comfort.
- 5-Convenient to operate.
- 6-Maximum stimulus--80.between 0-40,the patient feels ache and anesthesia, alive teeth nerve.
- 7-Between 40-80,with above-mentioned reaction,part of teeth nerve dead; 80,no above-mentioned reaction,dead teeth nerve.
- $8\mbox{-}\mbox{The device turns itself off in three minutes after operation (No wasted batteries).}$

This product has been applied layout-design patent, Counterfeit must be investigated; NO.ZL200730333925.9



IMPLANT SYSTEM



