



PROPHY BRUSH AND CUP



The product is used in conjunction with handpiece. It consists of metal sodium and a brush. The metal sodium is connected to the handpiece. The product is non-sterile and does not contain abrasive polishing materials. The brush head is made of nylon or bristle, and the stem is made of nickel-plated brass. The product is suitable for polishing, sanding and removing excess parts of restorations, or cleaning of implants. It is also used for descaling and polishing the tooth surface.

Prophy Brush(hard brush)



PB01A

PB01A,prophy brushes,cup-shaped,nylon,hard.

Item: PB01A

Specification: (Stick)22.8mm,
(Brush head)5.5mm

Packing: 100 pcs/box,100 boxes/ctn
Carton Size: 400X225X290mm

- Red
- Orange
- Yellow
- Green
- Blue
- Purple
- White



Prophy Brush(soft brush)



PB01B

PB01B,prophy brushes,cup-shaped,nylon,soft.

Item: PB01B

Specification: (Stick)22.8mm,
(Brush head)5.5mm

Packing: 100 pcs/box,100 boxes/ctn
Carton Size: 400X225X290mm

- Red
- Orange
- Yellow
- Green
- Blue
- Purple
- Transparent



Prophy Brush(bristles)



PB01C

PB01C,prophy brushes,cup-shaped,bristles.

Item: PB01C

Specification: (Stick)22.8mm,
(Brush head)5.5mm

Packing: 100 pcs/box,100 boxes/ctn
Carton Size: 400X225X290mm



Prophy Brush(hard brush)



PB02B

PB02B,prophy brushes,flat-shaped,nylon,soft.

Item: PB02B

Specification: (Stick)21.6mm,
(Brush head)5.5mm

Packing: 100 pcs/box,100 boxes/ctn
Carton Size: 400X225X290mm

- Red
- Orange
- Yellow
- Green
- Blue
- Purple
- Transparent



Prophy Brush(soft brush)



PB02C

PB02C,prophy brushes,flat-shaped,bristles.

Item: PB02C

Specification: (Stick)21.6mm,
(Brush head)5.5mm

Packing: 100 pcs/box,100 boxes/ctn
Carton Size: 400X225X290mm



AJH

HEALTHCARE MANAGEMENT PTE LTD

Prophy Brush(bristles)



■ PB03A,prophy brushes,tapered-shaped,nylon,hard.

Item: PB03A

Specification: (Stick)21.6mm,
(Brush head)7.5mm

Packing: 100 pcs/box,100 boxes/ctn

Carton Size: 400X225X290mm

Red
Orange
Yellow
Green
Blue
Purple
White



Prophy Brush(hard brush)



■ PB03B,prophy brushes,tapered-shaped,nylon,soft.

Item: PB03B

Specification: (Stick)21.6mm,
(Brush head)7.5mm

Packing: 100 pcs/box,100 boxes/ctn

Carton Size: 400X225X290mm

Red
Orange
Yellow
Green
Blue
Purple
Transparent



Prophy Brush(soft brush)



■ PB03C,prophy brushes,tapered-shaped,bristles.

Item: PB03C

Specification: (Stick)21.6mm,
(Brush head)7.5mm

Packing: 100 pcs/box,100 boxes/ctn

Carton Size: 400x225x290mm



Prophy Brush(bristles)



■ PC01A/PC01B,6-web prophy cup.

Item: PC01A(Hard)/PC01B(Soft)

Specification: (Stick)14.5mm,
(Silicone head)8mm

Packing: 100 pcs/box,100 boxes/ctn

Carton Size: 400X225X290mm

Red
Orange
Yellow
Green
Cyan
Purple



The product is used in conjunction with handpiece.It consists of metal sodium and a silicone cup.The metal sodium is connected to the handpiece.The product is non-sterile and does not contain abrasive polishing materials.The brush head is made of silicone,and the stem is made of nickel-plated brass.The product is suitable for polishing,sanding and removing excess parts of restorations,or cleaning of implants.It is also used for descaling and polishing the tooth surface.

Prophy Cup(hard/soft)



■ PC02A/PC02B,umbrella shape prophy cup.

Item: PC02A(Hard)/PC02B(Soft)

Specification: (Stick)14.5mm,
(Silicone head)8mm

Packing: 100 pcs/box,100 boxes/ctn

Carton Size: 400X225X290mm

Red
Orange
Yellow
Green
Cyan
Purple



Prophy Cup(hard/soft)



■ PC03A/PC03B,ribbed shape prophy cup.

Item: PC03A(Hard)/PC03B(Soft)

Specification: (Stick)14.5mm,
(Silicone head)8mm

Packing: 100 pcs/box,100 boxes/ctn

Carton Size: 400X225X290mm

Red
Orange
Yellow
Green
Cyan
Purple



Prophy Cup(hard/soft)



■ PC04A/PC04B,8-web prophy cup.

Item: PC04A(Hard)/PC04B(Soft)

Specification: (Stick)14.5mm,
(Silicone head)8mm

Packing: 100 pcs/box,100 boxes/ctn

Carton Size: 400X225X290mm

Red
Orange
Yellow
Green
Cyan
Purple



Prophy Cup(hard/soft)



■ PC01A/PC01B,6-web prophy cup.

Item: PC01A(Hard)/PC01B(Soft)

Specification: (Stick)14.5mm,
(Silicone head)8mm

Packing: 100 pcs/box,100 boxes/ctn

Carton Size: 400X225X290mm

Red
Orange
Yellow
Green
Cyan
Purple



■ PC02A/PC02B,umbrella shape prophy cup.

Item: PC02A(Hard)/PC02B(Soft)

Specification: (Stick)14.5mm,
(Silicone head)8mm

Packing: 100 pcs/box,100 boxes/ctn

Carton Size: 400X225X290mm

Red
Orange
Yellow
Green
Cyan
Purple



■ PC03A/PC03B,ribbed shape prophy cup.

Item: PC03A(Hard)/PC03B(Soft)

Specification: (Stick)14.5mm,
(Silicone head)8mm

Packing: 100 pcs/box,100 boxes/ctn

Carton Size: 400X225X290mm

Red
Orange
Yellow
Green
Cyan
Purple





■ PC04A/PC04B,8-web prophy cup.

Item: PC04A(Hard)/PC04B(Soft)
Specification: (Stick)14.5mm,
(Silicone head)8mm
Packing: 100 pcs/box,100 boxes/ctn
Carton Size: 400X225X290mm

- Red ■
- Orange ■
- Yellow ■
- Green ■
- Cyan ■
- Purple ■



COTTON ROLL DISPENSER

