

### Tray Base



TB01

Item: TB01

**Specification:** Square Head/Lower: L:100X78X28mm; M:90X69X28mm; S:81X60X28mm.  
Cusp Head/Upper: L:100X82X28mm; M:90X72X28mm; S:81X63X28mm.

**Packing:** 1 set/box,45 boxes/ctn

**Carton Size:** 440X320X250mm

**Weight:** 13 KGS

Blue 



### Teeth Whitening Tray



TT01

- Silicone medical EVA material thermoplastic dental tray, with high transparency, molded tray can maintain high transparency, especially suitable for use with cold light whitening instrument.

- Can be repeatedly molded, even if molded effect is not good, and it can be reshaped in hot water at 90 ° C.

- As a tray for tooth whitening and whitening glue, it can also be used as an anti-wear tooth tray.

Item: TT01

**Packing:** 2 pcs/unit,10 units/bag,50 bags/ctn

**Carton Size:** 460X350X260mm

**Weight:** 7.5 KGS



### Bite Tray



BT

Item: BT

**Packing:** 20 pcs/bag,50 bags/ctn

**Carton Size:** (BT01)520X350X375mm,(BT02)460X350X260mm, (BT03)430X430X290mm,(BT04)420X270X230mm, (BT05)280X260X230mm

**Weight:** (BT01)8.1 KGS,(BT02)5.5 KGS,(BT03)7.4 KGS,(BT04)3.6 KGS, (BT05)1.8 KGS



### Disposable Fluoride Foam Tray



DT01

- DT01, it is made of PE material, which is mainly used for washing teeth and protecting fluorine. The product is soft and flexible, which is comfortable for patients to bite. It can bite the upper and lower teeth at the same time, making the treatment process more convenient and quick.

Item: DT01

**Specification:** (L)193X65mm,(M)180X60mm,(S)160X55mm

**Packing:** (L)50 pcs/bag,50 bags/ctn,(M)50 pcs/bag,60 bags/ctn, (S)50 pcs/bag,80 bags/ctn

**Carton Size:** 710X360X655mm

**Weight:** (L)5.8 KGS,(M)5.5 KGS,(S)5.8 KGS

Green (S)   
Yellow (M)   
Blue (L) 



# MIXING BOWL AND SPATULA



### Hy-Drop

Patent Number: ZL 2018 3 0219158.7



■ Hy-drop consists of a medicine bottle, a dropper and a silicone pad. It is a small device for storing medicines for a short period of time. Mainly used to store oral medicine, such as camphor, formaldehyde cresol, 3% hydrogen peroxide, 2% sodium fluoride, etc.

**Item:** MW03  
**Specification:** 33X35X43mm  
**Packing:** 1 pc/box,144 boxes/ctn  
**Carton Size:** 440X440X250mm  
**Weight:** 4.8 KGS

- Light Yellow
- Light Green
- Light Blue
- Light Purple



### Disposable Dappen Dish

Patent Number: ZL 2015 2 0375947.0



■ Be put into using completely for one time,may save time and be guarded against cross-infection.Be used for etching agent and binder.

**Item:** DD01  
**Specification:** 23X32mm  
**Packing:** 100 pcs/box,20 boxes/ctn  
**Carton Size:** 490X170X190mm  
**Weight:** 4.5 KGS

- Yellow
- Pink
- Blue
- White



### Glass Dappen Dish



■ Assorted coloured glass dappen dishes,suitable to mix small dose.

**Item:** GD01  
**Specification:** 33X33mm  
**Packing:** 1 pc/box,12 boxes/unit,18 units/ctn  
**Carton Size:** 360X170X280mm  
**Weight:** 11 KGS

- Amber
- Green
- Blue
- White



### Prophy Paste Ring

Patent Number: ZL 2015 2 0376134.3



■ Put it on our fingers.Convenient for dentists to take medicine.

**Item:** DD030A  
**Packing:** 20 pcs/bag,100 bags/ctn  
**Carton Size:** 430X430X290mm  
**Weight:** 7.3 KGS

- Light Green
- Light Blue
- Light Purple



### Prophy Paste Ring

Patent Number: ZL 2015 2 0376134.3



■ Put it on our fingers.Convenient for dentists to take medicine.

**Item:** DD030B  
**Packing:** 20 pcs/bag,100 bags/ctn  
**Carton Size:** 430X430X290mm  
**Weight:** 4.5 KGS

- Light Green
- Light Blue
- Light Purple



### Disposable Dappen Dish

Patent Number: ZL 2015 2 0376134.3



■ New disposable dishes can be stucked on the hands for using more conveniently.

**Item:** DD02  
**Specification:** 32X13mm  
**Packing:** 100 pcs/bag,75 bags/ctn  
**Carton Size:** 520X520X420mm  
**Weight:** 7.4 KGS

- White



### Mixing Well



**Item:** MW01/MW02  
**Specification:** (MW01)51X51X10mm,(MW02)51X28X10mm  
**Packing:** (MW01)200 pcs/box,24 boxes/ctn  
(MW02)200 pcs/box,45 boxes/ctn  
**Carton Size:** (MW01)400X400X370mm,  
(MW02)440X340X320mm  
**Weight:** (MW01)12.3 KGS,(MW02)11.5 KGS

- White



### Mixing Rod



**Item:** MR01  
**Packing:** 100 pcs/bag,70 bags/ctn  
**Carton Size:** 290X260X230mm  
**Weight:** 12.8 KGS

- White



### Mixing Bowl



■ Flexible plastic products,can be used repeatedly after cleaning.Four size

Item	Capacity	Specification	Packing	Carton Size	Weight	Color
MB01	690ml	130X96mm	100 pcs/ctn	560X270X270mm	12.6 KGS	Blue,Clear Green,Clear Red
MB02	280ml	108X59mm	180 pcs/ctn	430X330X290mm	14.9 KGS	Blue,Clear Green,Clear Red
MB03	150ml	84X50mm	250 pcs/ctn	420X270X240mm	9.2 KGS	Blue,Clear Green,Clear Red
MB04	500ml	130X86mm	100 pcs/ctn	560X270X270mm	12 KGS	Green



### Silicone Small Bowl



■ Flexible silicone rubber products,five different models suitable for different work,washable for repeated use.Sucker cup bottom, with easy to put the bowls firmly on the smooth surface such as glass.

Item	Capacity	Specification	Packing	Carton Size	Weight
SB01 XL	100ml	70X55mm	1 pc/box,144 boxes/ctn	440X440X250mm	6 KGS
SB01 L	30ml	50X40mm	12 pcs/box,25 boxes/ctn	480X290X320mm	6 KGS
SB01 M	10ml	36X38mm	12 pcs/box,56 boxes/ctn	390X390X370mm	9 KGS
SB01 S	8ml	31X23mm	12 pcs/box,120 boxes/ctn	430X430X290mm	8.1 KGS
SB01 SS	5ml	25X25mm	12 pcs/box,144 boxes/ctn	430X430X290mm	7.5 KGS

### Spatula



Dentists can choose different types of spatula according to different medicines and mixing amounts; the deformed design of the handle and the texture of the matte increase friction and are convenient to use.

**Item:** SP01  
**Specification:** 195X32mm  
**Packing:** 50 pcs/bag,10 bags/ctn  
**Carton Size:** 420X270X240mm  
**Weight:** 10 KGS

- Pink
- Blue
- White



### Spatula



**Item:** SP02  
**Specification:** 185X25mm  
**Packing:** 50 pcs/bag,20 bags/ctn  
**Carton Size:** 420X270X240mm  
**Weight:** 9 KGS

- Pink
- Blue
- White



### Spatula



**Item:** SP03  
**Specification:** 185X33mm  
**Packing:** 25 pcs/bag,16 bags/ctn  
**Carton Size:** 420X270X240mm  
**Weight:** 10 KGS

- Pink
- Blue
- White



### Disposable Mouth Mirror



**Item:** MM01  
**Specification:** 160X20mm  
**Packing:** 1 pc/pouch,100 pouches/bag,20 bags/ctn  
**Carton Size:** 430X430X290mm  
**Weight:** 14 KGS

- Yellow
- Pink
- Blue
- White



### Disposable Mouth Mirror



**Item:** MM04  
**Specification:** 170X20mm  
**Packing:** 1 pc/pouch,100 pouches/bag,20 bags/ctn  
**Carton Size:** 430X430X290mm  
**Weight:** 15 KGS

- White



## BASE PLATE WAX

