# RETRACTOR








#### **Dental Retractor**

Patent Number: ZL 2019 3 0018928.6



Item: CR10M Packing:2 pcs/box,100boxes/ctn Carton Size:460X350X260mm Weight:5KGS



#### Advantage:







#### **Dental Retractor**



■ The parts of the retractor in contact with the two sides of the mouth of the patient use soft rubber to make protective design, which makes the patient more comfortable to use.

Item: CR10S

Packing: 2 pcs/box,100boxes/ctn



#### **Mouth Prop** -

Patent Number: ZL 2015 3 0162734.5



Easily avoiding patients fatigue by using mouth prop, make patients enjoy a relaxing treatment.

Packing: 2 pcs/bag,12 bags/box,32 boxes/ctn

Carton Size: 440X440X250mm

Weight: (L)13 KGS,(M)9.5 KGS,(S)7.7 KGS

Item No	Specification	
MP01L	39.5X39.5X20.4mm	
MP01M	34.5X36.5X17mm	
MP01S	33X28X17mm	





#### **Mouth Prop & Tongue Guard**

Patent Number: ZL 2014 2 0507326.9



**Cheek Retractor** 

28

Tongue guard block tongue moving so that no injury would happen during the whole treatment.



Weight: (L)11.5 KGS,(M)10.5 KGS,(S)9.5 KGS

It consists of C-shaped connecting rod and retractor on both sides. After special processing, the surface is white translucent, smooth, high strength and good toughness.

Item: CR09

Remark: Frosted Surface







#### **Cheek Retractor**



- This product consists of a semilunar retractor and a handle. During oral treatment, the medical staff can place the retractor on one side of the mouth corner of the patient and gently pull the open lip to provide a clear surgical field vision.
- Type: T type Color: Clear

color. Cicai				
Item No	Specification	Carton Size	Packing	Weight
CR01L	120X53.5X18mm	430X330X320mm	2 pcs/box,	6.3 KGS
CR01S	114X42X17mm	360X360X280mm	200 boxes/ctn	5.8 KGS



#### **Cheek Retractor**



- This product is composed of C-shaped connecting rod and retractor on both sides. The material has good strength and toughness. During oral treatment, put it at the corner of both sides of the patient, open the mouth and make it keep a certain degree of opening.
- Type: C type Color: Clear

Item No	Specification	Carton Size	Packing	Weight
CR02L	132X90X17mm	500X360X440mm	2 pcs/box, 200 boxes/ctn	8.2 KGS
CR02M	118X85X17mm			8 KGS
CR02S	97X70X15mm	430X335X290mm		5.5 KGS



#### **Cheek Retractor**



- This product is M type design.In oral treatment, medical staff can place it in both sides of the mouth position separately to open the mouth, block the tongue, and provide a clear surgical field vision.
- Type: M type Color: Clear

т.					
	Item No	Specification	Carton Size	Packing	Weight
	CR03L	126X57X51mm		2 pcs/box, 147 boxes/ctn	6 KGS
	CR03M	106X48X54mm	420X420X470mm	2 pcs/box,	8.6 KGS
	CR03S	100X44X40mm		246 boxes/ctn	7.5 KGS



#### **Cheek Retractor**



- This product is a circular and hollow design.In oral treatment, medical staff can place it on both sides of the patient's mouth to open the oral vestibule and provide a clear surgical vision.
- Type: O type

Color: Clear				
Item No	Specification	Carton Size	Packing	Weight
CR04L	84X80x22mm	440X440X250mm	10 pcs/bag,	6.5 KGS
CR04S	70X68X35mm	440,440,7230111111	100 bags/ctn	6.2 KGS



#### **Cheek Retractor**



- This product is composed of double-headed semilunar retractor and handle are formed.
- Type: Double head Color: Clear

li			
		CR05	
	-		

Item No	Specification	Carton Size	Packing	Weight
CR05L	167X55X24mm	400X400X370mm	2 pcs/box, 125 boxes/ctn	7.7 KGS
CR05S	145X45X20mm	400X400X370mm	2 pcs/box, 120 boxes/ctn	5.9 KGS



www.ajjhealthcare.com 29

### **Aligner Chewies**

Patent Number: ZL 2018 2 1528611.3



Mint flaver Juicy peach flaver

Grape flaver Green apple flaver

Lemon flaver Watermelon flaver

- Help aligner tray fit teeth closer to make sure reach preset orthodontics force.
- Bite for 5-10 minutes to avoid braces face.
- Gum bite and handle are designed separately,convenience and hygiene.
- Triangular prism shape uses wedge-force principle to fit aligner tray for adjacent anterior or premolar region.
- Orlinder is used to fit aligner tray in molar region with greater occlusal force, and it can also be used to exercise bite muscle.

Elliptic cylinder is mainly used to fit aligner tray in occlusal surface of posterior teeth.





tem: MP03

	Вох	
Packing 1	3 pcs/box,400 boxes/ctn	2 pcs/pouch,12 pouches/bag, 30 bags/ctn
Packing 2	2 pcs/bax, 400 baxes/ctn	3 pcs/pouch,12 pouches/bag, 30 bags/ctn

#### **Retainer Retrievers**

Patent Number: ZL 2019 3 0018952.X



- Retainer retrievers help patients remove appliances cleanly and efficiently.
- Helps patients remove appliances from the buccal side cleanly and efficiently.
- The ring handle is easy to hold by hand when operating.
- This handy tool features specially designed hooks which help to remove clear retainers or aligners with ease.

Item: MP0

Packing: 1 pc/unit,12 units/bag,30 bags/ctn Carton Size: 290X260X230mm

Weight: 2.1 KGS







#### **Dry Field System Tubing**

Patent Number: ZL 2016 2 1097848.1



White (L) Green (S)

- The Dry Field System completely eliminates saliva and fully exposes both arches making it ideal for all bonding applications, The system is easy to assemble and be inserted and used by a single operator.
- Set salivary absorption, retractor, tongue guard in one.

Item: CR08

Packing: (L)1 set/box,48 boxes/ctn,(S)1 set/box,50 boxes/ctn
Carton Size: (L)440X440X250mm,(S)460X350X260mm
Weight: (L)5.8 KGS,(S)5 KGS



### **Dry Field System Component**



Item: CR0802

Packing: 1 set/bo

Remark: Components consist tongue guard, connection arm, position block, 2 pcs of silicone tubing and Y-connector.

## **DISPOSABLE SALIVA**

## **EJECTORS**

