



BUR HOLDER BOX



Bur Holder Box



- BB01, it is made of plastic material, can place high speed burs and low speed burs.
- Bottom part can be autoclavable at 121°C.

Item: BB01
 Specification: 89X66mm
 Packing: 1 pc/box, 60 boxes/ctn
 Carton Size: 480X290X320mm
 Weight: 6 KGS
 Remarks: 60 slots

Blue
 White



Bur Holder Box



- BB02, it is made of plastic, specially designed for high-speed burs.
- Bottom part can be autoclavable at 121°C.

Item: BB02
 Specification: 74X60X36mm
 Packing: 1 pc/box, 100 boxes/ctn
 Carton Size: 430x270x240mm
 Weight: 4.5 KGS
 Remarks: 24 slots

Blue
 White



Bur Holder Box



- BB03, it is made of plastic material, can place high speed burs and low speed burs.
- Bottom part can be autoclavable at 121°C.

Item: BB03
 Specification: 89X91mm
 Packing: 1 pc/box, 60 boxes/ctn
 Carton Size: 480X390X320mm
 Weight: 7 KGS
 Remarks: 60 slots

Blue
 White



Bur Holder Box



Item: BB09
 Specification: 120X53mm
 Packing: 1 pc/box, 54 boxes/ctn
 Carton Size: 400X400X370mm
 Weight: 7 KGS
 Remarks: 166 slots

Blue



Bur Holder Cleanser



- BC01, it is made of autoclavable plastic materials, suitable for the placement of tools such as burs, root canals, and enlarged needles before oral treatment; with transparent materials, dentists can easily find the tools they need.

Item: BC01
 Specification: 83X90mm
 Packing: 1 pc/box, 48 boxes/ctn
 Carton Size: 390X390X270mm
 Weight: 7.5 KGS
 Remarks: 60 slots

Blue
 Clear



When the dentist is performing root canal treatment, the closedness of the apex often determines the success or failure of the final treatment. Sometimes the root canal is well prepared, but it is easy to overfill or underfill when filling. The main reason is that the size of gutta percha point is not enough, either too large or too small, and some dentists ignore the use of gutta-percha measuring ruler.


Gutta-Percha Measuring Ruler



- BF0517, promisee Dental developed the gutta-percha measuring ruler which made of Zirconia.
- Length range is 30mm, hole diameter is from 20-80, total 11 tapered holes, which can meet the requirements of gutta percha point accuracy Length and diameter.
- Suitable for accurate measurement of the diameter of the apex of the main apex of the gutta percha point, and can also be used as a tool for inspecting the wear of the root canal instrument.
- Ultra-high accuracy, the accuracy does not change under different environments.

Advantage:

- Health and environmental protection. Ceramic material itself has chemical stability, antioxidant, restlessness, acid and corrosion resistant.
- Stylish appearance, mellow color, pure jade texture, bring noble enjoyment.
- Smooth Surface, easy to clean.
- High hardness, high wear resistance.
- High temperature and pressure resistant. It is applicable to accurate tip diameter measurement and also can be used as a tool to test tip wear of root canal instruments.



Item: BF0517
Specification: 90X20X3.5mm
Packing: 1 pc/box



Bur Holder Frame



Item: BF01
Specification: 65X38X12mm
Packing: 20 pcs/bag, 36 bags/ctn
Carton Size: 440X340X300mm
Weight: 8.6 KGS
Remarks: 16 slots



Root Canal Frame



BF02, the product is made of PA material and has a hollow structure, which is convenient for taking instruments. It is mainly used for measuring the length of root canal files during dental treatment.

Item: BF02
Specification: 75X52.5X10mm
Packing: 3 pcs/box, 48 boxes/ctn
Carton Size: 265X170X180mm
Weight: 4.5 KGS



Root Canal Measuring Table




BF03A/BF03B, the product is made of PA material, which is autoclavable 134 ° C, and has a hollow structure design, which is convenient for handling equipment.

Item: BF03A/BF03B
Specification: (BF03A)89x20.6x32mm,
(BF03B)89X26X32mm
Packing: 1 pc/box, 250 boxes/ctn
Carton Size: (BF03A)320X320X250mm,
(BF03B)430x270x240mm
Weight: (BF03A)9 KGS, (BF03B)9.6 KGS




Microdot



BF04, the product is made of PA material, which is autoclavable 134 °, easy to do operation.

Item: BF04
Specification: 90X13X2.8mm
Packing: 3 pcs/package, 10 packages/box, 14 boxes/ctn
Carton Size: 240X155X107mm
Weight: 2 KGS




Root Canal Measuring Ruler



BF05, the product is made of PA material, which is autoclavable 134 °, easy to do operation.

Item: BF05
Specification: 42X23.3X23.9mm
Packing: 2 pcs/package, 24 packages/box, 6 boxes/ctn
Carton Size: 240X155X107mm
Weight: 1.5 KGS



Grinding Carriage

Patent Number: ZL 2015 2 0376130.



Item: BF06
Specification: 128X82X23mm
Packing: 1 pc/box, 200 boxes/ctn
Carton Size: 500X360X440mm
Weight: 14 KGS

