> K-Flex Files













Material: Stainless Steel Package: 6 pcs per pack (assorted/single size)

Size: ISO sizes Sterilization: ≤126°C Length: 21,25,28,31 mm

Features and Benefits:

More flexible than conventional endodontic files of comparable size

• Prepare root canals with minimal potential for ledging and perforation.



> Barbed Broaches













Material: Stainless Steel Package: 10 pcs per blister (assorted/single size)

Size: ISO sizes Sterilization: ≤126°C Length: 21,25,28 mm






Other **Instruments**



> U Files




> Paste Carriers



> Peeso Reamers



> Gates Drills





[→] HD-P



Material: Aluminium alloy handle; Ni-Ti working part or stainless steel working part.

Package: 1 piece per box

Advantages:

- ISO sizing and tapering, ISO color coded
- Conical tip and flat tip, ideal for both lateral and vertical compaction of gutta-percha
- Smooth surface, not stick the fillings when pressing
- Aluminum alloy handle, light weight and greater oxidation resistance



06

Endodontic Material

