

**FY0247  
Suction Liner**

**Product Characteristic:**

Simplicity  
Safety  
Economical

**Indication for Use:**

This product is used for collecting the biohazardous fluids which produced in ICU and operation room.

**Product Features:**

1. Double safety design to ensure no liquid is suctioned back to the tubing and suction unit, effectively prevent cross infection and protect the medical staff safe.
2. It's made of environment-friendly polyethylene. It's safe and non-toxic.
3. Easy to use, can be a single use, also can use with the stander.
4. Available in some specifications: with or without shutoff filter, with or without solidifier, 1000ml, 1500ml and 2000ml.


**Product Advantages:**

1. Soft and transparent (Convenient for observation liquid waste)
2. Double safety design (To prevent fluid into the suction systems by mistake)
3. It can be connected with relevant suction device or vacuum device, soft bag is seal and reliable inside and can be withstand the pressure of 60kg when fully filled with liquid (Prevent infections liquid overflow caused by medical personnel).



**Specification**

Item	Description
Stand	Four tandem stand 4 canisters ( a set )
	Twin tandem stand 2 canisters ( a set )
Hanger	Two tandem hanger stand
	Single hanger stand
	Single canister



**HYPODERMIC**





**FY0601**  
**Disposable Hypodermic Syringe**

1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml  
Three-part syringes/Two part syringes ( luer slip,  
special design to prevent plunger slipping out )  
Luer slip/luer lock  
Latex/latex free piston  
Sterilization: EO  
PE/Blister individual package, 100pcs/Inner Box  
W/needle, W/O needle



**FY0602**  
**Auto Self-destruction Syringe**

0.05mm, 0.1ml, 0.5ml, 1ml  
Fixed needle 27G-30G  
Sterilization: EO  
PE/Blister individual package, 200pcs/Inner Box



**FY0603**  
**Tubercle Bacillus Syringe**

1ml  
luer slip  
Latex/latex free piston  
PE/Blister individual package  
Sterilization: EO



**FY0604**  
**Insulin Syringe**

0.3ml, 0.5ml, 1ml  
Fixed needle 27G-31G  
PE/Blister individual package  
Sterilization: EO

**FY0617**  
**Disposable Destruction Syringe**

1ml, 2/3ml, 5ml, 10ml, 20ml  
After injection , when you pull the plunger,  
the needle will be back into the barrel.  
It's able to prevent it from hurting people  
by covering hypodermic needle with barrel.  
Luer lock  
Sterilization: EO



**FY0618**  
**Safety Clip Syringe**

1ml, 2/3ml, 5ml, 10ml, 20ml  
Safety syringe with protective cap  
Luer lock  
After injection, Push the clip over the needle and lockit  
onto the clip base  
Sterilization: EO



**FY0401**  
**Irrigation Syringe**

50ml, 60ml, 80ml  
Flat/Bulb/Ring  
PE/Blister individual package  
Sterilization: EO



**FY0501**  
**Disposable Infusion Set**

For infusion bags and bottles  
20/60 drops = 1±0.1ml  
Small Chamber  
Latex tube  
With Filter  
Different types of components and parts can be  
changed as customer's requirements  
For gravity use  
EO sterile  
PE/Blister individual package







### FY0502 Disposable Infusion Set

For infusion bags and bottles  
20/60 drops =  $1 \pm 0.1$ ml  
Small Chamber  
With Y site  
Different types of components and parts can be changed as customer's requirements  
For gravity use  
EO gas sterile  
PE/Blister individual package



### FY0503 Disposable Infusion Set

For infusion bags and bottles  
20/60 drops =  $1 \pm 0.1$ ml  
Big Chamber  
With Filter  
Different types of components and parts can be changed as customer's requirements  
For gravity use  
EO gas sterile  
PE/Blister individual package



### FY0504 Burette Infusion Set

Burette Capacity: 100ml, 150ml  
For infusion bags and bottles  
60 drops =  $1 \pm 0.1$ ml  
With Filter  
With Y Site  
Different types of components and parts can be changed as customer's requirements  
For gravity use  
EO gas sterile  
PE/Blister individual package



### FY0505 Disposable Blood Transfusion Set

For blood bags  
20 drops =  $1 \pm 0.1$ ml  
Different types of components and parts can be changed as customer's requirements  
For gravity use  
EO gas sterile  
PE/Blister individual package


### FY0605 Disposable Hypodermic Needle

Size: 18-30G  
Length of needle: 1/2", 5/8", 3/4", 1", etc.  
Blister individual package  
Sterilization: EO



### FY0606 Scalp Vein Set

Size: 18-30G  
Length of needle: 1/2", 5/8", 3/4", 1", etc.  
Tube length: 30cm  
Luer slip/luer lock  
PE/Blister individual package  
Sterilization: EO



### FY0616 Insulin Pen Needle

Ultra fine cannula for maximum patient comfort, patient point with triple beveling for easy and pain decreases  
The needle is fit for Humapen ( EliLily, USA )  
NovopenIII ( Novo Nordisk, Denmark )  
Sterile, non-pyrogenic  
Sterilization: EO



### FY0701 Disposable Blood Bag

For the collection of a blood volume: 450ml, etc.  
Is filled with 63ml of CPDA-1 ANTICOAGULANT SOLUTION  
Each 100ml CPDA-1 (U.S.P.) contains:  
Citric Acid (anhydrous)-----0.299g  
Sodium Citrate (dihydrate)-----2.63g  
Monobasic Sodium Phosphate (monohydrate)-----0.222g  
Dextrose (monohydrate)-----3.19g  
Adenine-----0.0275g  
Water for injection, a sufficient quantity to make-----100ml  
Needle: 16 gauge  
Steam sterilization  
Sterile, non-toxic, non-pyrogenic fluid path





**FY0508  
Heparin Cap**

Material: PC  
Sterilization: EO  
Individual Hard Blister Packing, 100pcs/Box, 5000pcs/ctn



**FY0510  
Medical Disposable  
Consumable Spiral Stoppers**

Material: ABS Sterilization : EO  
Color: Red, White, Blue, etc.  
Individual Hard Blister Packing, 100pcs/Box

**FY0620 Medical Connector**



**FY062001  
Needle Free Connector**

**Features:**  
1. Transparent housing material: Polycarbonate or Copolyester.  
2. Metal-free & MRI-compatible.  
3. Latex-free.  
4. Multiple insertion at least 100 times/day.  
5. Priming Volume: 0.09ml.  
6. Perfect flowing rate: 350 ml/min under one meter water pressure



**FY062002  
Y-port Valves**

Material: Polycarbonate or Copolyester  
Size: ID 4.01mm, 3.6mm



**FY062003  
Tube End Valves**

Material: Polycarbonate or Copolyester  
Size: ID 2.0mm, 2.5mm, 2.8mm, 3.7mm, 4.1mm



**FY062004  
Positive Pressure Needle Free Connector**

**Features:**  
1. Patent design of positive pressure avoids back flow of blood during removal of syringe. This will help to prevent blood clot at the catheter tip inside the blood vessel.  
2. The housing is impregnated with PC Ag+ which helps to lower infection.  
3. Protruding design of the stalk at the injection port follows infection control protocol.  
4. High quality spring ensures multiple insertion at the injection site without leaking.  
5. Two sealing rings at the upper and lower end of the stalk keeps the connector out air, fluid and outside substance.  
6. The straight flow of fluid path creates less turbulence and complies with proper transfusion protocol.  
7. The holding capacity of remaining fluid is less than 0.06ml.



**FY0511  
Micro-flow Regulator**

Flow range: 5ml/h-250ml/h  
Materials:  
Regulator: ABS  
Tube: PVC  
Other accessories: PC & ABS  
Extension Tube Length: 45cm is common, others are also available.  
Sterilization: EO  
Packing Details: Poly-bag per pcs, 200pcs/inner box, 20boxes/carton



**Manifolds**

**FY0507  
Three Way Stopcock**

Parts: Body(PC), Plas(PE), Screw Cap(ABS)  
Individual Hard Blister Packing, 50pcs/Box, 1000pcs/ctn  
Sterilization: EO





## FY0509 Extension Tube



### FY050901 Disposable Extension Tube

The length of extension tube can be changed as customer's requirements  
 Different types of components and parts can be changed as customer's requirements  
 Sterilization: EO  
 Packing: Individual blister Packing

#### Tube Length:

10-25cm, 50pcs/ctn  
 30-50cm, 50pcs/ctn  
 80cm, 50pcs/ctn  
 100-120cm, 40pcs/ctn



### LB0630 Central Venous Catheter

Advanced automatic drilling technology and equipment, a more accurate location and the pore diameter of the pipe, surface smooth and without cutting residues and injury.

Special custom polyurethane fine materials, delicate and smooth tube, welding head end, ultra smooth surface treatment and X line to show the position, make a catheter inserted into the blood vessels in more easy and safe.



Screw Extension Tube

Mini-Quantity Infusion Extension Tube (OD:2.0mm)

T-Extension Tube

## FY050902 Needleless Adapters for Infusion Set



1. This device is designed to meet the needs of general I.V. Therapy, Anaesthesia Cardiovascular,
2. ICU&CCU, Recovery&Oncology.
3. Safe for infusion IV medication including most high risk anti-neoplastics, or aspirating blood without risk of hemolysis.


#### Catheter Diameter:

1. One-Way Extension Tube with One Normal Needle Free Connector
2. One-Way Extension Tube with One Positive Pressure Needle Free Connector
3. Two-Way Extension Tube with Two Normal Needle Free Connectors
4. Two-Way Extension Tube with Two Positive Pressure Needle Free Connectors
5. Three-Way Extension Tube with Three Normal Needle Free Connectors
6. Three-Way Extension Tube with Three Positive Pressure Needle Free Connectors
7. Four-Way Extension Tube with Four Normal Needle Free Connectors
8. Four-Way Extension Tube with Four Positive Pressure Needle Free Connectors



### FY3007 Peripherally Inserted Central Catheter Kit

1. The connector is round and smooth with PASV valve.  
 A) Valve opens inward for sampling under negative pressure.  
 B) Valve opens outward for infusion under positive pressure.  
 C) Valve remains closed to prevent air embolism, blood reflux and clot under normal pressure
2. The catheter is made from premium silicone, which has excellent biocompatibility.
3. Adjustable length of catheter.
4. Detectable under X ray, the end tip enhances radiopacity.
5. The style of catheter ensures placement more convenient and smooth.
6. Variety of kit configuration, initial guide wire and positive, pressure needless connector is available for choice.



### FY0631 Anesthetic Products Kit

#### Standard Kits include:

1. Central Venous Catheter
2. Guide-wire
3. Vessel Dilator
4. Clamp
5. Fastener: Catheter Clamp
6. Introducer Needle
7. Introducer Syringe
8. Injection Needle
9. Injection Cap

#### Optional Compound Kits include:

1. Central Venous Catheter Standard Kit Accessories
2. 5ml Syringe
3. Surgical Gloves
4. Surgical Pledget
5. Surgery Sheet
6. Surgery Towel
7. Sterile Brush
8. Gauze Pad
9. Suture of Needle
10. Wound Dressing
11. Scalpel





### FY0632 Anesthetic Products Kit

#### Standard Kits include:

1. Drainage Catheter
2. Guide-wire
3. Vessel Dilator
4. Introducer Needle
5. Introducer Syringe
6. Injection Needle
7. Injection Cap
8. Connecting tube

#### Optional Compound

##### Kits include:

1. Drainage Catheter
2. 5ml Syringe
3. Surgical Gloves
4. Surgical Pledget
5. Surgery Sheet
6. Surgery Towel
7. Sterile Brush
8. Gauze Pad
9. Suture of Needle
10. Wound Dressing
11. Scalpel



### FY0635 Disposable Pressure Transducer

Disposable Pressure Transducer is made of medical PC and PVC material.

#### Operating Parameters:

Operating Pressure Range:  $-30\text{mmHg} \sim +300\text{mmHg}$

Operating Temperature Range:  $15^{\circ}\text{C} \sim 40^{\circ}\text{C}$

Storage Range:  $-25^{\circ}\text{C} \sim +70^{\circ}\text{C}$

Atmosphere Pressure Range:  $56.7\text{kPa} \sim 113\text{kPa}$

Zero Output:  $\pm 25\text{mmHg}$

Sensitivity:  $5 \mu\text{V/V/mmHg} \pm 2\%$

Input Impedance:  $1200 \Omega \sim 3200 \Omega$

Output Impedance:  $285 \Omega \sim 315 \Omega$

Zero Thermal Drift:  $\pm 0.3\text{mmHg}/^{\circ}\text{C}$

Zero Drift:  $\leq 2\text{mmHg}$

Natural Frequency:  $40\text{Hz}$  nominal for a standard kit  $> 200\text{Hz}$  for transducer alone

Sensitivity Thermal Drift:  $\pm 0.1\%/^{\circ}\text{C}$

Leak Current:  $< 2 \mu\text{A}$  at  $230\text{VRMS}$ ,  $50\text{Hz}$

Note: rated voltage  $6.00\text{VDC}$ ; ambient temperature  $25^{\circ}\text{C}$ .

#### Specification:

130cm pressure tubing with M/F connector	1pc
30cm pressure tubing with M/F connector	1pc
Disposable pressure transducer	1pc
Flush device 3 cc/h (30cc/h)	1pc
Male non-vented Cap	2pcs
150cm single IV set	1pc
3-way stopcock	2pcs



### FY0633 Blood Purification Products Kit

#### Standard Kits include:

1. Hemodialysis Catheter
2. Vessel Dilator
3. Introducer Needle
4. Syringe
5. Guide-wire
6. Adhesive Wound Dressings
7. Heparin Caps
8. Scalpel
9. Needle with Suture

#### Optional Compound

##### Kits include:

1. Hemodialysis Catheter
2. 5ml Syringe
3. Surgical Gloves
4. Surgical Pledget
5. Surgery Sheet
6. Surgery Towel
7. Sterile Brush
8. Gauze Pad



### FY0634 Blood Purification Products Kit

#### Characteristics:

1. The Catheter is made from premium silicone, which has excellent biocompatibility.
2. The dacron cuff can protect immigration of catheter and reduce the chance of CRBI (Catheter Related Bloodstream Infections).
3. The silicone material makes catheter more floppy, hence it can minimize the trauma to the vein and reduce the incident of clot.

#### Standard Kits include:

1. Silicone Catheter
2. Vessel Dilator
3. Valved Peelable Introducer
4. Introducer Needle
5. Guide-wire
6. Tunneler (Tocar)
7. Adhesive Wound Dressings
8. Heparin Caps
9. Gauze Sponges
10. Scalpel
11. Syringe
12. Needle with Suture



### FY0615 Epidural Set

#### Features:

Epidural needle: 16G/18\*80/90mm

Epidural catheter: 1.0/0.8mm, clear marked, length longer than 850mm with 3 holes or more holes, closed tip

Epidural filter: 0.22um, Luer lock

LOR syringe: 7ml, 10ml available

Catheter adaptor: TBA, secure the catheter completely

Single use to avoid the cross infection

Use in the clinic of epidural nerve or subarachnoid anesthesia





**FY050601  
Dental Needle**

Size: 27g x 38mm 30g x 38mm , etc.  
Sterilization : EO  
Packing: 100pcs/Box, 5000pcs/ctn




**FY050602  
Dental Needle**

Size: 18G-25G  
Packing: 100pcs/bag, 600bags/ctn



**FY0611 Spinal Needle**

Size: 18g~28g x 90mm  
Sterilization: EO  
Packing: 1pc/Bag, 50pcs/Box, 1000pcs/ctn



**FY0612 A.V. Needle**

Size: 17g, 18g, etc.  
Sterilization: EO  
Packing: 1pc/Bag, 10pcs/Box, 500pcs/ctn



**FY0607**



**FY0608**

**I.V. Cannula**

Normal/Safety  
Pen type/With wings/With injection port  
Size: 14G-26G  
PTFE/FEP/PU flex catheter  
Sterilization: EO  
Packing: Hard Blister Packing



**FY0609**

All the popular sizes are available in both a carbon or stainless steel sterile Blade shapes #9 to 17 are used in histology units and are all compatible with the associated #3 handles (Blade shape 11 can be also used for gardening with a plastic or metal handle). Blade shapes #18 to 36, #60 and #70 are used for major invasive procedures and are all compatible with the associated #4 handles  
Sterilization: Gamma  
Packing: 1pcs/Sealed Aluminum Foil Bag, 100pcs/box, 5000pcs/ctn



**FY0901  
Disposable Surgical Blade**

**FY2209 Disposable electrosurgical pencil ( Tips )**



63.5mm 63.5mm 63.5mm 63.5mm 63.5mm 63.5mm

**FY2209 Disposable electrosurgical pencil**

**FY220901**



Electrode Shape: Blade  
Handle Length: 158mm  
Hand control  
Cable Length: 3m  
Exposed electrode length: 40mm

**FY220902**



Electrode Shape: Loop  
Handle Length: 158mm  
Hand control  
Cable Length: 3m  
Electrode Exposed Length: 40mm

**FY220903**



Electrode Shape: Blade  
Handle Length: 158mm  
Foot control  
Cable Length: 3m  
Exposed electrode length: 40mm



### FY0902 Stitch Cutters

All sizes are available in both a carbon or stainless steel sterile  
Designed with a rounded bull nose tip, the Index family of Stitch  
Cutter Blades are ideal for  
cutting the wide range of modern sutures available today whilst  
minimising patient discomfort  
The range consists of a short standard stitch cutter with a larger  
grip. The long handle enables  
a more aseptic technique allowing good visibility whilst working  
further from the wound site  
Size: 65mm, 110mm  
Sterilization: Gamma  
Packing: 1pcs/Sealed Aluminum Foil Bag, 100pcs/box, 5000pcs/ctn



### FY0903 Surgical Handle

Index stainless steel metal handles are supplied  
non-sterile and individually wrapped  
Sterilization: Gamma  
Packing: 1pcs/Sealed Aluminum Foil Bag, 10pcs/box, 500pcs/ctn



### FY0904 Surgical Scalpels

Standard sheath guarded and blade protected sterile  
disposable scalpels offer convenience in a number of  
user environments: on the ward, in emergency or trauma  
units or even out in the field with paramedics  
Material: Carbon Steel or Stainless Steel  
Sterilization: Gamma  
Packing: 10pcs/Box, 500pcs/ctn



### FY0905 Steel Lancet

Low cost: simple, flat, steel lancet  
Sample volume: excellent blood flow  
High quality: hardened surgical steel  
Minimum discomfort: modern asymmetric design  
Sterile: sterilized by gamma radiation  
Size: 40 × 5.5 × 1.5mm  
Packing: 5pcs of packaged, 200pcs/Box, 20000pcs/ctn



### FY0906 Twist Lancets

Conventional twist-off cap design  
Sterilized by gamma-radiation  
Smooth tri-bevel needle tip for  
comfortable sampling experience  
Versatile with most lancing devices  
Size: 21G, 23G, 26G, 28G, 30G  
Packing: 100pcs/Box, 20000pcs/ctn



### FY0907 Twist Lancets (Exclusive to Roche Diagnostics)

5 twist lancets (28G) are used with  
Roche Accu-Chek Softclix Lancing Device and the  
needle tip is siliconized for easy and smooth motion  
during lancing  
Size: 28G, 30G  
Packing: 100pcs/Box, 20000pcs/ctn



### FY0910 Safety Lancet

Ergonomic design offers a more secure hold when  
to activate the device  
Tri-bevel needle tip ensures the sharpness and accuracy  
The spring buffer the puncture making the procedure  
a more comfortable experience  
The complete range of needle sizes from 18G to 32G  
ensure the accurate  
Blood samples for all popular meters  
The needle retracts immediately after use  
Sterilization : Gamma  
Size: 21G, 23G, 26G



### FY0911 Pressure Activated Safety Lancets

Needle is fully shielded before and after use  
Pre-loaded, pressure activated design  
Ensures a good blood sample every time  
Ultra sharp tri-bevel needle or bi-bevel blade,  
high speed penetration  
Variety of needle specifications to suit most  
capillary blood requirements  
Sterilization : Gamma  
Size: 21G, 23G, 26G



### FY0912 Button Activated Safety Lancets

Needle tip is sealed with cap through injection moulding  
Pre-loaded, button operated design  
Single-use only, sterilised by gamma irradiation,  
cannot be reused  
High quality, ultra sharp tri-bevel needle,  
high speed penetration  
Variety of needle specifications to suit most capillary  
blood requirements  
Sterilization : Gamma  
Size: 21G, 23G, 26G



### FY0908 Lancing Device ( Releasable ) FY0909 Lancing Device

The user can install the lancet with the device easily.  
And the device will work with the lancet  
Smoothly without any blockage. The lancet won't  
fall on its own when to fire the device. Retracting  
automatically after firing, the launcher is safe and  
reliable.





**FY0913 / FY0914 / FY0915 Medical Razor**

Carbon steel or stainless steel Single or double blade All razors come with the highest quality stainless steel blades or carbon steel , Comfort-coated for smooth shaves



**Surgical Suture**

Cat. NO.	Type
FY1001	Chromic Catgut
FY1002	Plain Catgut
FY1003	Silk Braided
FY1004	Nylon Monofilament
FY1005	Polyester Braided
FY1006	Polypropylene
FY1007	Polyglycolic Acid

Thread Length: 45cm, 60cm, 75cm, 90cm, 150cm etc.  
 Thread Diameter: 8/0,7/0,6/0,5/0,4/0,3/0,2/0,1/0,1,2,3 etc.  
 Needle Length: 6mm,8mm,12mm,18mm,35mm,40mm,50mm etc.  
 Needle Curvature: Straight,1/2circle,1/2circle(double),1/4circle,1/4circle(double),3/8circle,3/8circle(double),5/8circle,loop round  
 Cross Section: round bodied,round bodied(heavy),curved cutting, curved cutting(heavy),reverse cutting,reverse cutting(heavy) etc.



**FY1502 Splints**

Aluminum plate Macromolecule foam material Nontoxic and nonbenzenoid  
 Free of effect of temperature and climate Type available: Fold, Unfold and Roll  
 Normal size: 5cm×12cm, 5cm×15cm, 7cm×11cm  
 11cm×44cm, 11cm×91.5cm, 11cm×100cm



**FY1501 Thermal Accident Blanket**

Material: PET Film  
 Waterproof, use to offset hypothermic reactions, reflective surface provides high visibility, reflects sun and heat for cooling in hot climates, space-age design retains up to 90% of body heat,Assists in prevention of post-accident shock and first aid application  
 Size: 132×210cm, 140×220cm, 160×210cm or other size per customer demand  
 Color: Gold, Sliver, Gold+Sliver



**EXAMINATION**

