



# PORTABLE DENTAL UNIT



## S202 Dental Unit

### Complete function Portable and practical

- ◆ S202 portable dental unit (economy) is made of high quality alloy aluminum.
- ◆ The box body is strong, beautiful, portable, and internally integrated with oil-free compressors, gas storage tanks, water, gas-electric control systems and other components, It is used in combination with folding chairs, cold lamps, etc.
- ◆ Weight:26kg
- ◆ Volume :510X360X720mm
- ◆ source(split) independent air
- ◆ Portable dental unitS202



Simple installation  
Easy operation  
Easy to carry

- ◆ The coldlight lamp adoptsanLEDsource,portable and durable,is matched with an adjustable three-di-mensional lamp holder,easy to be folded and is convenient to carry.

- ◆ S202 portable dental unit (standard) is made of high quality engineering plastics.



- ◆ Portable folding stool and portable folding chair



- ▶ Portable folding chair is made of high-strengthaluminumalloystructure, light and firm, comfortable and safe, bearing capacity of 135kg.

