




Diagnostic Consumables

Test Tubes

Round bottom plastic test tubes



- PS & PP material available

Cat No	Volume (mL)	Size	Packagaing
MF901111203	3	Φ12 x 60mm	500pc/bag, 10bag/ctn
MF901111205	5	Φ12 x 75mm	500pc/bag, 10bag/ctn
MF901111305	5	Φ13 x 75mm	500pc/bag, 10bag/ctn
MF901111610	10	Φ16 x 100mm	500pc/bag, 10bag/ctn
MF902111612	12	Φ16 x 102mm	500pc/bag, 10bag/ctn

Test Tubes Plug Caps

Flange type test tubes stopper

- PE material
- Neutral, White, Red, Mauve, Blue, Green colors available



Cat No	Color	Size	Packagaing
MF87010012	Various	Φ12mm	1000pc/bag, 20bag/ctn
MF87010013	Various	Φ13mm	1000pc/bag, 20bag/ctn
MF87010016	Various	Φ16mm	1000pc/bag, 20bag/ctn

Sample Cups

Used for hematology and coagulation analysis of whole blood specimens

- PS material

Cat No	Compatibility	Size	Packagaing
MF91000211	Hitachi	Φ12 x 38mm	5000pc/ctn
MF91000212	Beckman CX	Φ13 x 24mm	5000pc/ctn
MF91000213	Beckman	Φ14 x 25mm	10000pc/ctn

MF91000212

MF91000211



MF91000213




Diagnostic Consumables

Microscope Glass Slides (Pathology Grade)




Non-interleaved
Dim.: 24.4mm x 76.2mm (1"x3")
Thickness: 1.0 - 1.2mm

Material Element: Glass
Storage Requirements:
Prevent stress, dampproof, Normal temperature
Packing: 50pcs/box, 50 boxes/ctn
36 x 18 x 16cm



Code	Description	Corner
MF99117101	Fine polished edges, plain	90 Degrees
MF99117105	Fine polished edges, Frosted one end one side	90 Degrees
MF99117107	Fine polished edges, Frosted one end both sides	90 Degrees
MF99147102	Cut edges, plain	90 Degrees
MF99147105	Cut edges, Frosted one end one side	90 Degrees
MF99147107	Cut edges, Frosted one end both sides	90 Degrees
MF99127101	Fine polished edges, plain	45 Degrees
MF99127105	Fine polished edges, Frosted one end one side	45 Degrees
MF99127107	Fine polished edges, Frosted one end both sides	45 Degrees
MF99137101	Bevelled edges, plain	45 Degrees
MF99137105	Bevelled edges, frosted one end one side	45 Degrees
MF99137107	Bevelled edges, frosted one end both sides	45 Degrees
MF99127103	Single cavity, fine polished edges	45 Degrees
MF99127104	Double cavity, fine polished edges	45 Degrees





Diagnostic Consumables

Microscope Glass Slides (Standard Grade)

Thickness Available

Nr. 1	0.13.0.16mm
Nr. 1.5	0.16.0.19mm
Nr. 2	0.19.0.23mm

Material Element: Glass



Storage Requirements:

prevent stress, dampproof, normal temperature

*Special Dimensions available upon request

Code	Dim.	Packing
MF98111818	18x18	200pcs/vacuum pack, 2000pcs/inner box, 40000pcs/ctn
MF98112020	20x20	200pcs/vacuum pack, 2000pcs/inner box, 40000pcs/ctn
MF98112222	22x22	200pcs/vacuum pack, 2000pcs/inner box, 40000pcs/ctn
MD98112424	24x24	200pcs/vacuum pack, 2000pcs/inner box, 40000pcs/ctn
MF98112232	22x32	100pcs/vacuum pack, 1000pcs/inner box, 20000pcs/ctn
MF98112432	24x32	100pcs/vacuum pack, 1000pcs/inner box, 20000pcs/ctn
MF98112240	22x40	100pcs/vacuum pack, 1000pcs/inner box, 20000pcs/ctn
MF98112440	24x40	100pcs/vacuum pack, 1000pcs/inner box, 20000pcs/ctn
MF98112250	22x50	100pcs/vacuum pack, 1000pcs/inner box, 20000pcs/ctn
MF98112450	24x50	100pcs/vacuum pack, 1000pcs/inner box, 20000pcs/ctn
MF98112260	22x60	100pcs/vacuum pack, 1000pcs/inner box, 20000pcs/ctn
MF98112460	24x60	100pcs/vacuum pack, 1000pcs/inner box, 20000pcs/ctn

Code	Description	Corner
MF99117101	Fine polished edge, plain	90 Degrees
MF99117105	Fine polished edges, Frosted one end one side	90 Degrees
MF99117107	Fine polished edges, Frosted one end both sides	90 Degrees
MF99117102	Cut edges, plain	90 Degrees
MF99117105	Cut edges, Frosted one end one side	90 Degrees
MF99117107	Cut edges, Frosted one end one side	90 Degrees
MF99127101	Fine polished edges, plain	45 Degrees
MF99127105	Fine polished edges, Frosted one end both sides	45 Degrees
MF99127107	Fine polished edges, Frosted one end both sides	45 Degrees
MF99137101	Bevelled edges, plain	45 Degrees
MF99137105	Bevelled edges, Frosted one end both sides	45 Degrees
MF99137107	Bevelled edges, Frosted one end both sides	45 Degrees
MF99127103	Single cavity, Fine polished edges	45 Degrees
MF99127104	Double cavity, Fine polished edges	45 Degrees





Diagnostic Consumables

Microscope Glass Slides (Economy Grade)



Code	Description	Corner
MF993271015	Rough ground edges, plain	45 Degrees
MF993171025	Cut edges, plain	90 Degrees
MF993271055	Rough ground edges, frosted one end one side	45 Degrees
MF993471055	Cut edges, frosted one end one side	90 Degrees
MF993271075	Rough ground edges, frosted one end both side	45 Degrees
MF993171075	Cut edges, frosted one end both side	90 Degrees





Diagnostic Consumables

Glass Covers Slides (Economy Grade)



Thickness Available

Nr. 1 0.13.0.16mm

Nr. 1.5 0.16.0.19mm

Nr. 2 0.19.0.23mm

Special Dimensions available upon request

Glass covers are pieces of glass, typically used with microscope when viewing of specimens.

- Economy and Standard Grade available
- Vacuum packaging available

Code	Dim.	Packing
MF98211818	18x18	100pcs/box, 1000pcs/vacuum pack, 50000pcs/ctn
MF98212020	20x20	100pcs/box, 1000pcs/vacuum pack, 50000pcs/ctn
MF98212222	22x22	100pcs/box, 1000pcs/vacuum pack, 50000pcs/ctn
MF98212424	24x24	100pcs/box, 1000pcs/vacuum pack, 50000pcs/ctn
MF98212232	22x32	100pcs/box, 1000pcs/vacuum pack, 50000pcs/ctn
MF98212432	24x32	100pcs/box, 1000pcs/vacuum pack, 50000pcs/ctn
MF98212240	22x40	100pcs/box, 1000pcs/vacuum pack, 40000pcs/ctn
MF98212440	24x40	100pcs/box, 1000pcs/vacuum pack, 40000pcs/ctn
MF98212250	22x50	100pcs/box, 1000pcs/vacuum pack, 30000pcs/ctn
MF98212450	24x50	100pcs/box, 1000pcs/vacuum pack, 30000pcs/ctn
MF98212260	22x60	100pcs/box, 1000pcs/vacuum pack, 30000pcs/ctn
MF98212460	24x60	100pcs/box, 1000pcs/vacuum pack, 30000pcs/ctn



Diagnostic Consumables

Embedding Cassettes

EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
MF69010101P	Pink	Polyoxymethylene (P.O.M)

Biopsy-Square holes with removable lid,
P.O.M. Material

Packing: 250pcs/dispenser box, 8x250pcs/ctn



EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
MF69010104W	White	Polyoxymethylene (P.O.M)

Biopsy-Fine Square holes with Removable lid,
P.O.M. Material

Packing: 250pcs/dispenser box, 8x250pcs/ctn



EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
MF69010102W	White	Polyoxymethylene (P.O.M)

Histology-Strip holes with removable lid,
P.O.M. Material

Packing: 250pcs/dispenser box, 8x250pcs/ctn



EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
MF6903(04)0105P	Pink	Polyoxymethylene (P.O.M)

Biopsy-Fine Square holes with Detachable/Hinged lid,
P.O.M. Material

Packing:

MF6903: Cassettes and lids separately packed,
4x500 cassettes + 2x1000 lids/ctn

MF6904: Assembled cassettes and lids,
375pcs/dispenser box, 8x375pcs/carton



EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
MF69020103G	Green	Polyoxymethylene (P.O.M)

Standard-Round holes without lid,
P.O.M. Material

Packing: 500pcs/dispenser box, 8x500pcs/ctn



EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
MF6903(04)0106B	Blue	Polyoxymethylene (P.O.M)

Histology-Strip holes with Detachable/Hinged lid,
P.O.M. Material

Packing:

MF6903: Cassettes and lids separately packed,
4x500 cassettes + 2x1000 lids/ctn

MF6904: Assembled cassettes and lids,
375pcs/dispenser box, 8x375pcs/carton





Diagnostic Consumables

Embedding Cassettes

EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
------	-------	----------

MF6903(04)0107G	Green	Polyoxymethylene (P.O.M)
-----------------	-------	--------------------------

Histology-Strip holes with Detachable/Hinged lid, P.O.M. Material

Packing:

MF6903: Cassettes and lids separately packed, 4x500 cassettes + 2x1000 lids/ctn

MF6904: Assembled cassettes and lids, 375pcs/dispenser box, 8x375pcs/carton



EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
------	-------	----------


MF6903(04)0110Y	Purple	Polyoxymethylene (P.O.M)
-----------------	--------	--------------------------

Histology-Strip holes with Detachable Lid P.O.M. Material

Packing:

MF6903: Cassettes and lids separately packed, 4x500 cassettes + 2x1000 lids/ctn

MF6904: Assembled cassettes and lids, 375pcs/dispenser box, 8x375pcs/carton



EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
------	-------	----------

MF69020108Y	Yellow	Polyoxymethylene (P.O.M)
-------------	--------	--------------------------

Standard-Round holes without lid P.O.M. Material

Packing: 500pcs/dispenser box, 8x500pcs/ctn



EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
------	-------	----------


MF6903(04)01110	Purple	Polyoxymethylene (P.O.M)
-----------------	--------	--------------------------

Biopsy-Fine Square Holes with Detachable Lid P.O.M. Material

Packing:

MF6903: Cassettes and lids separately packed, 4x500 cassettes + 2x1000 lids/ctn

MF6904: Assembled cassettes and lids, 375pcs/dispenser box, 8x375pcs/carton



EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
------	-------	----------

MF69020109M	Purple	Polyoxymethylene (P.O.M)
-------------	--------	--------------------------

Standard-Square holes without lid P.O.M. Material

Packing: 500pcs/dispenser box, 8x500pcs/ctn



EMBEDDING CASSETTES ☺☺

Suitable for most Automatic Embedding Cassette labeling instruments

CODE	COLOR	MATERIAL
------	-------	----------

MF6903(04)0112M	Purple	Polyoxymethylene (P.O.M)
-----------------	--------	--------------------------

Histology-Strip Holes with Detachable Lid P.O.M. Material

Packing:

MF6903: Cassettes and lids separately packed, 4x500 cassettes + 2x1000 lids/ctn

MF6904: Assembled cassettes and lids, 375pcs/dispenser box, 8x375pcs/carton

