

# CentriVet™

GK Blood Glucose & Ketone  
Monitoring System

vs Competition



CentriVet™ GK



Abbott  
Precision Xtra®



NovaVet Blood  
Glucose/Ketone Meter



Wellion®Vet Belua



KetoVet TD-4235

Calibrated for Animals?	✓	X	✓	✓	✓
Sample Volume	β-Ketone: 1.2 μL	β-Ketone: 1.5 μL	β-Ketone: 1.2 μL	β-Ketone: 0.8 μL	β-Ketone: 0.7 μL
	Glucose: 0.8 μL	Glucose: 0.6 μL	Glucose: 0.8 μL	Glucose: 0.8 μL	X
Test Time	β-Ketone: 10 seconds	β-Ketone: 10 seconds	β-Ketone: 10 seconds	β-Ketone: 8 seconds	β-Ketone: 5 seconds
	Glucose: 5 seconds	Glucose: 5 seconds	Glucose: 4 seconds	Glucose: 6 seconds	X
Memory	450 tests with date, test type, and time	450 tests	400 tests	500 tests for Glucose 100 tests for Ketone	150 tests
Multiple Species?	✓	X	X	✓	X
Test Strip Ejector	✓	X	X	✓	X
Meter Cost	\$	\$\$\$	\$\$	\$\$\$	\$
Strip Cost	\$	\$\$\$\$	\$\$\$	\$\$\$	\$\$\$

All trademarks are the property of their respective owners. Competitor specifications excerpted from user manuals and test strip package inserts that were current as of March 2016.

ACON Laboratories does not guarantee the accuracy or completeness of the information and data, makes no representations of warranties, express or implied, and assumes no liability associated with use or reliance thereof. Pricing determined based on currency rates as of March 2016. For informational purposes only.

© 2016 ACON Laboratories, Inc. San Diego, CA 92121