

CareX ObxTable



Specifications:

- Bed dimensions: L1525*W730*H900mm (± 20 mm)
- Gas spring-assisted backrest adjustment from 0- 75° $\pm 5^\circ$
- Sliding, retractable foot section cushion.
- Manual pelvic tilt (gynecological examination): 5°
- Integrated a pair of retractable stirrups.
- 5 mono block polyester ABS drawers (3 with bilateral opening 2 open from the front).
- Bed frame resistant to corrosion, anti-rust.
- Retractable foot stand which is from the front.
- Paper roll holder under the back section.

Explore more here:



Available in 8 colours:

#1C1B1A	#303f66	#AA9F96	#98968F
Black	Blue	Beige	Smoky Cashmere
#A09880	#CBB680	#9FB4BF	#69312E
Creamy Latte	Lemon Meringue	Blue Skies	Fine Wine

*Other colors, available on demand

✉ general@ajjhealthcare.com

☎ +65 6988 4251