

# SURGICAL POSITIONING PAD




**AJJ+**  
HEALTHCARE MANAGEMENT PTE LTD


✉ [general@ajjhealthcare.com](mailto:general@ajjhealthcare.com)

☎ +65 6908 4251

📍 8 Commonwealth Lane #02-04 Singapore 149555



## | 仰卧位手术体位 | Resupine Positioning



### 配置项目

#### ▶ Configuration Items

名称 name	型号 model	规格 specifications	数量 number	名称 name	型号 model	规格 specifications	数量 number
圆形头圈 Donut head pad	MT-2-2	20*7*5	1	腕关节固定带 Wrist Protectors	MT-12-2	32*5*1	2
手术组合床垫 Surgery combined pad	MT-27-1	50*23*1.5	1	斜形垫 slanting shaped pad	MT-18-3	30*22*5/1	2
	MT-27-2	50*50*1.5		半圆形垫 Universal Dome Shaped pad	MT-21-11	50*15*7	1
	MT-27-3	50*100*1.5		跟骨垫 Heel pad	MT-23-5	18*10*7	2
凹形体位垫 Contoured Arm/Leg pad	MT-22-1	50*15*4	2				

## | 俯卧位手术体位 | Prostrate Positioning



### 配置项目


#### ▶ Configuration Items

名称 name	型号 model	规格 specifications	数量 number
俯卧位头垫 Prostrate head pad	MT-1-1	28*24*14	1
俯卧位垫 Prostrate position pad	MT-14-2	62*47*17	1
凹形体位垫 Contoured Arm/Leg pad	MT-22-1	50*15*4	3
半圆形垫 Universal Dome Shaped pad	MT-21-3	50*15*10	1





## 侧卧位手术体位 | Lateral Positioning



### 配置项目

#### Configuration Items

名称 name	型号 model	规格 specifications	数量 number	名称 name	型号 model	规格 specifications	数量 number
开放式头圈 Horseshoe head pad	MT-2-13	20*7*5	1	体部固定带 Body protectors	MT-20-1	60*7*1	1
侧卧位垫 Lateral position pad	MT-13-1	70*50*17	1	床垫 Operation table pad	MT-27-3	50*100*1.5	1
半圆形垫 Universal Dome Shaped pad	MT-21-5	50*10*4	1	隧道垫 Tunnel pad	MT-26-1	64*46*25	1
凹形体位垫 Contoured Arm/Leg pad	MT-22-1	50*15*4	2	肘部固定垫 Elbow pads	MT-11-1	40*15*1.5	1

## 截石位手术体位 | Lithotomy Positioning



### 配置项目

#### Configuration Items

名称 name	型号 model	规格 specifications	数量 number	名称 name	型号 model	规格 specifications	数量 number
心形头垫 Heart shaped head pad	MT-2-6	29*28.5*7/3.3	1	臀垫 Hip pad	MT-17-4	31*32.5*5	1
通用方垫 Universal Shaped pad	MT-29-17	40*24*1.5	2	手术床垫 Operation table pad	MT-27-1	50*23*1.5	1
凹形体位垫 Contoured Arm/Leg pad	MT-22-1	50*15*4	2	手术床垫 Operation table pad	MT-27-3	50*100*1.5	1

| 俯卧位头垫 | Prostrate Head Pad 



XHT-1-1 使用示范


**用途功能说明**  
 Function Description

适用于俯卧位手术，用于头面部的支撑保护，方便使用麻醉管路，使病人保持呼吸畅通。(XHT-1-2也适用于眼科病人术后康复)

Applicable to the prone position surgery ,to support pratient' s facial part, easy to use anesthesia circuit, keep patient' s breathe smoothly. (XHT - 1-2 can also be applied to ophthalmic patients postoperative rehabilitation)




名称 name 俯卧位头垫 Prostrate Head Pad 型号 model MT-1-1 尺寸 size 28\*24\*14 (cm)



名称 name 俯卧位头垫 Prostrate Head Pad 型号 model MT-1-2 尺寸 size 25\*20\*10 (cm)




名称 name 俯卧位头垫 Prostrate Head Pad 型号 model MT-1-6 (儿童型) (pediatric) 尺寸 size 20\*16\*11 (cm)



名称 name 俯卧位头垫 Prostrate Head Pad 型号 model MT-1-3 尺寸 size 23\*20\*11 (cm)



名称 name 俯卧位头垫 Prostrate Head Pad 型号 model MT-1-4 尺寸 size 23\*20\*6 (cm)



名称 name 俯卧位头垫 Prostrate Head Pad 型号 model MT-1-5 (一对) (a pair) 尺寸 size 18\*22\*5 (cm)

名称 name 俯卧位头垫 Prostrate Head Pad 型号 model MT-1-7 (带镜面可调节) (With a mirror can be adjusted) 尺寸 size 23\*20\*14/11 (cm)

| 圆形头圈 | Donut Head Pad 



型号 model	尺寸 (cm) size (cm)
MT-2-1	20*7*7
MT-2-2	20*7*5
MT-2-3 (儿童型) (pediatric)	14*5*3
MT-2-4 (儿童型) (pediatric)	10*4*2.5
MT-2-5 (新生儿型) (Neonatal)	8*3*2
MT-2-7 (儿童型) (pediatric)	14*7*2.5
MT-2-9	18*10*3.5
MT-2-10	20*12*3.5
MT-2-11	22*15*3.5

名称 name	型号 model	尺寸 (cm) size (cm)
凹边圆形头圈 concave round head pad	MT-8-2	20*9*5 (-1.5)
心形头垫 heart shaped head pad	MT-2-6	29*28.5*7/3.3

**用途功能说明**  
 Function Description

适用于仰卧位手术，用于支撑、稳定、保护头枕部，避免压疮的产生，一般应用于脑外科、五官科手术。

Applied to supine surgery, to support and protect patient' s head and neck to avoid pressure sores, usually used in brain surgery, ENT surgery.





## 开放式头圈 | Horseshoe Head Pad

### 用途功能说明

#### Function Description

适用于仰卧、侧卧、俯卧位手术，用于支撑保护头面部，便于全麻插管的使用。

Applicable to the supine position, lateral position and prone position surgery, used to support and protect the head and face, easy to use anesthesia intubation.



型号 model	尺寸 (cm) size (cm)
MT-2-12	20*7*7



型号 model	尺寸 (cm) size (cm)
MT-2-13	20*7*5
MT-2-14 (儿童型) (pediatric)	14*5*3

型号 model	尺寸 (cm) size (cm)
MT-2-15 (儿童型) (pediatric)	10*4*2.5
MT-2-16 (新生儿型) (Neonatal)	8*3*2



型号 model	尺寸 (cm) size (cm)
MT-2-17	20*7*5
MT-2-18 (儿童型) (pediatric)	14*5*3

型号 model	尺寸 (cm) size (cm)
MT-2-19 (儿童型) (pediatric)	9*4*2.5

## 通用头垫 | Universal Flat Head Pad



型号 model	尺寸 (cm) size (cm)
MT-3-1	25*23*6
MT-3-2	25*18*3
MT-3-3	45*15*10

### 用途功能说明

#### Function Description

适用于仰卧位手术缓解头部压力，也可用于侧卧位手术同辅助支架配合使用，加强对腰部、膝部及会阴部等部位的保护。

Suitable for supine surgery to alleviate the pressure to the head, also can be used in lateral position operation with the auxiliary bracket used to strengthen the protection of the waist, knee and other parts of the perineum.



## 碗形头圈 | Bowl Shaped Head Pad



型号 model	尺寸 (cm) size (cm)
MT-4-1	8*2/20*5
MT-4-2 (儿童型) (Pediatric)	6*1.5/14*3

### 用途功能说明

#### Function Description

适用于仰卧位手术，用于支撑、稳定、保护头枕部，避免压疮的产生，一般应用于脑外科、五官科手术。

The operation is suitable for lying on the back, which is used for supporting, stabilizing and protecting the head rest pillow and avoiding the generation of pressure sore.

## 肩部保护垫 | Shoulder protective Pad



型号 model	尺寸 (cm) size (cm)
MT-5-1 (一对) (a pair)	16*15*2.5

### 用途功能说明

#### Function Description

保护肩部时，适用于头高脚低体位。

for the protection of the shoulder, and apply to reverse trendenbug position.

## 凹形枕垫 | Concave head Pad



型号 model	尺寸 (cm) size (cm)
MT-6-1	18*18*8/4

### 用途功能说明

#### Function Description

适用于仰卧位手术，稳定头枕部；也可用于侧卧位手术，肩部及膝部的保护。

Suitable for supine operation, stable occipital; can also be used for bending operation of lateral decubitus, shoulder hip and knee protection.

## 眼科头垫 | Ophthalmic Head Pad



MT-7-2 产品图

型号 model	尺寸 (cm) size (cm)
MT-7-2	22*19*6

### 用途功能说明

#### Function Description

用于稳定患者头部，适用于眼科、耳鼻喉科及颌面、整形外科等手术。

For the stability of the patient's head, apply to the eye, Department of ENT and maxillofacial, plastic surgery and other operations.




MT-7-3 产品图

型号 model	尺寸 (cm) size (cm)
MT-7-3	29*28.5*7/3.3



## | 甲状腺垫 | Thyroid Pad




MT-9-1 产品图

型号 model	尺寸 (cm) size (cm)
MT-9-1	45*35*4

### 用途功能说明

#### ▶ Function Description

适用于头颈部及甲状腺手术，充分暴露手术视野，加强对颈后部的支撑和保护，避免手术过程中的不适应症。



MT-9-2 产品图

型号 model	尺寸 (cm) size (cm)
MT-9-2	48*39*9

Suitable for head and neck and thyroid surgery, fully exposed the operation field of vision, to strengthen the support and protection of the back of the neck, to avoid the operation process of the indications.

## | 上肢垫 | Universal Armboard Pad



型号 model	尺寸 (cm) size (cm)
MT-10-1	40*16*1

MT-10-2	60*16*1
---------	---------

MT-10-3	60*12.5*2
---------	-----------

MT-10-4	40*12.5*2
---------	-----------

### 用途功能说明

#### ▶ Function Description

通用于手术床上臂支架，保护肘部和臂部。

General to the upper arm of the operating arm support, the protection of the elbow and arm.



## | 肘部护垫 | Elbow Pads



型号 model	尺寸 (cm) size (cm)
MT-11-1	40*15*1.5

### 用途功能说明

#### ▶ Function Description

保护肘部，避免损伤。  
Protect elbows to avoid injury.

## | 腕关节固定带 | Wrist Protectors



型号 model	尺寸 (cm) size (cm)
MT-12-1	25*5*1

MT-12-2	32*5*1
---------	--------

### 用途功能说明

#### ▶ Function Description

用于手术中，对腕关节的固定，避免腕关节的神经损伤，预防压疮。

For the operation, the fixation of the wrist joint, to avoid the nerve injury of the wrist joint, prevention of pressure ulcers.

## | 侧卧位垫 | Lateral Positioning Pad



型号 model	尺寸 (cm) size (cm)
MT-13-1 (海绵底托) (Sponge tap)	70*50*17

### 用途功能说明

#### ▶ Function Description

适用于所有侧卧位手术，减轻肩膀和上肢的压力并保护神经组织。

Applicable to all lateral decubitus operation, reduce the pressure of shoulder and upper limb and protect the nerve tissue.





## | 俯卧位垫 | Prostrate Pad



MT-14-1 产品图

**型号**  
model

**尺寸 (cm)**  
size (cm)

MT-14-1  
(海绵底托)  
(Sponge base) 63\*46\*18



MT-14-2 产品图

**型号**  
model

**尺寸 (cm)**  
size (cm)

MT-14-2  
(海绵底托)  
(Sponge base) 62\*47\*17

### 用途功能说明

#### ▶ Function Description

适用于俯卧位手术，减轻胸部压力，保持良好的呼吸循环，多用于颈椎、腰椎、脊柱手术。

Suitable for surgery in the prone position, reduce the chest pressure force, to maintain a good respiratory cycle, and more in the cervical spine, lumbar spine, spine surgery.



**型号**  
model

**尺寸 (cm)**  
size (cm)

MT-14-3  
电动  
(electric) 66\*(30-40)\*29

MT-14-4  
手动  
(electric) 66\*(30-40)\*22

### 用途功能说明

#### ▶ Function Description

适用于俯卧位手术，可调节式俯卧位垫，可配置与专业脊柱床使用。

Suitable for surgery in the prone position, adjustable prone pad can be configured with professional spine bed.



MT-14-5 产品图

**型号**  
model

**尺寸 (cm)**  
size (cm)

MT-14-5  
(一对)  
(a pair) 71/58\*20\*7.6

### 用途功能说明

#### ▶ Function Description

适用于俯卧位手术，可调节式俯卧位垫，可配置与专业脊柱床使用。

Suitable for surgery in the prone position, adjustable prone pad can be configured with professional spine bed.



MT-14-7/8/9 产品图


**型号**  
model

**尺寸 (cm)**  
size (cm)

MT-14-7  
(一对)  
(a pair) 30\*16\*13/9

MT-14-8  
(一对)  
(a pair) 40\*16\*13/9

MT-14-9  
(一对)  
(a pair) 50\*16\*13/9



MT-14-10 产品图

**型号**  
model

**尺寸 (cm)**  
size (cm)

MT-14-10  
(一对)  
(a pair) 50\*12\*10/4



MT-14-11 产品图

**型号**  
model

**尺寸 (cm)**  
size (cm)

MT-14-11  
(分体式)  
(split type) 675\*500\*11

### 用途功能说明

#### ▶ Function Description

适用于俯卧位手术，可调节式俯卧位垫，可配置与专业脊柱床使用。

Suitable for surgery in the prone position, adjustable prone pad can be configured with professional spine bed.

## | 胸髋垫 | Chest-Hipbone Pad



**型号**  
model

**尺寸 (cm)**  
size (cm)

MT-15-1 50\*20\*7



MT-15-2 7\*4\*4

**型号**  
model

**尺寸 (cm)**  
size (cm)

MT-15-3 36\*23\*3.8

### 用途功能说明

#### ▶ Function Description

适用于俯卧位手术，支撑保护胸髋部，保持良好的呼吸循环，避免产生压疮；也可用于仰卧位手术中，支撑肩部、胸部。

Suitable for surgery in the prone position, supporting and protecting the chest hip, maintain good respiratory circulation, avoid pressure sores; also can be used in the supine position during the operation, supporting shoulder and chest.

## | 骨科柱垫 | Orthopaedic Pad



**型号**  
model

**尺寸 (cm)**  
size (cm)

MT-16-1 27\*20\*1.5

### 用途功能说明

#### ▶ Function Description

配合骨科牵引床使用，置于会阴部，避免压疮产生。

With the Department of orthopedics traction bed, placed in the perineal, to avoid pressure ulcers.

| 臀垫 | Hip Pad



型号  
model


尺寸 ( cm )  
size ( cm )

MT-17-1 50\*40\*1.5

用途功能说明

▶ Function Description

保护骶尾部，避免压疮的产生。  
To protect sacrococcygeal region from pressure sores



型号  
model


尺寸 ( cm )  
size ( cm )

MT-17-3 40\*40\*1.8

用途功能说明

▶ Function Description

保护骶尾部，避免压疮的产生。  
To protect sacrococcygeal region from pressure sores



型号  
model

尺寸 ( cm )  
size ( cm )


MT-17-2 43\*35\*5

用途功能说明

▶ Function Description

保护骶尾部，避免压疮的产生，主要用于肛周预防保健、肛周手术前后、肛周及臀部的疾患护理。

Protection of the sacral tail, avoid the generation of pressure ulcers, mainly applicable to the anus week to prevent health care, the anus week operation, anus week and buttocks of the disease care.



型号  
model


尺寸 ( cm )  
size ( cm )

MT-17-4 31\*32.5\*5

用途功能说明

▶ Function Description

保护骶尾部，避免压疮的产生。  
To protect sacrococcygeal region from pressure sores



型号  
model

尺寸 ( cm )  
size ( cm )

MT-17-5 44\*42\*3

用途功能说明

▶ Function Description

保护骶尾部，避免压疮的产生。  
To protect sacrococcygeal region from pressure sores



型号  
model

尺寸 ( cm )  
size ( cm )

MT-17-6 50\*50\*3

用途功能说明

▶ Function Description

保护骶尾部，避免压疮的产生。  
To protect sacrococcygeal region from pressure sores

| 斜形垫 | Prostrate Pad(Slanting Shaped Pads)



型号  
model

尺寸 ( cm )  
size ( cm )

MT-18-1	30*18*7/3
MT-18-2	30*18*11/8
MT-18-3	30*22*5/1
MT-18-4	25*15*6/2
MT-18-5	50*12*10/4
MT-18-6	50*16*13/9
MT-18-7	25*21*9/1
MT-18-10	30*13*5/2

用途功能说明

▶ Function Description

适用于仰卧位、侧卧位和俯卧位手术，主要用于仰卧位手术时腰部、肝部、脾部的支撑和保护。

Applicable to the supine position, lateral position and prone position surgery, mainly used in the supine position during surgery, liver, spleen, waist support and protection.

| 柱形垫 | Pillar Shaped Pad



型号  
model

尺寸 ( cm )  
size ( cm )

MT-19-1	50*12*10
MT-19-2 (儿童型)	10*5*5
MT-19-3	40*12*10
MT-19-4	30*12*10
MT-19-5	20*12*10

用途功能说明

▶ Function Description

适用于仰卧位和俯卧位手术，用于支撑保护膝踝部、肩部；组合应用于仰卧位手术，支撑保护肩部和髋部；也可应用于侧卧位手术，对腋下的支撑和保护。

Applicable to the supine and prone position surgery, used to support and protect the knee and ankle, shoulder; combination in supine operation, supporting and protecting the shoulders and hips; also can be used in lateral position operation, support and protection of the armpit.

| 体部固定带 | Body Protectors



型号  
model

尺寸 ( cm )  
size ( cm )

MT-20-1 60\*7\*1

用途功能说明

▶ Function Description

固定身体和关节，避免普通固定带对身体的损伤。

To fix the body and joints, to avoid the damage to the body.



## 半圆形垫 | Universal Dome Shaped Pad



型号 model	尺寸 (cm) size (cm)
MT-21-1 (儿童型) (pediatric)	30*10*7
MT-21-2	50*10*7
MT-21-3	50*15*10
MT-21-4	50*10*10
MT-21-5	50*10*4
MT-21-6	30*10*4
MT-21-7	25*10*7
MT-21-8	20*10*7
MT-21-9	30*15*7
MT-21-10	40*15*7
MT-21-11	50*15*7
MT-21-12	40*10*7
MT-21-13	50*18*15
MT-21-14	20*15*7

### 用途功能说明 Function Description

适用于仰卧位和俯卧位手术，用于支撑保护膝踝部、肩部；组合应用于仰卧位手术，支撑保护胸部和髋部；也可应用于侧卧位手术，对腋下的支撑和保护。

Applicable to the supine and prone position surgery, used to support and protect the knee and ankle, shoulder; combination in supine surgery, chest and hip support; also can be used in lateral position operation, support and protection of the armpit.

## 凹形体位垫 | Contoured Arm/Leg Pad



型号 model	尺寸 (cm) size (cm)
MT-22-1	50*15*4
MT-22-2	60*15*4
MT-22-3	40*15*4

### 用途功能说明 Function Description

适用于俯卧位手术，支撑保护膝盖，也可用于上臂的支撑保护。

Suitable for surgery in the prone position, support and protect the knee, can also be used to support the protection of the upper arm.

## 跟骨垫 | Heel Pad




型号 model	尺寸 (cm) size (cm)
MT-23-1 (一对) (a pair)	14*11*5.5

### 用途功能说明 Function Description

保护踝部、足跟部，也可用于骨科跟骨牵引，避免形成压疮。

Apply to protect patient's ankle and heel, also can be used in orthopedics and skeletal traction to avoid pressure sores.



型号 model	尺寸 (cm) size (cm)
MT-23-2 (一对) (a pair)	12*9.5*4.5

### 用途功能说明 Function Description

保护踝部、足跟部，也可用于骨科跟骨牵引，避免形成压疮。

Apply to protect patient's ankle and heel, also can be used in orthopedics and skeletal traction to avoid pressure sores.



型号 model	尺寸 (cm) size (cm)
MT-23-5 (一对) (a pair)	18*10*7

### 用途功能说明 Function Description

保护踝部、足跟部，也可用于骨科跟骨牵引，避免形成压疮。

Apply to protect patient's ankle and heel, also can be used in orthopedics and skeletal traction to avoid pressure sores.



型号 model	尺寸 (cm) size (cm)
MT-23-6 (一对) (a pair)	15*8*4.5

### 用途功能说明 Function Description

保护踝部、足跟部，也可用于骨科跟骨牵引，避免形成压疮。

Apply to protect patient's ankle and heel, also can be used in orthopedics and skeletal traction to avoid pressure sores.



型号 model	尺寸 (cm) size (cm)
MT-23-3 (一对) (a pair)	45*17*1

### 用途功能说明 Function Description

保护踝部、足跟部，也可用于骨科跟骨牵引，避免形成压疮。

Apply to protect patient's ankle and heel, also can be used in orthopedics and skeletal traction to avoid pressure sores.

## 脚踝固定带 | Ankle Protectors



型号 model	尺寸 (cm) size (cm)
MT-12-2	32*5*1


### 用途功能说明 Function Description

用于手术中，对腕关节的固定，避免腕关节的神经损伤，预防压疮。

For the operation, the fixation of the wrist joint, to avoid the nerve injury of the wrist joint, prevention of pressure ulcers.



## | 隧道垫 | Tunnel Pad




型号 model	尺寸 (cm) size (cm)
MT-26-1	64*46*25

### 用途功能说明

#### ▶ Function Description

适用于侧卧位手术，保护下肢体位。  
Apply to protect lower limb in lateral surgery.

## | 手术床垫 | Operation Table Pads



型号 model	尺寸 (cm) size (cm)
MT-27-1	50*23*1.5
MT-27-2	50*50*1.5
MT-27-3	100*50*1.5
MT-27-4 (带凹槽) (With cutout)	50*50*1.5
MT-27-5 (带凹槽) (With cutout)	100*50*1.5


型号 model	尺寸 (cm) size (cm)
MT-27-6 (一体) (One)	180*50*1.5
MT-27-7 (一体) (One)	170*50*1.5
MT-27-8 (一体) (One)	150*50*1.5
MT-27-9 (一体) (One)	120*50*1.5

### 用途功能说明

#### ▶ Function Description

配合手术床使用，保护皮肤避免压疮和神经损伤。  
Used on operation tables to protect the skin from pressure sores and avoid nerve damage.

## | 通用方垫 | Universal Square Pad



型号 model	尺寸 (cm) size (cm)
MT-29-1	50*40*1.5
MT-29-2	50*30*1.5
MT-29-3	50*20*1.5
MT-29-4	50*80*1.5
MT-29-5	50*70*1.5
MT-29-6	15*15*1.5
MT-29-7	15*75*1.5
MT-29-8	15*60*1.5
MT-29-9	15*50*1.5
MT-29-10	15*40*1.5

型号 model	尺寸 (cm) size (cm)
MT-29-11	15*30*1.5
MT-29-12	16*60*1
MT-29-13	16*50*1
MT-29-14	16*40*1
MT-29-15	16*30*1
MT-29-16	30*20*1.5
MT-29-17	40*24*1.5
MT-29-18	30*20*2.5
MT-29-19	50*20*2.8
MT-29-20	40*40*1.5
MT-29-21	40*30*1.5

## | 多用途垫 | Multipurpose Pad



型号 model	尺寸 (cm) size (cm)
MT-28-1	10*4*2
MT-28-2	33*7*5
MT-28-3	25*12.5*5
MT-28-4	33*12*6
MT-28-5	38*12.5*6.3
MT-28-6	50*12.5*7
MT-28-7	36*12.5*7
MT-28-8	17*8*4.5

型号 model	尺寸 (cm) size (cm)
MT-28-9	25*13*3
MT-28-10	35*12*3
MT-28-11	46*11*3
MT-28-12	36*23*4
MT-28-13	20*14*2
MT-28-15	13*5*5
MT-28-29	50*30*2

## | 积雪包 | Snow Pack



### 使用方法： Using method:

积雪包垫在患者身体下，摆成理想体位后，将积雪包气孔连接中心负压吸引接口或电动负压引流瓶的接口，抽取负压，15秒钟即可固定患者体位

Put the snow pack under the body of patient, place it in an ideal position, connection the vents of snow pack with the interface of center negative pressure suction or the interface of electric negative pressure suction drainage bottle extract negative pressure, the position can be fixed after 15 seconds.



型号 model	尺寸 (cm) size (cm)
MT-A-1	100*73
MT-A-2	100*90



### 用途功能说明 Function Description

适用于侧卧位手术  
Suitable for lateral position surgery

型号	尺寸 (cm)
MT-A-3	40*50

### 用途功能说明 Function Description

固定腰部或保护足部  
Fix the waist or protect the foot








型号 model	尺寸 ( cm ) size ( cm )
MT-A-4	165*72
MT-A-5	180*72

### 用途功能说明 Function Description

全身垫  
Body Pad



型号 model	尺寸 ( cm ) size ( cm )
MT-B-1	100*73

### 用途功能说明 Function Description

适用于侧卧位手术  
Suitable for lateral  
position surgery



型号 model	尺寸 ( cm ) size ( cm )
MT-C-1	72*88
MT-F-1	100*88

### 用途功能说明 Function Description

甲状腺手术时，快速固定颈部位置  
Thyroid gland operation,  
fast fixed neck position



型号 model	尺寸 ( cm ) size ( cm )
MT-E-1	100*73
MT-E-2	100*90

### 用途功能说明 Function Description

适用于侧卧位手术  
Suitable for lateral  
position surgery

## 注意事项

1. 轻拿轻放；
2. 避免接触坚硬和锐力物体；
3. 不要用腐蚀性强和含碘的消毒清洁剂清洗垫子表面；
4. 平时存放要放平整，避免阳光照射，防尘；
5. 为加强舒适性，建议在侧卧位及俯卧位手术时，在体位垫上平铺一层手术单巾；
6. 避免用力将手术体位垫推挤于患者身体下方，保证垫子同身体接触面的平整；
7. 使用产品前请仔细阅读各部位产品的使用方法说明；
8. 建议如果使用时间过长（特别是俯卧位手术），每三小时应将易受压部位放松按摩 5 分钟。

### Matters needing attention

1. Should be used gently to avoid the drop hit.
2. Keep away from sharp objects.
3. Do not clean the pad by strong corrosive detergent or iodine detergent.
4. Pads should be stored properly away from direct solar radiation and dust.
5. To be more comfortable, a towel or a sheet can be covered on the pad in lateral and prostrate surgery.
6. Do not foist the pad under patient body toughly. Make sure the touching surface between patient body and the pad smooth.
7. Read the operation and function description carefully before using the positioning pads.
8. Release the pad 5 minutes after 3 hours continuous use.

## 禁忌症

### contraindications

- 1) 有透气性要求的体表损伤部位禁用；  
1) 1. Do not directly contact with injured skin.
- 2) 对聚氨酯材料接触性过敏的患者禁用；  
2) People who is allergic to polyurethane materials are prohibited to use the pad.
- 3) 极度肥胖且需应用俯卧位手术者禁用。  
3) Extremely obese people are prohibited to use the pad in prostrate surgery.

## 清洁及消毒方法

### Cleaning and disinfection methods

手术体位垫是不吸水的，对于垫子表面的污物，可用清水或清洁剂清洗，也可以用酒精消毒；

- ▶ Pad can be cleaned and disinfected by water, detergent and alcohol.

不建议用含碘消毒剂擦洗；

- ▶ Do not clean the pad by iodine detergent.

不能高温高压消毒和熏蒸；

- ▶ Do not sterilize and fumigate the pads by high temperature and high pressure.

不要在消毒液中长时间浸泡；

- ▶ Do not soak the pads in disinfectant for a long time.

避免用紫外线照射消毒。

- ▶ Avoid long time exposure to the ultraviolet rays.